



AFRL-RH-WP-TR-2020-0063

**COG-PACK™ DASHBOARD
USERS MANUAL**

**Mr. Allen W. Dukes
2Lt. Matt T. Rommel
Airman Biosciences Division**

**Mr. Courtney J. Downs
Mr. Scott J. Duberstein
Mr. Ethan B. Blackford
Ball Aerospace & Technologies**

**May 2020
Interim Report**

DISTRIBUTION A. Approved for public release:

**AIR FORCE RESEARCH LABORATORY
711TH HUMAN PERFORMANCE WING,
AIRMAN SYSTEMS DIRECTORATE,
WRIGHT-PATTERSON AIR FORCE BASE, OH 45433
AIR FORCE MATERIEL COMMAND
UNITED STATES AIR FORCE**

NOTICE AND SIGNATURE PAGE

Using Government drawings, specifications, or other data included in this document for any purpose other than Government procurement does not in any way obligate the U.S. Government. The fact that the Government formulated or supplied the drawings, specifications, or other data does not license the holder or any other person or corporation; or convey any rights or permission to manufacture, use, or sell any patented invention that may relate to them.

This report was cleared for public release by 88th Air Base Wing Public Affairs Office and is available to the general public including foreign nationals. Copies may be obtained from the Defense Technical Information Center (DTIC) (<http://www.dtic.mil>).

AFRL-RH-WP-TR-2020-0063 HAS BEEN REVIEWED AND IS APPROVED FOR PUBLICATION.

DOUGLAS L. FRANCK, DR-III, DAF
Chief, Performance Optimization Branch
Airman Biosciences Division

MCKINLEY.RICHARD.A.126223
RD.A.1262239710
Digitally signed by
MCKINLEY.RICHARD.A.126223
9710
Date: 2020.09.03 10:08:15
-04'00'

R. ANDY MCKINLEY, PhD
Core Research Area Lead
Performance Optimization Branch
Airman Biosciences Division

This report is published in the interest of scientific and technical information exchange, and its publication does not constitute the Government's approval or disapproval of its ideas or findings.

REPORT DOCUMENTATION PAGE				<i>Form Approved</i> OMB No. 0704-0188	
The public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number. PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.					
1. REPORT DATE (DD-MM-YY) 04 08 20		2. REPORT TYPE Interim		3. DATES COVERED (From - To) January 2019 to August 2020	
4. TITLE AND SUBTITLE COG-PACK™ Dashboard Users Manual				5a. CONTRACT NUMBER FA8650-16-C-6610	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER 62202F	
6. AUTHOR(S) COURTNEY J. DOWNS+, ALLEN W. DUKES*, ETHAN B. BLACKFORD+, SCOTT J. DUBERSTEIN+, MATT T. ROMMEL*				5d. PROJECT NUMBER 5329	
				5e. TASK NUMBER 08	
				5f. WORK UNIT NUMBER H0KG (53290819)	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) +Ball Aerospace & Technologies 2675 Presidential Drive, Fairborn OH 45324				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES) *Air Force Materiel Command Air Force Research Laboratory 711 Human Performance Wing Airman Systems Directorate Airman Biosciences Division Performance Optimization Branch Cognitive Neuroscience Section Wright-Patterson AFB, OH 45433				10. SPONSORING/MONITORING AGENCY ACRONYM(S) 711HPW/RHBC	
				11. SPONSORING/MONITORING AGENCY REPORT NUMBER(S) AFRL-RH-WP-TR-2020-0063	
12. DISTRIBUTION/AVAILABILITY STATEMENT DISTRIBUTION A. Approved for public release,					
13. SUPPLEMENTARY NOTES 88ABW-2020-2718, cleared 28 August 2020					
14. ABSTRACT The design of the COG Pack Dashboard allows a user to set up, configure, and manage multiple Sensors by providing a standard look and feel for each Sensor. Generally, when manufacturers develop hardware sensors, they provide a custom application designed to only work with their particular device. From an integration standpoint, this provides a high barrier of entry for researchers and users when an experiment requires multiple disparate sensors. Each independent application often provides unique functionality for actions such as saving files, managing sensor state connectivity, and observing the data-streams in real-time. One of the objectives of the COG Pack Dashboard is to simplify these steps and consolidate them into a single interface. Additionally, by leveraging a database back-end for data storage, COG Pack automatically manages data timestamps and organization. By not requiring a user to execute multiple mouse actions and keyboard inputs for file management, COG Pack streamlines this process and reduces the potential for 'fat-fingered' input errors. The goal of this document is to provide a user or researcher with the instructions they need to work with the COG Pack Dashboard. This document contains descriptions of the processes and procedures required to interact with the supported Sensors and Signals.					
15. SUBJECT TERMS COG Pack™, Dashboard, User Interface, UJI, Manual					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT: SAR	18. NUMBER OF PAGES 18	19a. NAME OF RESPONSIBLE PERSON (Monitor) Armando Soto 19b. TELEPHONE NUMBER (Include Area Code) N/A
a. REPORT U	b. ABSTRACT U	c. THIS PAGE			



COG PACK™ DASHBOARD

COURTNEY J. DOWNS, ALLEN W. DUKES, ETHAN B.
BLACKFORD, SCOTT J. DUBERSTEIN, MATT T.
ROMMEL

User's Manual

TABLE OF CONTENTS

1.0	OVERVIEW	1
2.0	STARTING THE COG PACK DASHBOARD	2
3.0	ASSIGNING SENSORS TO PARTICIPANTS.....	2
3.1	HSI PHYSIO Configuration	3
3.2	Profusa Lumee Configuration	3
3.3	Intelligent Optical Systems (IOS) MASES Configuration	4
3.4	Smart Eye Pro Configuration.....	5
4.0	ATTACHING ADDITIONAL SIGNAL PROCESSING STEPS.....	8
4.1	Attaching a Threshold Calculator.....	8
4.2	Attaching a “Core Temperature Estimator”	10
5.0	GRAPHS AND VISUALIZATIONS.....	11
6.0	LOADING AND SAVING CONFIGURATIONS	12
7.0	CONCLUSION	13

LIST OF FIGURES

Figure 1 - COG Pack Desktop Shortcut.....	2
Figure 2 - COG Pack Initialization Screen	2
Figure 3 - Configuring an HSI PHYSIO Sensor.....	3
Figure 4 - Configuring a Profusa Lumee Sensor	4
Figure 5 - Configuring an IOS MASES Sensor	5
Figure 6 - Configuring a Smart Eye Sensor.....	7
Figure 7 – Assigning Signal Processing.....	8
Figure 8 – Assigning a Threshold Calculator	9
Figure 9 – Assigning Threshold Ranges to a Signal	9
Figure 10 – Assigning a Core Temperature Estimator	10
Figure 11 – Configuring the Starting Temperature for Core Temperature Estimation	11
Figure 12 – Graphical Display of Signals	11
Figure 13 – Dashboard Configuration Save and Recall.....	12
Figure 14 – Loading a Dashboard Configuration	13

1.0 OVERVIEW

The design of the COG Pack Dashboard allows a user to set up, configure, and manage multiple Sensors by providing a standard look and feel for each Sensor. Generally, when manufacturers develop hardware sensors, they provide a custom application designed to only work with their particular device. From an integration standpoint, this provides a high barrier of entry for researchers and users when an experiment requires multiple disparate sensors. Each independent application often provides unique functionality for actions such as saving files, managing sensor state connectivity, and observing the data-streams in real-time. One of the objectives of the COG Pack Dashboard is to simplify these steps and consolidate them into a single interface. Additionally, by leveraging a database back-end for data storage, COG Pack automatically manages data timestamps and organization. By not requiring a user to execute multiple mouse actions and keyboard inputs for file management, COG Pack streamlines this process and reduces the potential for 'fat-fingered' input errors.

The goal of this document is to provide a user or researcher with the instructions they need to work with the COG Pack Dashboard. This document contains descriptions of the processes and procedures required to interact with the supported Sensors and Signals.

2.0 STARTING THE COG PACK DASHBOARD

The COG Pack system requires two main software components for execution: the back-end web service and the front-end Dashboard UI. These are the necessary components required, regardless of the configuration. In the most commonly deployed configuration, a single icon on the desktop loads both components. (Figure 1.)



Figure 1 - COG Pack Desktop Shortcut

3.0 ASSIGNING SENSORS TO PARTICIPANTS

After loading the COG Pack processes, the initialization screen is presented. (Figure 2.) This screen presents the first step in the process required for set up and configuration, allowing a user to associate particular Sensors and Signals with a participant.




Figure 2 - COG Pack Initialization Screen

To begin the association, complete the following steps:

1. Add a participant by clicking the “Add Participant” button on the left.
 - a. By default, the system generates a Participant ID
2. Click the generated Participant ID in the list.

3. Select the desired type of Sensor to assign to this participant by selecting the name of the Sensor from the array of available Sensors found in the “Add Sensor(s)” section.
4. Enter the values and settings required for each [Sensor](#).

Proceed to the following sections for examples showing how to configure the individual [Sensors](#).

3.1 HSI PHYSIO Configuration

Figure 3 displays the card for the configuration of an HSI PHYSIO Sensor. To assign an HSI PHYSIO Sensor to a participant, first label it with a Sensor ID by entering a value into the Sensor ID text field. Then select all of the Signals to capture. These selected Signals are then available to the Signal processing pipeline for viewing, logging, and further analysis. For this [Sensor](#), COG Pack has provided the ability to toggle between live and simulated data. To configure this option, toggle the button labeled “Emulated” to tell the system if this device is live or simulated. Finally, the user has the option to attach additional Signal processing algorithm(s). (See also [Attaching Additional Signal Processing Steps](#))

Figure 3 - Configuring an HSI PHYSIO Sensor

3.2 Profusa Lumee Configuration

To add a Profusa Lumee Sensor to a participant, locate the 4 character alpha-numeric ID on the Sensor, and enter that value into the Sensor ID text field. Then select all of the signals to capture. Finally, attach any additional Signal processing algorithms. (See also [Attaching](#)

[Additional Signal Processing Steps](#)) Figure 4 demonstrates the configuration of a Profusa Lumee Sensor.

The screenshot shows the COG Pack software interface. At the top, there is a header bar with the COG Pack logo and name. Below the header, there is a sidebar on the left with a list of configuration options: 'Add Participant', 'Study Configuration', 'Admin Settings', and a list of sensor IDs (0000001, 0000002, 0000003). The main area is titled 'Add Sensor(s)' and contains several buttons for different sensor types: SmartEye, PHYSIO, ThresholdCalculator, CoreTempEstimator, MASES, and RegionCalculator. The 'Lumee' button is selected. Below the buttons, there is a section titled 'Lumee' with the following configuration options: 'Enter Sensor ID:' with a text box containing 'E858', 'Assign Signal Processing:' with a dropdown menu showing 'Signal Processing', and 'Select Desired Signals:' with a grid of checkboxes. The checkboxes are arranged in two rows: Temperature, Oxygen Concentration, Accelerometer, Motion Classification, Diagnostic Bits, Battery, Baseline, Intensity, Gain, and LED Brightness. All checkboxes are checked. At the bottom of the 'Lumee' section, there is an 'Initialize' button.

Figure 4 - Configuring a Profusa Lumee Sensor

3.3 Intelligent Optical Systems (IOS) MASES Configuration

Adding an IOS MASES Sensor follows the same process as described for other Sensors. The IOS MASES Sensor requires only a few configuration options. First, a Sensor ID to uniquely describe this instance, then the IP Address of the interface connected to the IOS MASES Sensor. Additionally, specify the port number that the IOS MASES device uses to send data. COG Pack also supports emulation for this Sensor and can be toggled with the “Emulated” button. Finally, select all of the Signals to capture visually by clicking their respective check boxes. Figure 5 displays the configuration of an IOS MASES Sensor.

 **COG Pack**

Add Participant

Study Configuration

Admin Settings

0000001

0000002

Add Sensor(s)

SmartEye
PHYSIO
ThresholdCalculator
CoreTempEstimator
MASES
RegionCalculator

Lumee

Mases

Enter Sensor ID:

Enter Sensor Address:

Assign Signal Processing:

☐ Emulated

Select Desired Signals:

☒ Temperature
☒ Oxygen Concentration
☒ Carbon Dioxide Concentration
☒ Relative Humidity
☒ Pressure

Figure 5 - Configuring an IOS MASES Sensor

3.4 Smart Eye Pro Configuration

Currently, Smart Eye Pro is the most complicated Sensor to configure in COG Pack. To start, fill out all of the listed fields below with the settings formatted as described in Table 1.

Table 1. Smart Eye Settings, Descriptions, and Data Types

Setting Name	Description	Data Type
Sensor ID	Unique Identifier for this Smart Eye Instance	String
IP Address	IP Address of the COG Pack system connected to Smart Eye	String
Machine Name	Name of the COG Pack system	String
Port Number	Port Number used by Smart Eye to connect to COG Pack	Integer
World Model File	Full path to *.sew file used by Smart Eye	String
Minimum Fixation Dwell Time	For AdvancedFixations calculations, define the minimum number of milliseconds required for a Fixation	Double
Vergence Angle	For AdvancedFixations calculations, define the Visual Angle, in degrees, describing the maximum dispersion of points possible for a Fixation	Double
Monitor Size in Pixels	A comma-separated list of 4 Doubles that define the pixel space of a single monitor. DEPRECATED – replaced by values in the World Model File	DEPRECATED
Monitor Size in Meters	A comma-separated list of 4 Doubles that define the physical dimensions of a single monitor. DEPRECATED – replaced by values in the World Model File	DEPRECATED
Sample Rate	The frequency at which Smart Eye is configured. Most commonly 60 or 120 Hz.	Integer
Playback File	Full path to a *.sme file recorded previously by COG Pack software to replay previously recorded binary data. (Leave empty for live data capture)	String

These are the settings required for the configuration of the interface between COG Pack and Smart Eye. For additional settings of the visual elements of COG Pack, see the report titled “COG PACK Eye/Head/Gaze/Heatmap 3D Display Documentation”. Figure 6 shows the configuration of a Smart Eye Sensor.


COG Pack

Add Participant

Study Configuration

Admin Settings

0000001

0000002

0000003

0000004

Add Sensor(s)

SmartEye

PHYSIO

ThresholdCalculator

CoreTempEstimator

MASES

RegionCalculator

Lumee

SmartEye

☐ Binary Logging

Binary Logging Directory
 No file chosen

Sensor ID

IP Address

Machine Name

Port Number

World Model File

Minimum Fixation Dwell Time

Vergence Angle

Monitor Size in Pixels

Monitor Size in Meters

Sample Rate

Playback File

Assign Signal Processing:

Signal Processing

☒ Select All

Signals

☒ System Performance	☒ Gaze Intersection	☒ Blink	☒ Perclose
☒ Head Position	☒ Head Orientation	☒ Eye Orientation	☒ Eye Position
☒ Fixation	☒ Gaze Direction Quality	☒ Head Position Quality	☒ Advanced Fixations
☒ Raw Gaze	☒ Filtered Gaze	☒ Raw Estimated Gaze	☒ Filtered Estimated Gaze
☒ Closest World Intersection	☒ All World Intersections	☒ Filtered Closest World Intersection	
☒ Filtered All World Intersections	☒ Estimated Closest World Intersection		
☒ Estimated All World Intersections	☒ Robust Closest World Intersection		
☒ Robust All World Intersections	☒ Head Closest World Intersection		
☒ Head All World Intersections	☒ Calibration Gaze Intersection	☒ Tagged Gaze Intersection	
☒ Left Closest World Intersection	☒ Left All World Intersections		
☒ Right Closest World Intersection	☒ Right All World Intersections		
☒ Filtered Left Closest World Intersection	☒ Filtered Left All World Intersections		
☒ Filtered Right Closest World Intersection	☒ Filtered Right All World Intersections		
☒ Estimated Left Closest World Intersection	☒ Estimated Left All World Intersections		
☒ Estimated Right Closest World Intersection	☒ Estimated Right All World Intersections		

Figure 6 - Configuring a Smart Eye Sensor

4.0 ATTACHING ADDITIONAL SIGNAL PROCESSING STEPS

Each of the above cards for a Sensor's configuration has an option to add additional signal processing algorithm(s). To view all of the available algorithms for the selected Sensor, select the drop-down directly under the title **Assign Signal Processing (Figure 7)**.

The screenshot displays the COG Pack software interface. At the top, there is a header bar with a gear icon and the text "COG Pack". Below this, on the left, is a sidebar with buttons for "Add Participant", "Study Configuration", "Admin Settings", and a text field containing "0000001". The main area is titled "Add Sensor(s)" and contains several sensor cards: "MASES", "SmartEye", "ThresholdCalculator", "RegionCalculator", "CoreTempEstimator", and "PHYSIO". The "Physio" card is selected and expanded, showing a sub-header "Physio". Below this, there is a section "Enter Sensor ID:" with a text field containing "Mario" and a toggle switch labeled "Emulated". To the right of this is a section "Assign Signal Processing:" with a dropdown menu. The dropdown menu is open, showing three options: "Signal Processing", "CoreTempEstimator", and "ThresholdCalculator". Below the dropdown is a section "Select Desired Signals:" with a grid of checkboxes for various signals: Heart Rate, Battery, IBI, Respiration Rate, SpO2, Movement, Chest Temperature, Underarm Temperature, and Skin Temperature. At the bottom of the "Physio" card is an "Initialize" button.

Figure 7 – Assigning Signal Processing

4.1 Attaching a Threshold Calculator

To attach "Threshold Calculator" to a [Sensor](#) you must perform the following actions:

1. Select the Signal to Threshold for the current [Sensor](#).
2. Select the "Threshold Calculator" option for Signal processing as shown in **Figure 8**.
3. Fill in the fields with the desired range labels and values.
4. Click the button labeled "Attach".


COG Pack

Add Participant

Study Configuration

Admin Settings

0000001

0000002

0000003

0000004

0000005

Add Sensor(s)

SmartEye
PHYSIO
ThresholdCalculator
CoreTempEstimator
MASES
RegionCalculator

Lumee

Physio

Enter Sensor ID: ☐ Emulated

Assign Signal Processing:

ThresholdCalculator X

CoreTempEstimator

Select Desired Signals:

☐ Heart Rate

☐ Battery

☐ IBI

☐ Respiration Rate

☐ SpO2

☐ Movement

☐ Chest Temperature

☐ Underarm Temperature

☐ Skin Temperature

Thresholding

Range Labels	Range Values	
New Label	-	Attach

Initialize

Figure 8 – Assigning a Threshold Calculator

For example, to set the following values for the Signal “Respiration Rate,” enter the values in the same order as shown in **Figure 9**.

Low: < 60

Medium: >= 60 and < 120

High: >= 120

Thresholding

Range Labels	Range Values	
Low	60	
Medium	60 120	
High	120	
New Label	-	Attach

Figure 9 – Assigning Threshold Ranges to a Signal

4.2 Attaching a “Core Temperature Estimator”

The process for attaching a “Core Temperature Estimator” follows a similar method:

1. Select the **Signal** (Heart Rate) used to derive the Core Temperature Estimate (Figure 10.)
2. Select the “Core Temperature Estimator” option for **Signal** processing.
3. Fill in the value for “Starting Core Temperature” as seen in Figure 11.
4. Click the button labeled “Attach”.

The screenshot displays the COG Pack software interface. At the top, there is a header with the COG Pack logo and name. Below this, on the left, is a sidebar with buttons for 'Add Participant', 'Study Configuration', 'Admin Settings', and a participant ID '0000001'. The main area is titled 'Add Sensor(s)' and contains several sensor type buttons: 'PHYSIO', 'SmartEye', 'ThresholdCalculator', 'CoreTempEstimator', 'RegionCalculator', and 'Lumee'. The 'PHYSIO' button is selected, leading to a 'Physio' configuration panel. This panel includes an 'Enter Sensor ID:' field with the value 'Mario', an 'Emulated' toggle switch, and an 'Assign Signal Processing:' dropdown menu set to 'CoreTempEstimator'. Below these are 'Select Desired Signals:' checkboxes for Heart Rate, Battery, IBI, Respiration Rate, SpO2, Movement, Chest Temperature, Underarm Temperature, and Skin Temperature. A 'Core Temperature Estimation' section contains a 'Starting Core Temperature' dropdown menu set to 'Core Temp', with Fahrenheit and Celsius unit buttons below it. An 'Attach' button is located to the right of this section. At the bottom of the main panel is an 'Initialize' button.

Figure 10 – Assigning a Core Temperature Estimator

Core Temperature Estimation

Starting Core Temperature

37.1

Attach

°F

°C

Figure 11 – Configuring the Starting Temperature for Core Temperature Estimation

5.0 GRAPHS AND VISUALIZATIONS

When **Sensor** and **Signal** configuration is complete, click the button at the bottom of the screen labeled **Initialize**. This action will start the communication process between the front-end UI and back-end **COG_Pack.Communcations.exe** processes. All of the settings and configurations will be sent to the backend, and the **Sensors** will begin connecting. The UI will transition to a screen that renders all of the graphs configured in the previous screen. **Figure 12** is an example of a single participant and the output of the configured **Sensors** and **Signals**. This example contains the output of various Threshold **Signals** derived from the **Signals** produced by a single HSI PHYSIO **Sensor** and 2 separate Profusa Lumee **Sensors**.




Figure 12 – Graphical Display of Signals


The different colored sections on these graphs and gauges represent the configured values for each [Signal's](#) respective "Threshold Calculator" range values. In the COG Pack Threshold [Signal](#) definition, these values can change over time, and the graphs and gauges react accordingly.

6.0 LOADING AND SAVING CONFIGURATIONS

Another very important feature of this application is its ability to save the current setup and suit of sensors that you have configured for a group of participants into a single .JSON file. This file captures the settings and setup that you have entered into the application and can represent configuration for a specific trial.

For example, if a researcher wishes to use COG Pack for a study and it always requires the same 5 participants with the same set of sensors then you can just set it up in the application one time and hit save. Then, you will be able to load that saved file into the application later and it will restore the setup without you having to go through and reenter all of that data every time.

Therefore, on the main screen of the COG Pack Dashboard is an option on the left labeled **Save/Load (Figure 13)**.



The screenshot displays the COG Pack Dashboard interface. At the top, there is a header bar with the COG Pack logo and the text "COG Pack". Below the header, the dashboard is divided into several sections. On the left, there is a sidebar with a tab labeled "Add Participant" and a section titled "Study Configuration" containing a "Save/Load" tab. The "Save/Load" tab is active, showing a list of four participant IDs: 0000001, 0000002, 0000003, and 0000004. To the right of the sidebar, there is a main content area with two sections: "Upload Sensor Configuration File" and "Save Current Sensor Configurations". The "Upload Sensor Configuration File" section has a "Choose File" button and the text "No file chosen". The "Save Current Sensor Configurations" section has a "Save Filename:" input field and a "Save" button. At the bottom of the dashboard, there is an "Initialize" button.

Figure 13 – Dashboard Configuration Save and Recall

Loading a Configuration

The process for loading an existing configuration for the Dashboard is the following:

1. Navigate to the Save/Load tab on the left
2. Click the button labeled "Choose File"

3. Navigate to the .JSON file that you've previously saved
4. Click the open button as seen in **Figure 14** below.
5. Click the Initialize button.

After following the above steps and clicking the initialize button at the bottom of the screen you will see the screen displaying all of the visualizations and graphs associated with the sensors that were saved to the save file.




Figure 14 – Loading a Dashboard Configuration

[Saving a new Configuration](#)

To save your existing dashboard configuration, do the following:

1. Configure the sensors for your participants
2. Navigate to the Save/Load page in the application
3. Enter a name for the configuration file that will be generated
4. Click the save button

7.0 CONCLUSION

The COG Pack dashboard is a centralized place for multiple sensor configurations and visualizations. Our unique approach of using “Reflection Based” user interface generation allows this application to be very versatile when it comes to adding even more disparate sensors to its existing library of integrated sensors. The choice that we made to use Electron along with the React.js web framework allows this application to also be able to run on many different operating systems. We are constantly working on improving the design and user experience of the Dashboard to meet the needs and suggestions of our user base.