



Fragment Analysis Process for the Joint Trauma Analysis and Prevention of Injury in Combat Program

Karen Pizzolato

ARL-TR-5586

July 2011

NOTICES

Disclaimers

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

Citation of manufacturer's or trade names does not constitute an official endorsement or approval of the use thereof.

Destroy this report when it is no longer needed. Do not return it to the originator.

Army Research Laboratory

Aberdeen Proving Ground, MD 21005-5068

ARL-TR-5586**July 2011**

Fragment Analysis Process for the Joint Trauma Analysis and Prevention of Injury in Combat Program

Karen Pizzolato
Survivability/Lethality Analysis Directorate, ARL

REPORT DOCUMENTATION PAGE				Form Approved OMB No. 0704-0188	
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number. PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.					
1. REPORT DATE (DD-MM-YYYY) July 2011		2. REPORT TYPE Final		3. DATES COVERED (From - To) 2010	
4. TITLE AND SUBTITLE Fragment Analysis Process for the Joint Trauma Analysis and Prevention of Injury in Combat Program				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S) Karen Pizzolato				5d. PROJECT NUMBER AH80	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) U.S. Army Research Laboratory ATTN: RDRL-SLB-W Aberdeen Proving Ground, MD 21005-5068				8. PERFORMING ORGANIZATION REPORT NUMBER ARL-TR-5586	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution is unlimited.					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT Fragments from the threat and surrounding environment often become embedded in Service Members during an event in-theater. Analyzing these fragments is crucial to the characterization and identification of threats. Fragment analysis is one part of the Joint Trauma Analysis and Prevention of Injury in Combat (JTAPIC) Program. JTAPIC is a partnership program that was established at the U.S. Army Medical Research and Materiel Command. JTAPIC fragments are removed during autopsy by the Armed Forces Medical Examiner Systems from Service Members who were killed in action. These fragments are analyzed by the U.S. Army Research Laboratory (ARL). The analysis procedure includes the sterilization, scanning, and documentation of the physical properties of the fragments. ARL determines the elemental composition of fragments using qualitative and quantitative analysis procedures. The analysis results are used in event recreations and modeling and simulation. The results also assist JTAPIC partners with test designs and in understanding personnel vulnerabilities to threats.					
15. SUBJECT TERMS fragment, JTAPIC					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT UU	18. NUMBER OF PAGES 14	19a. NAME OF RESPONSIBLE PERSON Karen Pizzolato
a. REPORT UNCLASSIFIED	b. ABSTRACT UNCLASSIFIED	c. THIS PAGE UNCLASSIFIED			19b. TELEPHONE NUMBER (Include area code) (410) 278-4102

Standard Form 298 (Rev. 8/98)
Prescribed by ANSI Std. Z39.18

Contents

List of Figures	iv
1. Background	1
2. Procedure	1
3. Results	3
4. Conclusions	6
Distribution List	7

List of Figures

Figure 1. An example of the physical properties and description of a fragment recorded in the JTAPIC Fragment and Material Database.	3
Figure 2. A photograph taken of the fragment using an L-shaped forensic evidence ruler.....	4
Figure 3. A 3-D model of the fragment.	4
Figure 4. An SEM micrograph provides a high-resolution image of the surface of the fragment.	5
Figure 5. An EDS spectrum provides the elemental composition of the fragment.	5

1. Background

Based on DOD 6025.21E,¹ the Joint Trauma Analysis and Prevention of Injury in Combat (JTAPIC) Program was established at the U.S. Army Medical Research and Materiel Command. JTAPIC is a partnership program that links the operational, medical, and material communities with the purpose of collecting, integrating, and analyzing operational and injury data. Fragment analysis is one project that was established as a part of the JTAPIC Program.

Fragments embedded in U.S. Service Members during an event in-theater are removed during autopsies by the Armed Forces Medical Examiner Systems (AFMES). These fragments are transferred to the Warfighter Survivability Branch (WSB) of the U.S. Army Research Laboratory (ARL) where they are analyzed to determine their elemental composition and origin. The results of fragment analyses are used to characterize and identify threats to assist in determining enemy techniques, tactics, and procedures (TTPs). The results are used by JTAPIC partners to support event recreations, modeling and simulation, assist with armor and personal protective equipment test designs, and assist in understanding personnel vulnerabilities to threats. All data from the JTAPIC Fragment Analysis Project is stored in the JTAPIC Fragment and Material Database. This data includes fragment photographs, physical properties, three-dimensional (3-D) scans, and elemental analysis results. The fragment processing and analysis procedures used by ARL are described herein, and examples of the analysis results are included.

2. Procedure

JTAPIC fragments are removed by a medical examiner at the AFMES during autopsy. All fragments are considered forensic evidence and are therefore accompanied by U.S. Department of the Army Form 4173. A forensic investigator from AFMES transfers all fragments to a WSB evidence custodian for transport back to ARL. The following ARL organizations are involved in the fragment analysis process: the Safety and Health Physics Branch of the Associate Director for Laboratory Operations (ADLO) Directorate and the Multifunctional Materials Branch of ARL's Weapons and Materials Research Directorate (WMRD).

¹ U.S. Department of Defense (DOD). Directive 6025.21E; *Medical Research for Prevention, Mitigation, and Treatment of Blast Injuries*; 5 July 2006.

The first step in the analysis process is to ensure that all fragments are safe for handling by ARL personnel. Prior to direct handling, all fragments are scanned for radiation. The radiation scanning is conducted by ADLO personnel through the use of a hand-held geiger counter. Once the fragments have been determined to be free of radiation, they can be logged into the ARL Fragment Evidence Logbook, which is used to track the date and time of fragment transfers. Next, photographs are taken to record the condition of the fragments upon receipt by ARL. At this time, metallic and plastic fragments are separated into two groups for sterilization.

WMRD performs sterilization to remove any potential biological hazard. Two different sterilization procedures are used. Metallic fragments are sterilized using an autoclave method. A VSR model AS12 autoclave is used to heat the metallic fragments at 100 °C for 100 min at 17 pounds per square inch (psi) of pressure. Fragments that appear to be plastic are soaked in a 10% formalin solution for 24 h. This method ensures that the plastic will not melt while providing the same level of sterilization as the autoclave method.

After sterilization is completed, the fragments can be safely handled by the fragment analysts. Each fragment is photographed, and its physical properties are recorded. The mass, length, width, and depth of each fragment are measured, and a light microscope is used to record a basic physical description of the fragments. A second set of photographs is taken of all fragments using an L-shaped forensic evidence ruler to preserve the image of the fragment evidence.

Three-dimensional scans are taken of selected fragments. These scans, which are performed using a NextEngine HD 3-D scanner, provide fragment analysts with the exact 3-D geometry and shape factor of the fragment of interest. Two scans are recorded for each fragment. The second scan is taken in an orthogonal position in relation to the first scan. These two scans are merged using the Geomagic computer program to create a 360° 3-D representation of the fragment.

Once all physical properties have been recorded, qualitative elemental analysis is performed using Scanning Electron Microscopy-Energy Dispersive X-ray Spectroscopy (SEM-EDS). Particular cases of interest also undergo quantitative analyzing using Inductively Coupled Plasma-Atomic Emission Spectroscopy (ICP-AES). A Hitachi S-4700 Field Emission Scanning Electron Microscope (SEM) with EDS capabilities is used for qualitative elemental analysis. Several representative areas are identified on each fragment, and a low magnification (35×) micrograph is taken of these areas. Qualitative elemental analysis is performed using EDS. This technique employs an electron beam to penetrate the surface of the fragment. The electron beam excites the surface atoms of the fragment, producing X-rays that possess energies characteristic of the elements found in the sample. A detector is used to convert the X-ray energies into voltage signals, which are sent to a pulse-processor for detection. The results are displayed as an elemental spectrum. Based on the results of the qualitative analysis and any input from the JTAPIC partners, selected fragments will undergo quantitative analysis using ICP-AES. ICP-AES is used to determine the exact elemental percentages of the trace concentrations of metals in a fragment sample, which can be used for metal alloy identification.

After the completion of all analyses, the results are recorded in the JTAPIC Fragment and Material Database. This is a database that was created by WSB to provide fragment analysts with an easily accessible storage tool. The database contains all fragment photographs, physical properties data, 3-D scans, and elemental analysis results. It also contains chain-of-custody information and medical data received from the AFMES that relates to each case.

3. Results

The following is an example of the results from one fragment analysis. The fragment was determined to be 1.33 g, with dimensions of $13.93 \times 14.42 \times 3.01$ mm. Use of a light microscope revealed a silver colored fragment with brown matter and multiple tan and blue fibers adhering to the surface of the fragment. The physical properties, recovery location, and description of the fragment were recorded in the JTAPIC Fragment and Material Database (figure 1).

Mass: 1.33 g	Description: Silver colored with reddish-brown encrusted matter and ridges. One end is more jagged and irregular whereas other is more square with rounded edges. Many tan and reddish fibers on convex side along with a small blue fiber near a small hole (located near area where fragment folded onto itself). Near jagged end on reverse side are several tan fibers and fabric with spherical adhesions along concave side.
Dimension: $13.93 \times 14.42 \times 3.01$ mm	
Density: 2.20 g/mL	
Shape: Irregular	
Recovery Location: Lateral thigh	
Predominant Materials: Copper and Iron	

Figure 1. An example of the physical properties and description of a fragment recorded in the JTAPIC Fragment and Material Database.

A photograph of the fragment was taken with an L-shaped forensic evidence ruler (figure 2). A 3-D scan (figure 3) preserved the fragment's shape and provided the necessary geometry needed for model and simulation. Scanning was also used to determine the fragment's density, which was 2.20 grams per milliliter (g/mL).

A high-resolution (35 $\times$) image of the surface of the fragment was captured using SEM-EDS (figure 4). SEM-EDS was also used to determine the elemental composition of this fragment, which was found to be predominately copper and iron (figure 5).

The results of this fragment analysis were used to make determinations about the fragment identification and origin to assist analysts in threat characterization and determining enemy TTPs.



Figure 2. A photograph taken of the fragment using an L-shaped forensic evidence ruler.




Figure 3. A 3-D model of the fragment.




Figure 4. An SEM micrograph provides a high-resolution image of the surface of the fragment.




Figure 5. An EDS spectrum provides the elemental composition of the fragment.

4. Conclusions

ARL has created a procedure for processing and analyzing fragments for the JTAPIC Program. The results from ARL's fragment analyses are provided to all JTAPIC partners. The partners combine fragment results with operational and medical information to support event recreations, assist with test designs, help determine enemy TTPS, and assist in understanding personnel vulnerabilities to threats.

NO. OF
COPIES ORGANIZATION

1 (PDF only)	DEFENSE TECHNICAL INFORMATION CTR DTIC OCA 8725 JOHN J KINGMAN RD STE 0944 FORT BELVOIR VA 22060-6218
1	DIRECTOR US ARMY RESEARCH LAB IMNE ALC HRR 2800 POWDER MILL RD ADELPHI MD 20783-1197
1	DIRECTOR US ARMY RESEARCH LAB RDRL CIO LL 2800 POWDER MILL RD ADELPHI MD 20783-1197
1	DIRECTOR US ARMY RESEARCH LAB RDRL CIO MT 2800 POWDER MILL RD ADELPHI MD 20783-1197
1	DIRECTOR US ARMY RESEARCH LAB RDRL D 2800 POWDER MILL RD ADELPHI MD 20783-1197

NO. OF
COPIES ORGANIZATION

1 ASST SECY ARMY
(CD ACQSTN LOGISTICS & TECH
only) SAAL ZS RM 2E533
103 ARMY PENTAGON
WASHINGTON DC 20310-0103

1 DIRECTOR FORCE DEV
(CD DAPR FDZ
only) RM 3A522
460 ARMY PENTAGON
WASHINGTON DC 20310-0460

1 US ARMY TRADOC ANL CTR
ATRC W
A KEINTZ
WSMR NM 88002-5502

1 USARL
RDRL SLE
R FLORES
WSMR NM 88002-5513

ABERDEEN PROVING GROUND

1 US ARMY DEV TEST COM
TEDT TMT
314 LONGS CORNER RD
APG MD 21005-5055

1 US ARMY EVALUATION CTR
TEAE SV
P A THOMPSON
4120 SUSQUEHANNA AVE
APG MD 21005-3013

12 DIR USARL
RDRL SL
J BEILFUSS
J FEENEY
J FRANZ
M STARKS
P TANENBAUM
RDRL SLB A
R DIBELKA
G MANNIX
RDRL SLB D
R GROTE
RDRL SLB E
M PERRY
RDRL SLB G
P MERGLER

NO. OF
COPIES ORGANIZATION

RDRL SLB S
S SNEAD
RDRL SLB W
L ROACH