



Efficacy of an Extroral Suction System in Reducing Aerosol Contamination

Benjamin Wang, DDS
Capt, USAF, DC
AADOCR/CADR Annual Meeting
23 Mar 2022

Disclaimer



- The views expressed are those of the authors and do not reflect the official views or policy of the Uniformed Services University, Department of Defense, or its Components. The authors do not have any financial interest in the companies whose materials are discussed in this abstract. The views of pH Dental are not necessarily the official views of, or endorsed by, the U.S. Government, the Department of Defense, or the Department of the Air Force. No Federal endorsement of pH Dental is intended.

Co-Contributors



- Lori Henrichs, MS, MT (ASCP)
- Kelli Arricale, RDH
- Wen Lien, DMD, MS, MS
Col, USAF, DC
- Kraig Vandewalle, DMD, MS
Col (Ret.), USAF, DC

Outline



- Background
- Objective
- Materials and Methods
- Results
- Discussion

Background



- Ongoing worldwide COVID-19 pandemic
 - Predominantly respiratory mode of infection
- Routine production of contaminated aerosols during dental procedures
 - Transmit HIV, TB, SARS, measles, influenza, Legionnaire's disease, hepatitis (Cottone et al. 1996, Gruninger et al. 1992)
- Significant concerns amongst the public regarding the safety of dental procedures

Background



- Infection control practices in dentistry
 - Standard Precautions
 - ✦ Immunization, hand washing, personal protective equipment, sterilization/disinfection
 - Rubber dam
 - ✦ Eliminates 98% of microbial contamination (Harrel 1998)
 - Evacuators or low-volume suction
 - ✦ Neither reduced aerosols during ultrasonic scaling effectively (Holloman et al. 2015)



Background



- Previous studies support that extraoral scavenging (EOS) devices and air purifiers reduce dental aerosol contamination
- New paradigm: high-volume extraoral suction (EOS) with high-efficiency particulate air (HEPA) filtration
 - Zhao et al. (2020) reported HEPA air purifier filters removed 83% of aerosols in-vitro, compared to 54% by fine filters
 - Shahdad et al. (2020) found that EOS device usage produced 33% and 76% reduction in mean intensity contamination for the operator and assistant
 - Graetz et al. (2021) further substantiated that an EOS device significantly reduced the number of generated particles during dental treatment
- Limited previous studies have tested the clinical effectiveness of newer EOS devices with HEPA filtration in reducing aerosol contamination

PAX2000



- PAX2000 Extraoral Dental Suction System
 - Unique, multi-stage carbon filtration system
 - Medical grade H13 HEPA filtration
 - UV-C bulb light system
 - Kills collected microorganisms >0.3 microns at 99.95% efficiency



Objective



- To evaluate the efficacy of a chairside, extraoral scavenging device (PAX 2000 Extraoral Dental Suction System) on the reduction of aerosol contamination during dental treatment



Null Hypotheses



- There will be no difference in microbial load, composition, and spatial distribution in the dental operatory with or without the use of the extra-oral suction device during ultrasonic cleaning of a dental patient




Materials and Methods




- Ultrasonic scaling was performed on 40 patients treatment-planned for initial or supportive periodontal therapy
 - 2 test groups of 20 patients each; randomized exposure of EOS device during ultrasonic cleaning
 - Exempt human study: Implied consent was obtained from subjects using an informational informed consent letter

Materials and Methods

- Sample collection involved both passive and active sampling
 - Passive sampling utilized petri dishes with blood agar/R2A agar
 - ✦ 3 time points
 - Before treatment (30 mins)
 - During treatment
 - After treatment (30 mins)
 - ✦ 3 locations
 - Patient chest
 - Instrument tray table
 - Countertop by foot of dental chair (control)





Materials and Methods



Materials and Methods

- Colony forming units were incubated for a week and then quantified per location and timepoint after a week-long incubation period



Data Analysis/Power Analysis



- Outcome data (microbial load, composition, and spatial distribution in the dental operator) were calculated all locations and timepoints
- Data were analyzed using Mann-Whitney U tests ($\alpha=0.05$)
- 87% power was achieved by sample size of 20 per test group assuming a large effect size for a 2-sided test

Results



- All collected dental unit water line samples contained 0 CFU/mL
- Lowest mean CFUs were found in samples collected before and after treatment and from samples placed at the control location
- Use of the EOS suction device reduced the number of CFUs during treatment at all locations
 - Statistically significant ($p = 0.018$) at the patient chest area where the highest microbial load was present during all time points

Time	Location	Colony Forming Units Median (IQR)		
		No Suction Device	Suction Device [#]	P value
Before	Delivery Unit	1 (3)	1 (2)	0.097
	Patient Chest	1 (2)	1 (2)	0.217
	Patient Foot	1 (2)	0 (1)	0.147
During [#]	Delivery Unit	4 (3)	3 (3.75)	0.122
	Patient Chest	6 (10.75)	4 (5)	0.018*
	Patient Foot	3 (3)	2 (3)	0.067
After	Delivery Unit	2 (2)	2 (3)	0.842
	Patient Chest	2 (2)	2 (3)	0.156
	Patient Foot	1 (1.75)	2 (3)	0.142

[#]Suction device was only active during treatment
^{*}Indicates significant difference

Conclusion



- The EOS suction system may serve to reduce aerosol contamination in the clinical dental setting, especially in proximity to the patient's chest region, where most aerosols were generated



Questions?

