

AD A 120835

12

**REPORT DRXTH-TE-CR-82179**

**SURFACE SAMPLING TECHNIQUES**

Bruce E. Goodwin  
James R. Aronson  
Robert P. O'Neil  
Margaret A. Randel  
Emmett M. Smith

**ARTHUR D. LITTLE, INC.  
CAMBRIDGE, MA 02140**

**SEPTEMBER, 1982**

**FINAL REPORT  
Volume II  
Certification Testing Data**

**Distribution Unlimited  
Cleared for Public Release**

*prepared for*

U.S. Army Toxic and Hazardous Materials Agency,  
Aberdeen Proving Ground, Maryland 21010

**DTIC  
ELECTE  
OCT 29 1982  
S H**

**Best Available Copy**

copy

#### DISCLAIMER

The views, opinions and/or findings contained in this report are those of the authors and should not be construed as an official Department of the Army position, policy or decision, unless so designated by other documents.

***REPORT DRXTH-TE-CR-82179***

**SURFACE SAMPLING TECHNIQUES**

Bruce E. Goodwin  
James R. Aronson  
Robert P. O'Neil  
Margaret A. Randel  
Emmett M. Smith

***ARTHUR D. LITTLE, INC.  
CAMBRIDGE, MA 02140***

***SEPTEMBER, 1982***

**FINAL REPORT  
Volume II  
Certification Testing Data**

**Distribution Unlimited  
Cleared for Public Release**

***prepared for***

**U.S. Army Toxic and Hazardous Materials Agency,  
Aberdeen Proving Ground, Maryland 21010**



# TABLE OF CONTENTS :

|                                                                                                          | Page |
|----------------------------------------------------------------------------------------------------------|------|
| Abbreviations                                                                                            | ii   |
| List of Tables                                                                                           | iii  |
| List of Figures                                                                                          | xi   |
| > Semiquantitative Certification Data;                                                                   | 5    |
| Quantitative Method for the Determination of<br>DNP, RDX, TNB, 2,4-DNT, TNT, Tetryl and DPA on Surfaces; |      |
| Quantitative Method for the Determination of<br>2,6-DNT and NG on Surfaces;                              |      |
| Quantitative Method for the Determination of<br>PETN on Surfaces; o.v. 6                                 |      |
| Quantitative Concentration Data.                                                                         |      |



|                    |                                     |
|--------------------|-------------------------------------|
| Accession For      |                                     |
| NTIS GRA&I         | <input checked="" type="checkbox"/> |
| DTIC TAB           | <input type="checkbox"/>            |
| Unannounced        | <input type="checkbox"/>            |
| Justification      |                                     |
| By                 |                                     |
| Distribution/      |                                     |
| Availability Codes |                                     |
| Avail and/or       |                                     |
| Dist Special       |                                     |
| A                  |                                     |

# ABBREVIATIONS

| <u>Abbreviation</u> | <u>Chemical Name</u>                |
|---------------------|-------------------------------------|
| DNP                 | 2,4-Dinitrophenol                   |
| RDX                 | Cyclotrimethylenetrinitramine       |
| TNB                 | 1,3,5-Trinitrobenzne                |
| 2,4-DNT             | 2,4-Dinitrotoluene                  |
| 2,4,6-TNT           | 2,4,6-Trinitrotoluene               |
| Tetryl              | 2,4,6-Trinitrophenylmethylnitramine |
| DPA                 | Diphenylamine                       |
| 2,6-DNT             | 2,6-Dinitrotoluene                  |
| NG                  | Nitroglycerine                      |
| PETN                | Pentaerythritetetranitrate          |

LIST OF TABLES

| <u>TABLE NO.</u>                     |                                                                            | <u>PAGE</u> |
|--------------------------------------|----------------------------------------------------------------------------|-------------|
| Semiquantitative Certification Data: |                                                                            |             |
| I-1                                  | Semiquantitative Certification Testing<br>Analytical Method Conditions     | 1           |
| I-2                                  | Semiquantitative Certification Testing<br>Statistical Data Summary         | 2           |
| I-3                                  | DNP - Found Concentrations                                                 | 3           |
| I-4                                  | DNP - Analysis of Target-Found Concentration Points                        | 4           |
| I-5                                  | RDX - Found Concentrations                                                 | 6           |
| I-6                                  | RDX - Analysis of Target-Found Concentration Points                        | 7           |
| I-7                                  | TNB - Found Concentrations-HPLC Method                                     | 9           |
| I-8                                  | TNB - Analysis of Target-Found Concentration Points-<br>HPLC Method        | 10          |
| I-9                                  | TNB - Found Concentrations-GC/ECD Method                                   | 12          |
| I-10                                 | TNB - Analysis of Target-Found Concentration Points-<br>GC/ECD Method      | 13          |
| I-11                                 | 2,4 - DNT - Found Concentrations-HPLC Method                               | 15          |
| I-12                                 | 2,4 - DNT - Analysis of Target-Found Concentration<br>Points-HPLC Method   | 16          |
| I-13                                 | 2,4 - DNT - Found Concentrations-GC/ECD Method                             | 18          |
| I-14                                 | 2,4 - DNT - Analysis of Target-Found Concentration<br>Points-GC/ECD Method | 19          |
| I-15                                 | TNT - Found Concentrations-HPLC Method                                     | 21          |
| I-16                                 | TNT - Analysis of Target-Found Concentration Points-<br>HPLC Method        | 22          |

LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                        | <u>PAGE</u> |
|------------------|------------------------------------------------------------------------|-------------|
| I-17             | TNT - Found Concentrations-GC/ECD Method                               | 24          |
| I-18             | TNT - Analysis of Target Found Concentration<br>Points-GC/ECD Method   | 25          |
| I-19             | Tetryl - Found Concentrations                                          | 27          |
| I-20             | Tetryl - Analysis of Target-Found Concentration<br>Points              | 28          |
| I-21             | DPA - Found Concentrations                                             | 30          |
| I-22             | DPA - Analysis of Target-Found Concentration Points                    | 31          |
| I-23             | 2,6-DNT-Found Concentrations-HPLC Method                               | 33          |
| I-24             | 2,6-DNT-Analysis of Target-Found Concentration<br>Points-HPLC Method   | 34          |
| I-25             | 2,6-DNT-Found Concentrations-GC/ECD Method                             | 36          |
| I-26             | 2,6-DNT-Analysis of Target-Found Concentration<br>Points-GC/ECD Method | 37          |
| I-27             | NG-Found Concentrations                                                | 39          |
| I-28             | NG-Analysis of Target-Found Concentration Points                       | 40          |
| I-29             | PETN-Found Concentrations                                              | 42          |
| I-30             | PETN-Analysis of Target-Found Concentration Points                     | 43          |

LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                             | <u>PAGE</u> |
|------------------|-----------------------------------------------------------------------------|-------------|
|                  | Quantitative Certification Data                                             |             |
| II-1             | Quantitative Certification Testing<br>Statistical Data Summary              | 45          |
| II-2             | DNP on Metal-Target vs. Found Concentrations                                | 58          |
| II-3             | DNP on Metal-Analysis of Target-Found Concentration<br>Points               | 59          |
| II-4             | DNP on Metal-Inaccuracy and Imprecision Data                                | 61          |
| II-5             | DNP on Concrete-Target vs. Found Concentrations                             | 64          |
| II-6             | DNP on Concrete-Analysis of Target-Found Concentration<br>Points            | 65          |
| II-7             | DNP on Concrete-Inaccuracy and Imprecision Data                             | 67          |
| II-8             | DNP on Concrete (3 days) - Target vs. Found<br>Concentrations               | 70          |
| II-9             | DNP on Concrete (3 days) - Analysis of Target-Found<br>Concentration Points | 71          |
| II-10            | DNP on Concrete (3 days) - Inaccuracy and Imprecision<br>Data               | 73          |
| II-11            | DNP on Brick - Target vs. Found Concentrations                              | 76          |
| II-12            | DNP on Brick - Analysis of Target-Found Concentration<br>Points             | 77          |
| II-13            | DNP on Brick - Inaccuracy and Imprecision Data                              | 79          |
| II-14            | DNP on Transite - Target vs. Found Concentrations                           | 82          |
| II-15            | DNP on Transite - Analysis of Target-Found<br>Concentration Points          | 83          |
| II-16            | DNP on Transite - Inaccuracy and Imprecision Data                           | 85          |
| II-17            | DNP on Transite (3 days) - Target vs. Found Concen-<br>trations             | 88          |
| II-18            | DNP on Transite (3 days) - Analysis of Target-Found<br>Concentration Points | 89          |



LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                          | <u>PAGE</u> |
|------------------|--------------------------------------------------------------------------|-------------|
| II-19            | DNP on Transite (3 days) - Inaccuracy and Imprecision Data               | 91          |
| II-20            | RDX on Metal - Target vs. Found Concentrations                           | 94          |
| II-21            | RDX on Metal - Analysis of Target-Found Concentration Points             | 95          |
| II-22            | RDX on Metal - Inaccuracy and Imprecision Data                           | 97          |
| II-23            | RDX on Concrete - Target vs. Found Concentrations                        | 100         |
| II-24            | RDX on Concrete - Analysis of Target-Found Concentration Points          | 101         |
| II-25            | RDX on Concrete - Inaccuracy and Imprecision Data                        | 103         |
| II-26            | RDX on Brick - Target vs. Found Concentrations                           | 106         |
| II-27            | RDX on Brick - Analysis of Target-Found Concentration Points             | 107         |
| II-28            | RDX on Brick - Inaccuracy and Imprecision Data                           | 109         |
| II-29            | RDX on Transite - Target vs. Found Concentrations                        | 112         |
| II-30            | RDX on Transite - Analysis of Target-Found Concentration Points          | 113         |
| II-31            | RDX on Transite - Inaccuracy and Imprecision Data                        | 115         |
| II-32            | RDX on Transite (3 days) - Target vs. Found Concentrations               | 118         |
| II-33            | RDX on Transite (3 days) - Analysis of Target-Found Concentration Points | 119         |
| II-34            | RDX on Transite (3 days) - Inaccuracy and Imprecision Data               | 121         |
| II-35            | TNB on Metal - Target vs. Found Concentrations                           | 124         |
| II-36            | TNB on Metal - Analysis of Target-Found Concentration Points             | 125         |

LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                          | <u>PAGE</u> |
|------------------|--------------------------------------------------------------------------|-------------|
| II-37            | TNB on Metal - Inaccuracy and Imprecision Data                           | 127         |
| II-38            | TNB on Concrete - Target vs. Found Concentrations                        | 130         |
| II-39            | TNB on Concrete - Analysis of Target-Found Concentration Points          | 131         |
| II-40            | TNB on Concrete - Inaccuracy and Imprecision Data                        | 133         |
| II-41            | TNB on Brick - Target vs. Found Concentrations                           | 136         |
| II-42            | TNB on Brick - Analysis of Target-Found Concentration Points             | 137         |
| II-43            | TNB on Brick - Inaccuracy and Imprecision Data                           | 139         |
| II-44            | TNB on Transite - Target vs. Found Concentrations                        | 142         |
| II-45            | TNB on Transite - Analysis of Target-Found Concentration Points          | 143         |
| II-46            | TNB on Transite - Inaccuracy and Imprecision Data                        | 145         |
| II-47            | TNB on Transite (3 days) - Target vs. Found Concentrations               | 148         |
| II-48            | TNB on Transite (3 days) - Analysis of Target-Found Concentration Points | 149         |
| II-49            | TNB on Transite (3 days) - Inaccuracy and Imprecision Data               | 151         |
| II-50            | 2,4-DNT on Metal - Target vs. Found Concentrations                       | 154         |
| II-51            | 2,4-DNT on Metal - Analysis of Target-Found Concentration Points         | 155         |
| II-52            | 2,4-DNT on Metal - Inaccuracy and Imprecision Data                       | 157         |
| II-53            | 2,4-DNT on Concrete - Target vs. Found Concentrations                    | 160         |
| II-54            | 2,4-DNT on Concrete - Analysis of Target Found Concentration Points      | 161         |
| II-55            | 2,4-DNT on Concrete-Inaccuracy and Imprecision Data                      | 163         |

LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                            | <u>PAGE</u> |
|------------------|----------------------------------------------------------------------------|-------------|
| II-56            | 2,4-DNT on Brick-Target vs. Found Concentrations                           | 166         |
| II-57            | 2,4-DNT on Brick-Analysis of Target-Found Concentration Points             | 167         |
| II-58            | 2,4-DNT on Brick-Inaccuracy and Imprecision Data                           | 169         |
| II-59            | 2,4-DNT on Transite-Target vs. Found Concentrations                        | 172         |
| II-60            | 2,4-DNT on Transite-Analysis of Target-Found Concentration Points          | 173         |
| II-61            | 2,4-DNT on Transite-Inaccuracy and Imprecision Data                        | 175         |
| II-62            | 2,4-DNT on Transite (3 days)-Target vs. Found Concentrations               | 178         |
| II-63            | 2,4-DNT on Transite (3 days)-Analysis of Target-Found Concentration Points | 179         |
| II-64            | 2,4-DNT on Transite (3 days)-Inaccuracy and Imprecision Data               | 181         |
| II-65            | TNT on Metal - Target vs. Found Concentrations                             | 184         |
| II-66            | TNT on Metal - Analysis of Target-Found Concentration Points               | 185         |
| II-67            | TNT on Metal - Inaccuracy and Imprecision Data                             | 187         |
| II-68            | TNT on Concrete - Target vs. Found Concentrations                          | 190         |
| II-69            | TNT on Concrete - Analysis of Target-Found Concentration Points            | 191         |
| II-70            | TNT on Concrete - Inaccuracy and Imprecision Data                          | 193         |
| II-71            | TNT on Brick - Target vs. Found Concentrations                             | 196         |
| II-72            | TNT on Brick - Analysis of Target-Found Concentration Points               | 197         |
| II-73            | TNT on Brick - Inaccuracy and Imprecision Data                             | 199         |
| II-74            | TNT on Transite - Target vs. Found Concentrations                          | 202         |

LIST OF TABLES

(Continued)

| <u>TABLE NO.</u> |                                                                          | <u>PAGE</u> |
|------------------|--------------------------------------------------------------------------|-------------|
| II-75            | TNT on Transite - Analysis of Target-Found Concentration Points          | 203         |
| II-76            | TNT on Transite - Inaccuracy and Imprecision Data                        | 205         |
| II-77            | TNT on Transite (3 days) - Target vs. Found Concentrations               | 208         |
| II-78            | TNT on Transite (3 days) - Analysis of Target-Found Concentration Points | 209         |
| II-79            | TNT on Transite (3 days) - Inaccuracy and Imprecision Data               | 211         |
| II-80            | Tetryl on Metal - Target vs. Found Concentrations                        | 214         |
| II-81            | Tetryl on Metal - Analysis of Target-Found Concentration Points          | 215         |
| II-82            | Tetryl on Metal - Inaccuracy and Imprecision Data                        | 217         |
| II-83            | Tetryl on Concrete - Target vs. Found Concentrations                     | 220         |
| II-84            | Tetryl on Concrete - Analysis of Target-Found Concentration Points       | 221         |
| II-85            | Tetryl on Concrete - Inaccuracy and Imprecision Data                     | 223         |
| II-86            | Tetryl on Brick - Target vs. Found Concentrations                        | 226         |
| II-87            | Tetryl on Brick - Analysis of Target-Found Concentration Points          | 227         |
| II-88            | Tetryl on Brick - Inaccuracy and Imprecision Data                        | 229         |
| II-89            | Tetryl on Transite - Target vs. Found Concentrations                     | 232         |
| II-90            | Tetryl on Transite - Analysis of Target-Found Concentration Points       | 233         |
| II-91            | Tetryl on Transite - Inaccuracy and Imprecision Data                     | 235         |
| II-92            | DPA on Metal - Target vs. Found Concentrations                           | 238         |
| II-93            | DPA on Metal - Analysis of Target-Found Concentration Points             | 239         |

LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                          | <u>PAGE</u> |
|------------------|--------------------------------------------------------------------------|-------------|
| II-94            | DPA on Metal - Inaccuracy and Imprecision Data                           | 241         |
| II-95            | DPA on Concrete - Target vs. Found Concentrations                        | 244         |
| II-96            | DPA on Concrete - Analysis of Target-Found Concentration Points          | 245         |
| II-97            | DPA on Concrete - Inaccuracy and Imprecision Data                        | 247         |
| II-98            | DPA on Concrete (3 days) - Target vs. Found Concentrations               | 250         |
| II-99            | DPA on Concrete (3 days) - Analysis of Target-Found Concentration Points | 251         |
| II-100           | DPA on Concrete (3 days) - Inaccuracy and Imprecision Data               | 253         |
| II-101           | DPA on Brick - Target vs. Found Concentrations                           | 256         |
| II-102           | DPA on Brick - Analysis of Target-Found Concentration Points             | 257         |
| II-103           | DPA on Brick - Inaccuracy and Imprecision Data                           | 259         |
| II-104           | DPA on Transite - Target vs. Found Concentrations                        | 262         |
| II-105           | DPA on Transite - Analysis of Target-Found Concentration Points          | 263         |
| II-106           | DPA on Transite - Inaccuracy and Imprecision Data                        | 265         |
| II-107           | DPA on Transite (3 days) - Target vs. Found Concentrations               | 268         |
| II-108           | DPA on Transite (3 days) - Analysis of Target-Found Concentration Points | 269         |
| II-109           | DPA on Transite (3 days) - Inaccuracy and Imprecision Data               | 271         |
| II-110           | 2,6-DNT on Metal - Target vs. Found Concentrations                       | 274         |
| II-111           | 2,6-DNT on Metal - Analysis of Target-Found Concentration Points         | 275         |

LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                     | <u>PAGE</u> |
|------------------|---------------------------------------------------------------------|-------------|
| II-112           | 2,6-DNT on Metal - Inaccuracy and Imprecision Data                  | 277         |
| II-113           | 2,6-DNT on Concrete - Target vs. Found Concentrations               | 280         |
| II-114           | 2,6-DNT on Concrete - Analysis of Target-Found Concentration Points | 281         |
| II-115           | 2,6-DNT on Concrete - Inaccuracy and Imprecision Data               | 283         |
| II-116           | 2,6-DNT on Brick - Target vs. Found Concentrations                  | 286         |
| II-117           | 2,6-DNT on Brick - Analysis of Target-Found Concentration Points    | 287         |
| II-118           | 2,6-DNT on Brick - Inaccuracy and Imprecision Data                  | 289         |
| II-119           | 2,6-DNT on Transite - Target vs. Found Concentrations               | 292         |
| II-120           | 2,6-DNT on Transite - Analysis of Target-Found Concentration Points | 293         |
| II-121           | 2,6-DNT on Transite - Inaccuracy and Imprecision Data               | 295         |
| II-122           | NG on Metal - Target vs. Found Concentrations                       | 298         |
| II-123           | NG on Metal - Analysis of Target-Found Concentration Points         | 299         |
| II-124           | NG on Metal - Inaccuracy and Imprecision Data                       | 301         |
| II-125           | NG on Concrete Target vs. Found Concentrations                      | 304         |
| II-126           | NG on Concrete - Analysis of Target-Found Concentration Points      | 305         |
| II-127           | NG on Concrete - Inaccuracy and Imprecision Data                    | 307         |
| II-128           | NG on Brick - Target vs. Found Concentrations                       | 310         |
| II-129           | NG on Brick - Analysis of Target-Found Concentration Points         | 311         |

LIST OF TABLES  
(Continued)

| <u>TABLE NO.</u> |                                                                  | <u>PAGE</u> |
|------------------|------------------------------------------------------------------|-------------|
| II-130           | NG on Brick - Inaccuracy and Imprecision Data                    | 313         |
| II-131           | NG on Transite - Target vs. Found Concentrations                 | 316         |
| II-132           | NG on Transite - Analysis of Target-Found Concentration Points   | 317         |
| II-133           | NG on Transite - Inaccuracy and Imprecision Data                 | 319         |
| II-134           | PETN on Metal - Target vs. Found Concentrations                  | 322         |
| II-135           | PETN on Metal - Analysis of Target-Found Concentration Points    | 324         |
| II-136           | PETN on Concrete - Target vs. Found Concentrations               | 326         |
| II-137           | PETN on Concrete - Analysis of Target-Found Concentration Points | 328         |
| II-138           | PETN on Brick - Target vs. Found Concentrations                  | 33-         |
| II-139           | PETN on Brick - Analysis of Target-Found Concentration Points    | 332         |
| II-140           | PETN on Transite - Target vs. Found Concentrations               | 334         |
| II-141           | PETN on Transite - Analysis of Target-Found Concentration Points | 336         |

## LIST OF FIGURES

| <u>FIGURE NO.</u>                    |                                                                      | <u>PAGE</u> |
|--------------------------------------|----------------------------------------------------------------------|-------------|
| Semiquantitative Certification Data: |                                                                      |             |
| I-1                                  | DNP - Graph of Target-Found Concentration Points                     | 5           |
| I-2                                  | RDX - Graph of Target-Found Concentration Points                     | 8           |
| I-3                                  | TNB - Graph of Target-Found Concentration Points-<br>HPLC Method     | 11          |
| I-4                                  | TNB - Graph of Target-Found Concentration Points-<br>GC/ECD Method   | 14          |
| I-5                                  | 2,4-DNT-Graph of Target-Found Concentration<br>Points-HPLC Method    | 17          |
| I-6                                  | 2,4-DNT-Graph of Target-Found Concentration<br>Points-GC/ECD Method  | 20          |
| I-7                                  | TNT - Graph of Target-Found Concentration<br>Points-HPLC Method      | 25          |
| I-8                                  | TNT - Graph of Target-Found Concentration<br>Points-GC/ECD Method    | 26          |
| I-9                                  | Tetryl - Graph of Target-Found Concentration Points                  | 29          |
| I-10                                 | DPA - Graph of Target-Found Concentration Points                     | 32          |
| I-11                                 | 2,6-DNT-Graph of Target-Found Concentration Points-<br>HPLC Method   | 35          |
| I-12                                 | 2,6-DNT-Graph of Target-Found Concentration Points-<br>GC/ECD Method | 38          |
| I-13                                 | NG - Graph of Target-Found Concentration Points                      | 41          |
| I-14                                 | PETN - Graph of Target-Found Concentration Points                    | 44          |



## LIST OF FIGURES

| <u>FIGURE NO.</u> |                                                                       | <u>PAGE</u> |
|-------------------|-----------------------------------------------------------------------|-------------|
|                   | Quantitative Certification Data                                       |             |
| II-1              | DNP on Metal - Graph of Target-Found Concentration Points             | 60          |
| II-2              | DNP on Metal - Graph of Inaccuracy                                    | 62          |
| II-3              | DNP on Metal - Graph of Imprecision                                   | 63          |
| II-4              | DNP on Concrete - Graph of Target-Found Concentration Points          | 66          |
| II-5              | DNP on Concrete - Graph of Inaccuracy                                 | 68          |
| II-6              | DNP on Concrete - Graph of Imprecision                                | 69          |
| II-7              | DNP on Concrete (3 days) - Graph of Target-Found Concentration Points | 72          |
| II-8              | DNP on Concrete (3 days) - Graph of Inaccuracy                        | 74          |
| II-9              | DNP on Concrete (3 days) - Graph of Imprecision                       | 75          |
| II-10             | DNP on Brick - Graph of Target-Found Concentration Points             | 78          |
| II-11             | DNP on Brick - Graph of Inaccuracy                                    | 80          |
| II-12             | DNP on Brick - Graph of Imprecision                                   | 81          |
| II-13             | DNP on Transite - Graph of Target-Found Concentration Points          | 84          |
| II-14             | DNP on Transite - Graph of Inaccuracy                                 | 86          |
| II-15             | DNP on Transite - Graph of Imprecision                                | 87          |
| II-16             | DNP on Transite (3 days) - Graph of Target-Found Concentration Points | 90          |
| II-17             | DNP on Transite (3 days) - Graph of Inaccuracy                        | 92          |
| II-18             | DNP on Transite (3 days) - Graph of Imprecision                       | 93          |
| II-19             | RDX on Metal - Graph of Target-Found Concentration Points             | 96          |

LIST OF FIGURES  
(Continued)

| <u>FIGURE NO.</u> |                                                                       | <u>PAGE</u> |
|-------------------|-----------------------------------------------------------------------|-------------|
| II-20             | RDX on Metal - Graph of Inaccuracy                                    | 98          |
| II-21             | RDX on Metal - Graph of Imprecision                                   | 99          |
| II-22             | RDX on Concrete - Graph of Target-Found Concentration Points          | 102         |
| II-23             | RDX on Concrete - Graph of Inaccuracy                                 | 104         |
| II-24             | RDX on Concrete - Graph of Imprecision                                | 105         |
| II-25             | RDX on Brick - Graph of Target-Found Concentration Points             | 108         |
| II-26             | RDX on Brick - Graph of Inaccuracy                                    | 110         |
| II-27             | RDX on Brick - Graph of Imprecision                                   | 111         |
| II-28             | RDX on Transite - Graph of Target-Found Concentration Points          | 114         |
| II-29             | RDX on Transite - Graph of Inaccuracy                                 | 116         |
| II-30             | RDX on Transite - Graph of Imprecision                                | 117         |
| II-31             | RDX on Transite (3 days) - Graph of Target-Found Concentration Points | 120         |
| II-32             | RDX on Transite (3 days) - Graph of Inaccuracy                        | 122         |
| II-33             | RDX on Transite (3 days) - Graph of Imprecision                       | 123         |
| II-34             | TNB on Metal - Graph of Target-Found Concentration Points             | 126         |
| II-35             | TNB on Metal - Graph of Inaccuracy                                    | 128         |
| II-36             | TNB on Metal - Graph of Imprecision                                   | 129         |
| II-37             | TNB on Concrete - Graph of Target-Found Concentration Points          | 132         |
| II-38             | TNB on Concrete - Graph of Inaccuracy                                 | 134         |
| II-39             | TNB on Concrete - Graph of Imprecision                                | 135         |

LIST OF FIGURES  
(Continued)

| <u>FIGURE NO.</u> |                                                                       | <u>PAGE</u> |
|-------------------|-----------------------------------------------------------------------|-------------|
| II-40             | TNB on Brick - Graph of Target-Found Concentration Points             | 138         |
| II-41             | TNB on Brick - Graph of Inaccuracy                                    | 140         |
| II-42             | TNB on Brick - Graph of Imprecision                                   | 141         |
| II-43             | TNB on Transite - Graph of Target-Found Concentration Points          | 144         |
| II-44             | TNB on Transite - Graph of Inaccuracy                                 | 146         |
| II-45             | TNB on Transite - Graph of Imprecision                                | 147         |
| II-46             | TNB on Transite (3 days) - Graph of Target-Found Concentration Points | 150         |
| II-47             | TNB on Transite (3 days) - Graph of Inaccuracy                        | 152         |
| II-48             | TNB on Transite (3 days) - Graph of Imprecision                       | 153         |
| II-49             | 2,4-DNT on Metal - Graph of Target-Found Concentration Points         | 156         |
| II-50             | 2,4-DNT on Metal - Graph of Inaccuracy                                | 158         |
| II-51             | 2,4-DNT on Metal - Graph of Imprecision                               | 159         |
| II-52             | 2,4-DNT on Concrete - Graph of Target-Concentration Points            | 162         |
| II-53             | 2,4-DNT on Concrete - Graph of Inaccuracy                             | 164         |
| II-54             | 2,4-DNT on Concrete - Graph of Imprecision                            | 165         |
| II-55             | 2,4-DNT on Brick - Graph of Target-Found Concentration Points         | 168         |
| II-56             | 2,4-DNT on Brick - Graph of Inaccuracy                                | 170         |
| II-57             | 2,4-DNT on Brick - Graph of Imprecision                               | 171         |
| II-58             | 2,4-DNT on Transite - Graph of Target-Found Concentration Points      | 174         |

LIST OF FIGURES  
(Continued)

| <u>FIGURE NO.</u> |                                                                               | <u>PAGE</u> |
|-------------------|-------------------------------------------------------------------------------|-------------|
| II-59             | 2,4-DNT on Transite - Graph of Inaccuracy                                     | 176         |
| II-60             | 2,4-DNT on Transite - Graph of Imprecision                                    | 177         |
| II-61             | 2,4-DNT on Transite (3 days) - Graph of Target-<br>Found Concentration Points | 180         |
| II-62             | 2,4-DNT on Transite (3 days) - Graph of Inaccuracy                            | 182         |
| II-63             | 2,4-DNT on Transite (3 days) - Graph of Imprecision                           | 183         |
| II-64             | TNT on Metal - Graph of Target-Found Concentration<br>Points                  | 186         |
| II-65             | TNT on Metal - Graph of Inaccuracy                                            | 188         |
| II-66             | TNT on Metal - Graph of Imprecision                                           | 189         |
| II-67             | TNT on Concrete - Graph of Target-Found Concen-<br>tration Points             | 192         |
| II-68             | TNT on Concrete - Graph of Inaccuracy                                         | 194         |
| II-69             | TNT on Concrete - Graph of Imprecision                                        | 195         |
| II-70             | TNT on Brick - Graph of Target-Found Concen-<br>tration Points                | 198         |
| II-71             | TNT on Brick - Graph of Inaccuracy                                            | 200         |
| II-72             | TNT on Brick - Graph of Imprecision                                           | 201         |
| II-73             | TNT on Transite - Graph of Target-Found Concen-<br>tration Points             | 204         |
| II-74             | TNT on Transite - Graph of Inaccuracy                                         | 206         |
| II-75             | TNT on Transite - Graph of Imprecision                                        | 207         |
| II-76             | TNT on Transite (3 days) - Graph of Target-Found<br>Concentration Points      | 210         |
| II-77             | TNT on Transite (3 days) - Graph of Inaccuracy                                | 212         |
| II-78             | TNT on Transite (3 days) - Graph of Imprecision                               | 213         |

LIST OF FIGURES  
(Continued)

| <u>FIGURE NO.</u> |                                                                       | <u>PAGE</u> |
|-------------------|-----------------------------------------------------------------------|-------------|
| II-79             | Tetryl on Metal - Graph of Target-Found Concentration Points          | 216         |
| II-80             | Tetryl on Metal - Graph of Inaccuracy                                 | 218         |
| II-81             | Tetryl on Metal - Graph of Imprecision                                | 219         |
| II-82             | Tetryl on Concrete - Graph of Target-Found Concentration Points       | 222         |
| II-83             | Tetryl on Concrete - Graph of Inaccuracy                              | 224         |
| II-84             | Tetryl on Concrete - Graph of Imprecision                             | 225         |
| II-85             | Tetryl on Brick - Graph of Target-Found Concentration Points          | 228         |
| II-86             | Tetryl on Brick - Graph of Inaccuracy                                 | 230         |
| II-87             | Tetryl on Brick - Graph of Imprecision                                | 231         |
| II-88             | Tetryl on Transite - Graph of Target-Found Concentration Points       | 234         |
| II-89             | Tetryl on Transite - Graph of Inaccuracy                              | 236         |
| II-90             | Tetryl on Transite - Graph of Imprecision                             | 237         |
| II-91             | DPA on Metal - Graph of Target-Found Concentration Points             | 240         |
| II-92             | DPA on Metal - Graph of Inaccuracy                                    | 242         |
| II-93             | DPA on Metal - Graph of Imprecision                                   | 243         |
| II-94             | DPA on Concrete - Graph of Target-Found Concentration Points          | 246         |
| II-95             | DPA on Concrete - Graph of Inaccuracy                                 | 248         |
| II-96             | DPA on Concrete - Graph of Imprecision                                | 249         |
| II-97             | DPA on Concrete (3 days) - Graph of Target-Found Concentration Points | 252         |

LIST OF FIGURES  
(Continued)

| <u>FIGURE NO.</u> |                                                                       | <u>PAGE</u> |
|-------------------|-----------------------------------------------------------------------|-------------|
| II-98             | DPA on Concrete (3 days) - Graph of Inaccuracy                        | 254         |
| II-99             | DPA on Concrete (3 days) - Graph of Imprecision                       | 255         |
| II-100            | DPA on Brick - Graph of Target-Found Concentration Points             | 258         |
| II-101            | DPA on Brick - Graph of Inaccuracy                                    | 260         |
| II-102            | DPA on Brick - Graph of Imprecision                                   | 261         |
| II-103            | DPA on Transite - Graph of Target-Found Concentration Points          | 264         |
| II-104            | DPA on Transite - Graph of Inaccuracy                                 | 266         |
| II-105            | DPA on Transite - Graph of Imprecision                                | 267         |
| II-106            | DPA on Transite (3 days) - Graph of Target-Found Concentration Points | 270         |
| II-107            | DPA on Transite (3 days) - Graph of Inaccuracy                        | 272         |
| II-108            | DPA on Transite (3 days) - Graph of Imprecision                       | 273         |
| II-109            | 2,6-DNT on Metal - Graph of Target-Found Concentration Points         | 276         |
| II-110            | 2,6-DNT on Metal - Graph of Inaccuracy                                | 278         |
| II-111            | 2,6-DNT on Metal - Graph of Imprecision                               | 279         |
| II-112            | 2,6-DNT on Concrete - Graph of Target-Found Concentration Points      | 282         |
| II-113            | 2,6-DNT on Concrete - Graph of Inaccuracy                             | 284         |
| II-114            | 2,6-DNT on Concrete - Graph of Imprecision                            | 285         |
| II-115            | 2,6-DNT on Brick - Graph of Target-Found Concentration Points         | 288         |
| II-116            | 2,6-DNT on Brick - Graph of Inaccuracy                                | 290         |
| II-117            | 2,6-DNT on Brick - Graph of Imprecision                               | 291         |

LIST OF FIGURES  
(Continued)

| <u>FIGURE NO.</u> |                                                                  | <u>PAGE</u> |
|-------------------|------------------------------------------------------------------|-------------|
| II-118            | 2,6-DNT on Transite - Graph of Target-Found Concentration Points | 294         |
| II-119            | 2,6-DNT on Transite - Graph of Inaccuracy                        | 296         |
| II-120            | 2,6-DNT on Transite - Graph of Imprecision                       | 297         |
| II-121            | NG on Metal - Graph of Target-Found Concentration Points         | 300         |
| II-122            | NG on Metal - Graph of Inaccuracy                                | 302         |
| II-123            | NG on Metal - Graph of Imprecision                               | 303         |
| II-124            | NG on Concrete - Graph of Target-Found Concentration Points      | 306         |
| II-125            | NG on Concrete - Graph of Inaccuracy                             | 308         |
| II-126            | NG on Concrete - Graph of Imprecision                            | 309         |
| II-127            | NG on Brick - Graph of Target-Found Concentration Points         | 312         |
| II-128            | NG on Brick - Graph of Inaccuracy                                | 314         |
| II-129            | NG on Brick - Graph of Imprecision                               | 315         |
| II-130            | NG on Transite - Graph of Target-Found Concentration Points      | 318         |
| II-131            | NG on Transite - Graph of Inaccuracy                             | 320         |
| II-132            | NG on Transite - Graph of Imprecision                            | 321         |
| II-133            | PETN on Metal - Graph of Target-Found Concentration Points       | 325         |
| II-134            | PETN on Concrete - Graph of Target-Found Concentration Points    | 329         |
| II-135            | PETN on Brick - Graph of Target-Found Concentration Points       | 333         |
| II-136            | PETN on Transite - Graph of Target-Found Concentration Points    | 337         |

Table I-1. Semiquantitative Certification Testing!  
Analytical Method Conditions

| Analyte   | Instrument | Column                                                           | Precolumn                          | Program                                                                                              | Temp of Solvent                                                                                      | Retention Time | Solvent System/<br>Carrier Gas                                                           | Flow Rate  |
|-----------|------------|------------------------------------------------------------------|------------------------------------|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|----------------|------------------------------------------------------------------------------------------|------------|
| PEIN      | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Isocratic                                                                                            | Isocratic                                                                                            | 6.2 min        | 65/35 CH <sub>3</sub> CH <sub>3</sub> /H <sub>2</sub> O                                  | 1 mL/min   |
| 2,6-DNT   | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Isocratic                                                                                            | Isocratic                                                                                            | 24.6 min       | 35/65 CH <sub>3</sub> CN/<br>.005M t-butyl                                               | 1 mL/min   |
|           | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Isocratic                                                                                            | Isocratic                                                                                            | 30.5 min       | ammonium hydroxide                                                                       | 1 mL/min   |
|           | GC         | 3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column | None                               | 100C for 6 min<br>15C/min to 165C                                                                    | 100C for 6 min<br>15C/min to 165C                                                                    | 19.90 min      | 5% Methane/argon                                                                         | 30 mL/min  |
| 1,3,5-TNB | GC         | 3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column | None                               | hold 8 min<br>15C/min to 200 C                                                                       | hold 8 min<br>15C/min to 200 C                                                                       | 16.20 min      | 5% Methane/argon                                                                         | 30 mL/min  |
|           | GC         | 3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column | None                               | hold 6 min                                                                                           | hold 6 min                                                                                           | 30.95 min      | 5% Methane/argon                                                                         | 30 mL/min  |
|           | GC         | 3% OV-225 on 100/120 Gas ChromQ 1/8" x 2 mm ID x 6' glass column | None                               | "                                                                                                    | "                                                                                                    | 28.14 min      | 5% Methane/argon                                                                         | 30 mL/min  |
| DPA       | HPIC       | U Bondpack C18 4mm x 30 cm                                       | None                               | Isocratic                                                                                            | Isocratic                                                                                            | 2.30 sec       | 90/10 methanol/H <sub>2</sub> O                                                          | 1 mL/min   |
| Tetryl    | HPIC       | U Bondpack C18 4mm x 30 cm                                       | None                               | Isocratic                                                                                            | Isocratic                                                                                            | 416 sec        | 60/40 methanol/H <sub>2</sub> O                                                          | 1 mL/min   |
| DNT       | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | 12.5 min       | 0.08M acetic acid<br>adjusted to pH 3.1<br>with ammonium<br>hydroxide/CH <sub>3</sub> CN | 1.0 mL/min |
| RDX       | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | 13.5 min       | 0.08M acetic acid<br>adjusted to pH 3.1<br>with ammonium<br>hydroxide/CH <sub>3</sub> CN | 1.0 mL/min |
| TNB       | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | 17.2 min       | 0.08M acetic acid<br>adjusted to pH 3.1<br>with ammonium<br>hydroxide/CH <sub>3</sub> CN | 1.0 mL/min |
| 2,4-DNT   | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | 24.3 min       | 0.08M acetic acid<br>adjusted to pH 3.1<br>with ammonium<br>hydroxide/CH <sub>3</sub> CN | 1.0 mL/min |
| TNT       | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | 25.6 min       | 0.08M acetic acid<br>adjusted to pH 3.1<br>with ammonium<br>hydroxide/CH <sub>3</sub> CN | 1.0 mL/min |
| Tetryl    | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | 26.7 min       | 0.08M acetic acid<br>adjusted to pH 3.1<br>with ammonium<br>hydroxide/CH <sub>3</sub> CN | 1.0 mL/min |
| DPA       | HPIC       | Spherisorb ODS, 5µ, 250 x 4.6 mm                                 | Pellicular LC-18, 40µ, 50 x 4.6 mm | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | Initial: 30% CH <sub>3</sub> CN<br>Final: 50% CH <sub>3</sub> CN<br>Time: 35 min<br>Gradient: linear | 39.2 min       | 0.08M acetic acid<br>adjusted to pH 3.1<br>with ammonium<br>hydroxide/CH <sub>3</sub> CN | 1.0 mL/min |



Table I-2. Semiquantitative Certification Testing  
Statistical Data Summary

| Analyte   | Detector     | Attn.                 | Chart<br>Speed | Inj.<br>Volume | USAT/AMA<br>Method | Det.<br>Limit | Corr.<br>Coeff. | Slope | Int.   | MR<br>Reference |
|-----------|--------------|-----------------------|----------------|----------------|--------------------|---------------|-----------------|-------|--------|-----------------|
| PEIN      | UV at 230 nm | 0.01 AUFS             | 0.1 in/min     | 70 uL          | none               | 1.77 ug/mL    | 0.996           | 1.046 | 0.450  | 9               |
| 2,6-DNT   | UV at 230 nm | 0.01 AUFS             | 0.1 in/min     | 70 uL          | MRI Method         | 0.26 ug/mL    | 0.998           | 1.061 | 9.888  | 8               |
| NI        | UV at 230 nm | 0.01 AUFS             | 0.1 in/min     | 70 uL          | MRI Method         | 4.54 ug/mL    | 0.999           | 0.989 | 2119.5 | 8               |
| 2,4-DNT   | ECD @300C    | 1x10 <sup>-11</sup> g | 0.5 cm/min     | 1 uL           | 246 TNT-MA-02      | 0.11 ug/mL    | 0.967           | 0.743 | 0.034  | 5               |
| 2,6-DNT   | ECD @300C    | 1x10 <sup>-11</sup> g | 0.5 cm/min     | 1 uL           | 246 TNT-MA-02      | 0.09 ug/mL    | 0.989           | 1.069 | -0.012 | 5               |
| 1,3,5-TNB | ECD @300C    | 1x10 <sup>-11</sup> g | 0.5 cm/min     | 1 uL           | 246 TNT-MA-02      | 0.12 ug/mL    | 0.999           | 1.058 | 0.007  | 5               |
| 2,4,6-TNT | ECD @300C    | 1x10 <sup>-11</sup> g | 0.5 cm/min     | 1 uL           | 246 TNT-MA-07      | 0.12 ug/mL    | 0.979           | 1.054 | -0.009 | 5               |
| DPA       | UV at 254 nm | 0.1 AUFS              | 0.5 cm/min     | 200 uL         | DPA-MA-01          | 12.25 ug/mL   | 0.997           | 0.796 | -0.634 | 3               |
| Tetryl    | UV at 254 nm | 0.1 AUFS              | 0.5 cm/min     | 200 uL         | TETRYL-MA-02       | 0.08 ug/mL    | 0.995           | 0.862 | -0.004 | 5               |
| DMP       | UV at 254 nm | 0.01 AUFS             | 0.1 in/min     | 100 uL         | MRI Method         | 30 ng/mL      | 0.998           | 0.971 | 5.397  | 10              |
| BDX       | UV at 254 nm | 0.01 AUFS             | 0.1 in/min     | 100 uL         | MIR Method         | 70 ng/mL      | 0.998           | 1.009 | 9.698  | 10              |
| TNB       | UV at 254 nm | 0.01 AUFS             | 0.1 in/min     | 100 uL         | MRI Method         | 25 ng/mL      | 0.999           | 0.994 | 0.674  | 10              |
| 2,4-DNT   | UV at 254 nm | 0.01 AUFS             | 0.1 in/min     | 100 uL         | MRI Method         | 24 ng/mL      | 0.999           | 0.990 | 1.585  | 10              |
| TNT       | UV at 254 nm | 0.01 AUFS             | 0.1 in/min     | 100 uL         | MRI Method         | 24 ng/mL      | 0.999           | 1.012 | 0.269  | 10              |
| Tetryl    | UV at 254 nm | 0.01 AUFS             | 0.1 in/min     | 100 uL         | MRI Method         | 25 ng/mL      | 0.999           | 1.050 | -5.410 | 10              |
| DPA       | UV at 254 nm | 0.01 AUFS             | 0.1 in/min     | 100 uL         | MRI Method         | 51 ng/mL      | 0.999           | 0.993 | 7.745  | 10              |

Table I-3. DNP - Found Concentrations

| ARMY CERTIFICATION |                      |
|--------------------|----------------------|
| DNP                | FOUND CONC           |
| (1)                | 26 24 (2) 26 28      |
| (3)                | 26 29 (4) 52 52      |
| (5)                | 52 52 (6) 52 53      |
| (7)                | 104 104 (8) 104 126  |
| (9)                | 104 111 (10) 260 254 |
| (11)               | 260 261 (12) 260 264 |
| (13)               | 520 495 (14) 520 533 |
| (15)               | 520 498              |

Table I-4. DNP - Analysis of Target-Found Concentration Points

ARMY CERTIFICATION  
DNP  
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 192.4 SD 189.231830605

FOUND CONC  
MEAN= 192.266666667 SD= 184.077804461

NO. RUNS : TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 5.39670981661  
SLOPE= 0.971257571986  
USE FOR ACCURACY  
R= 0.998451979989  
MEAN SQR DEV OF POINTS FROM REGRESSION= 112.090610730  
ST ERROR EST= 10.6249993288  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 13  
TWO TAIL P LEVEL IS .1  
t= 1.77093170942  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U  
REPLICATES= 2  
y(c)= 19.6057408586  
x(d)= 29.1259198183

ARMY CERTIFICATION  
DNP  
FOUND CONC

533.00

$y(c)$

$x(d)$

0.00

VERTICAL AXIS TIC INTERVAL = 53.3  
TARGET CONC

520

416

312

208

104

Figure I-1. DNP - Graph of Target-Found Concentration Points

Table I-5. RDX - Found Concentrations

| ARMY CERTIFICATION |                          |
|--------------------|--------------------------|
| RDX                | FOUND CONC               |
| (1)                | 50 44 (2) 50 44          |
| (3)                | 50 60 (4) 100 95         |
| (5)                | 100 107 (6) 100 110      |
| (7)                | 200 235 (8) 200 193      |
| (9)                | 200 271 (10) 500 500     |
| (11)               | 500 516 (12) 500 521     |
| (13)               | 1000 1000 (14) 1000 1066 |
| (15)               | 1000 984                 |

Table I-6. RDX - Analysis of Target-Found Concentration Points

ARMY CERTIFICATION  
 RDX  
 ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
 MEAN= 378 SD= 363.907366349

FOUND CONC  
 MEAN= 383.066666667 SD= 368.11963635

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15  
 MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 9.69795037756  
 SLOPE= 1.00910463862  
 USE FOR ACCURACY  
 R= 0.997557791953  
 MEAN SQR DEV OF POINTS FROM REGRESSION= 711.942079499  
 ST ERROR EST= 26.6822427749  
 USE FOR PRECISION  
 T FOR CONFIDENCE BAND  
 D.F.= 13  
 TWO TAIL P LEVEL IS .1  
 t= 1.77093170942  
 X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U  
 REPLICATES= 2  
 u(c)= 45.3006040026  
 x(d)= 70.3255003599

ARMY CERTIFICATION

RDX

FOUND CONC

1066.00




Figure I-2. RDX - Graph of Target-Found Concentration Points

Table I-7. TNB - Found Concentrations - HPLC Method

**ARMY CERTIFICATION**

**TNB**

**FOUND CONC**

**(1) 25 25 (2) 25 26**

**(3) 25 28 (4) 50 49**

**(5) 50 50 (6) 50 52**

**(7) 100 96 (8) 100 99**

**(9) 100 104 (10) 250 243**

**(11) 250 253 (12) 250 250**

**(13) 500 486 (14) 500 518**

**(15) 500 491**



Table I-8. TNE - Analysis of Target-Found  
Concentration Points - HPLC Method

ARMY CERTIFICATION  
TNE  
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 185 SD= 181.953683274

FOUND CONC  
MEAN= 184.666666667 SD= 181.09889432

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.674487594391  
SLOPE= 0.99455231931  
USE FOR ACCURACY  
R= 0.999246618194  
MEAN SQR DEV OF POINTS FROM REGRESSION= 53.1983838461  
ST ERROR EST= 7.2937166826  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 13  
TWO TAIL P LEVEL IS .1  
t= 1.77893178942  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U  
REPLICATES= 2  
y(c)= 10.4285310439  
x(d)= 19.5511771834

ARMY CERTIFICATION

TNB

FOUND CONC

518.00




Figure 1-3. TNB - Graph of Target-Found Concentration Points - HPLC Method

Table I-9. TNB-Found Concentrations - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

1,3,5 TNB

FOUND CONC

|      |      |      |      |      |      |
|------|------|------|------|------|------|
| (1)  | 0.12 | 0.09 | (2)  | 0.12 | 0.1  |
| (3)  | 0.12 | 0.09 | (4)  | 0.12 | 0.17 |
| (5)  | 0.23 | 0.22 | (6)  | 0.23 | 0.31 |
| (7)  | 0.23 | 0.25 | (8)  | 0.46 | 0.53 |
| (9)  | 0.46 | 0.55 | (10) | 0.46 | 0.48 |
| (11) | 0.46 | 0.49 | (12) | 1.16 | 1.23 |
| (13) | 1.16 | 1.27 | (14) | 1.16 | 1.26 |
| (15) | 2.31 | 2.42 | (16) | 2.31 | 2.44 |
| (17) | 2.31 | 2.46 | (18) | 3    | 0    |
| (19) | 0    | 0    | (20) | 0    | 0    |
| (21) | 0    | 0    |      |      |      |

Table I-10. TNB - Analysis of Target-Found Concentration Points - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

1.3.5 TNB

ANALYSIS OF 21 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.639047619048 SD= 0.791573779812

FOUND CONC

MEAN= 0.68380952381 SD= 0.838793634874

NB. RUNS 1 TOTAL X-Y ALL RUNS 21 NB. CONCENTR 21  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.60709275380941

SLOPE= 1.05894575030

R= 0.99332439997

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 3.884665285E-4  
COMPUTE T

D.F. = 19

ENTER 2 TAIL P LEVEL (USUALLY .1. EACH CONFID BAND IS

.05 SO TOTAL P= .1)

.1

t= 1.729129474

REPLICATES ON UNKNOWN SAMPLE 3

Y(CC)= 0.0420524309347

X(CC)= 0.0657723484472

ARMY CERTIFICATION GC/ECD METHOD  
 1,3,5 TNB  
 FOUND CONC




Figure I-4. TNB - Graph of Target-Found Concentration Points - GC/ECD Method

Table I-11. 2,4-DNT - Found Concentrations - HPLC Method

ARMY CERTIFICATION

24DNT

FOUND CONC

(1) 24 26 (2) 24 26

(3) 24 27 (4) 48 50

(5) 48 48 (6) 48 50

(7) 96 91 (8) 96 94

(9) 96 100 (10) 240 238

(11) 240 241 (12) 240 240

(13) 480 471 (14) 480 492

(15) 480 468

Table I-12. 2,4-DNT - Analysis of Target-Found  
Concentration Points - HPLC Method

ARMY CERTIFICATION  
24DNT  
ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 177.6 SD= 174.675535943

FOUND CONC  
MEAN= 177.466666667 SD= 173.068882516

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 1.58495145631  
SLOPE= 0.990324973031  
USE FOR ACCURACY  
R= 0.999518474424  
MEAN SQR DEV OF POINTS FROM REGRESSION= 31.037567837  
ST ERROR EST= 5.37293170935  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F. = 13  
TWO TAIL P LEVEL IS .1  
t= 1.77093170942  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U  
REPLICATES= 2  
y(c)= 9.03775232583  
x(d)= 15.0116760434

ARMY CERTIFICATION

24DNT  
FOUND CONC

492.00

0

y(c)

0.00

480

384

288

192

96

TARGET CONC

VERTICAL AXIS TIC INTERVAL = 49.2

Figure I-5. 2,4-DNT - Graph of Target-Found Concentration  
Points - HPLC Method



Table I-13. 2,4-DNT-Found Concentrations - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

2,4 DNT

FOUND CONC

|      |      |      |      |      |      |
|------|------|------|------|------|------|
| (1)  | 0.03 | 0.05 | (2)  | 0.03 | 0.04 |
| (3)  | 0.03 | 0.06 | (4)  | 0.06 | 0.06 |
| (5)  | 0.06 | 0.08 | (6)  | 0.06 | 0.08 |
| (7)  | 0.11 | 0.1  | (8)  | 0.11 | 0.11 |
| (9)  | 0.11 | 0.14 | (10) | 0.11 | 0.21 |
| (11) | 0.28 | 0.2  | (12) | 0.28 | 0.27 |
| (13) | 0.28 | 0.31 | (14) | 0.28 | 0.31 |
| (15) | 0.56 | 0.43 | (16) | 0.56 | 0.42 |
| (17) | 0.56 | 0.42 | (18) | 0.03 | 0.07 |
| (19) | 0    | 0    | (20) | 0    | 0    |
| (21) | 0    | 0.01 | (22) | 0    | 0    |

Table I-14. 2,4-DNT - Analysis of Target-Found Concentration Points - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD  
 2,4 DNT  
 ANALYSIS OF 22 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
 MEAN= 0.160909090909 SD= 0.188677328259

FOUND CONC  
 MEAN= 0.153181818182 SD= 0.144992162798

NO. RUNS 1 TOTAL X-Y ALL RUNS 22 NO. CONCENTR 22  
 MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.0336084831092  
 SLOPE= 0.7431112192  
 R= 0.367805515185  
 MEAN SQ. DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.0014326811139  
 COMPUTE 1

D.F. = 20  
 ENTER 2 TAIL P LEVEL (USUALLY .1, EACH CONFID BAND IS  
 .05 SO TOTAL P= .1)  
 .1  
 t= 1.72471716621  
 REPLICATES ON UNKNOWN SAMPLE 3  
 y(c)= 0.0755821824279  
 x(d)= 0.110789286561

Figure 8 is a graph showing the relationship between Target Concentration (X-axis) and TIC Interval (Y-axis). The X-axis ranges from 0.00 to 0.56, and the Y-axis ranges from 0.13 to 0.43. Two parallel lines are plotted, with data points marked by '0'. A vertical line is drawn at X=0.112, intersecting the lower line at Y=0.13 and the upper line at Y=0.224. The vertical distance between these two points is labeled 'VERTICAL AXIS TIC INTERVAL = 0.043'.

# itle

Table I-15. TNT - Found Concentrations - HPLC Method

| ARMY CERTIFICATION |                      |
|--------------------|----------------------|
| TNT                |                      |
| FOUND CONC         |                      |
| (1)                | 24 25 (2) 24 26      |
| (3)                | 24 26 (4) 48 48      |
| (5)                | 48 49 (6) 48 50      |
| (7)                | 96 89 (8) 96 98      |
| (9)                | 96 102 (10) 244 242  |
| (11)               | 244 248 (12) 244 250 |
| (13)               | 488 490 (14) 488 507 |
| (15)               | 488 487              |

Table I-16. TNT - Analysis of Target-Found Concentration  
Points - HPLC Method

ARMY CERTIFICATION

TNT

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 180 SD= 177.950233014

FOUND CONC  
MEAN= 182.46666667 SD= 180.195475869

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.269286848564  
SLOPE= 1.01220766566  
USE FOR ACCURACY  
R= 0.999595522042  
MEAN SQR DEV OF POINTS FROM REGRESSION= 28.2819574468  
ST ERROR EST= 5.31807836035

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLET C/U U

REPLICATES= 2

y(c)= 7.3789069858

x(d)= 14.0139688505




Figure I-7. TNT - Graph of Target-Found Concentration Points - HPLC Method

Table I-17. TNT - Found Concentrations - GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD:

2,4,6 TNT

FOUND CONC

(1) 0.04 0.02 (2) 0.04 0.03

(3) 0.04 0.1 (4) 0.04 0.00

(5) 0.00 0.12 (5) 0.00 0

(7) 0.00 0.06 (8) 0.00 0.15

(9) 0.15 0.11 (10) 0.15 0.14

(11) 0.15 0.00 (12) 0.15 0.00

(13) 0.37 0.51 (14) 0.37 0.32

(15) 0.37 0.34 (16) 0.37 0.34

(17) 0.75 0.92 (18) 0.75 0.77

(19) 0.75 0.73 (20) 0.75 0.74

(21) 0 0 (22) 0 0

(23) 0 0 (24) 0 0

Table I-18. TNT - Analysis of Target Found Concentration  
Points - GC/ECD Method

# ARMY CERTIFICATION GC/ECD METHOD

2,4,6 TNT

ANALYSIS OF 24 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.23166666667 SD= 0.266338566272

FOUND CONC

MEAN= 0.235 SD= 0.28672058121

N0. RUNS 1 TOTAL X-Y ALL RUNS 24 N0. CONCENTR 24  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.00911473269162

SLOPE= 1.053732605

R= 0.978826300535

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.00360103498674  
COMPUTE T

D.F.= 22

ENTER 2 TAIL P LEVEL (USUALLY .1) CONFID BANG IS

.05 SO TOTAL P= .1)

.1

t= 1.71713909197

REPLICATES ON UNKNOWN SAMPLE 3

y(c)= 0.0566954873641

x(d)= 0.122913743297



2,4,6 THT  
FOUND CONC




Table I-19. Tetryl-Found Concentrations

|                                 |  |
|---------------------------------|--|
| ARMY CERTIFICATION              |  |
| TETRYL                          |  |
| FOUND CONC                      |  |
| RUN 2                           |  |
| (14) 0 0 (15) 0 0               |  |
| (16) 0 0 (17) 0 0               |  |
| (18) 0.05 0.042 (19) 0.05 0.039 |  |
| (20) 0.05 0.039 (21) 0.05 0.036 |  |
| (22) 0.1 0.086 (23) 0.1 0.086   |  |
| (24) 0.1 0.083 (25) 0.1 0.083   |  |
| (26) 0.2 0.172 (27) 0.2 0.179   |  |
| (28) 0.2 0.166 (29) 0.2 0.256   |  |
| (30) 0.5 0.367 (31) 0.5 0.383   |  |
| (32) 0.5 0.364 (33) 0.5 0.37    |  |
| (34) 1 0.874 (35) 1 0.888       |  |
| (36) 1 0.872 (37) 1 0.885       |  |

Table I-20. Tetryl - Analysis of Target-Found Concentration Points

# ARMY CERTIFICATION

TETRYL

ANALYSIS OF 24 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.308333333333 SD= 0.356817657326

FOUND CONC

MEAN= 0.26125 SD= 0.309147108643

N0. RUNS 1 TOTAL X-Y ALL RUNS 24 N0. CONCENTR 24  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.00450455321571

SLOPE= 0.861906659078

R= 0.994812845018

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.00103307217364  
COMPUTE 1

D.F. = 22

ENTER 2 TAIL P LEVEL (USUALLY .1, EACH CONFID BAND IS  
.05 SO TOTAL P= .1)

.1

t= 1.71713909197

REPLICATES ON UNKNOWN SAMPLE 3

y(c)= 0.0307393706507

x(d)= 0.00103307514615

ARMY CERTIFICATION  
TETRYL  
FOUND CONC




Figure I-9. Tetryl - Graph of Target-Found Concentration Points

Table I-21. DPA - Found Concentrations

| ARMY CERTIFICATION |                          |
|--------------------|--------------------------|
| DPA                | FOUND CONC               |
| (1)                | 5 4.31 (2) 5 4.3         |
| (3)                | 10 7.5 (4) 10 7.6        |
| (5)                | 20 15.93 (6) 20 16.96    |
| (7)                | 50 35.01 (8) 50 35.1     |
| (9)                | 100 79.74 (10) 100 81.62 |

Table I-22. DPA- Analysis of Target-Found Concentration Points

ARMY CERTIFICATION

DPA ANALYSIS OF 10 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 37 SD= 37.0585122925

FOUND CONC  
MEAN= 28.807 SD= 29.5748620923

NO. RUNS 1 TOTAL X-Y ALL RUNS 10 NO. CONCENTR 10  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.633624595469  
SLOPE= 0.795692556634

USE FOR ACCURACY

R= 0.997835330176

MEAN SQR DEV OF POINTS FROM REGRESSION= 5.82586021445

ST ERROR EST= 2.4136818793

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 8

TWO TAIL P LEVEL IS .1

t= 1.85954407943

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U

REPLICATES= 1

Y(C)= 4.30510937847

X(D)= 12.2547388662

ARMY CERTIFICATION

DPA

FOUND CONC

81.52




Figure I-10. DPA -- Graph of Target-Found Concentration Points

Table I-23. 2,6-DNT - Found Concentrations -  
HPLC Method

**ARMY CERTIFICATION**

**26DNT**

**FOUND CONC**

**(1) 0 0 (2) 0 0**

**(3) 258 338 (4) 258 302**

**(5) 515 460 (6) 515 512**

**(7) 1030 1246 (8) 1030 1222**

**(9) 2575 2660 (10) 2575 2560**

**(11) 5150 5760 (12) 5150 5200**



Table I-24. 2,6-DNT - Analysis of Target-Found  
Concentration Points - HPLC Method

ARMY CERTIFICATION  
26DNT  
ANALYSIS OF 12 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 1588 SD= 1878.84635784

FOUND CONC  
MEAN= 1695 SD= 1998.49698869

NO. RUNS 1 TOTAL X-Y ALL RUNS 12 NO. CONCENTR 12  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 9.88878232945  
SLOPE= 1.86115315974  
USE FOR ACCURACY

R= 0.997621596909  
MEAN SQR DEV OF POINTS FROM REGRESSION= 20873.6483813  
ST ERROR EST= 144.477155223  
USE FOR PRECISION

T FOR CONFIDENCE BAND  
D.F.= 10

TWO TAIL P LEVEL IS .1  
t= 1.81245868646  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
y(c)= 290.49849589  
x(d)= 524.784244692 NG/2ML

ARMY CERTIFICATION  
26DNT  
FOUND CONC

5760.00




Figure I-11. 2,6-DNT - Graph of Target-Found  
Concentration Points - HPLC Method

Table I-25. 2,6-DNT-Found Concentrations -  
GC/ECD Method

ARMY CERTIFICATION GC/ECD METHOD

2,6 DNT

FOUND CONC

(1) 0.04 0.03 (2) 0.04 0.06

(3) 0.04 0 (4) 0.04 0

(5) 0.07 0.03 (6) 0.07 0.03

(7) 0.07 0.08 (8) 0.07 0.04

(9) 0.14 0.16 (10) 0.14 0.21

(11) 0.14 0.12 (12) 0.14 0.06

(13) 0.36 0.45 (14) 0.36 0.34

(15) 0.36 0.39 (16) 0.36 0.41

(17) 0.72 0.84 (18) 0.72 0.73

(19) 0.72 0.71 (20) 0.72 0.71

(21) 0 0 (22) 0 0

(23) 0 0 (24) 0 0

Table I-26. 2,6-DNT - Analysis of Target-Found  
Concentration Points

# ARMY CERTIFICATION GC/ECD METHOD

2,6 DNT

ANALYSIS OF 24 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 0.22166666667 SD= 0.256831304246

FOUND CONC

MEAN= 0.225 SD= 0.277567069842

N0. RUNS 1 TOTAL X-Y ALL RUNS 24 N0. CONCENTR 24  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.0119009535528

SLOPE= 1.06872610625

R= 0.988886469517

MEAN SQR DEV OF POINTS FROM REGRESSION (ST ERROR EST)= 0.00178834051461  
COMPUTE 1

D.F.= 22

ENTER 2 TAIL P LEVEL (USUALLY .1, EACH CONFID BAND IS  
.05 SO TOTAL P= .1)

.1

t= 1.71713909197

REPLICATES ON UNKNOWN SAMPLE 3

y(c)= 0.0343436272298

x(d)= 0.0054574529401

ARMY CERTIFICATION GC/ECD METHOD  
2,6-DNT  
FOUND CONC




Figure I-12. 2,6-DNT - Graph of Target-Found  
Concentration Points - GC/ECD Method

Table I-27. NG-Found Concentrations

| ARMY CERTIFICATION                    |            |
|---------------------------------------|------------|
| NG                                    | FOUND CONC |
| (1) 0 0 (2) 0 0                       |            |
| (3) 7600 10120 (4) 7600 10740         |            |
| (5) 15360 13940 (6) 15360 17760       |            |
| (7) 30720 37600 (8) 30720 34260       |            |
| (9) 76000 78600 (10) 76000 78800      |            |
| (11) 153600 153900 (12) 153600 152560 |            |

Table I-28. NG - Analysis of Target-Found  
Concentration Points

ARMY CERTIFICATION

NG

ANALYSIS OF 12 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 47360 SD= 56038.9682599

FOUND CONC

MEAN= 49043.3333333 SD= 55472.3674591

NO. RUNS 1 TOTAL X-Y ALL RUNS 12 NO. CONCENTR 12  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 2199.55036995

SLOPE= 0.989100147031

USE FOR ACCURACY

R= 0.999202923622

MEAN SQR DEV OF POINTS FROM REGRESSION= 5393900.17082

ST ERROR EST= 2322.4771626

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 10

TWO TAIL P LEVEL IS .1

t= 1.81245868646

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 6710.21417725

X(D)= 9075.83941875 NG/2ML

ARMY CERTIFICATION

NG  
FOUND CONC  
153980.00




Figure I-13. NG - Graph of Target-Found Concentration Points



Table I-29. PETN-Found Concentrations

| ARMY CERTIFICATION |                        |
|--------------------|------------------------|
| PETN               | FOUND CONC             |
| (1)                | 1.2 1.5 (2) 1.2 1.06   |
| (3)                | 1.2 1.73 (4) 1.2 2.32  |
| (5)                | 2.4 2.59 (6) 2.4 2.07  |
| (7)                | 2.4 3.52 (8) 2.4 3.42  |
| (9)                | 4.8 5.24 (10) 4.8 5.42 |
| (11)               | 4.8 5.98 (12) 4.8 7.07 |
| (13)               | 12 12.22 (14) 12 11.6  |
| (15)               | 12 13.2 (16) 12 12.56  |
| (17)               | 24 25.03 (18) 24 24.19 |
| (19)               | 24 26.41 (20) 24 27.34 |

Table I-30. PETH - Analysis of Target-Found  
Concentration Points

ARMY CERTIFICATION  
PETH  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 8.88 SD= 8.65682572188

FOUND CONC  
MEAN= 9.7235 SD= 9.0930327724

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.437787216828  
SLOPE= 1.04568837648  
USE FOR ACCURACY  
R= 0.99353506421  
MEAN SQR DEV OF POINTS FROM REGRESSION= 0.779368060054  
ST ERROR EST= 0.882818248596  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18  
TWO TAIL P LEVEL IS .1  
t= 1.73486096488  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U U  
REPLICATES= 3  
Y(C)= 1.36703361921  
X(D)= 1.76926827162

# ARMY CERTIFICATION

PETN  
FOUND CONC

27.34




Figure I-14. PETN - Graph of Target-Found Concentration Points

TABLE II-1  
QUANTITATIVE CERTIFICATION TESTING  
STATISTICAL DATA SUMMARY

|           | Metal                             |       | Concrete                          |                          | Brick                             |       | Transite                          |                          |
|-----------|-----------------------------------|-------|-----------------------------------|--------------------------|-----------------------------------|-------|-----------------------------------|--------------------------|
|           | D.L.<br>$\mu\text{g}/\text{cm}^2$ | % Rec | D.L.<br>$\mu\text{g}/\text{cm}^2$ | % Rec                    | D.L.<br>$\mu\text{g}/\text{cm}^2$ | % Rec | D.L.<br>$\mu\text{g}/\text{cm}^2$ | % Rec                    |
| DNP       | 0.33                              | 96%   | 1.74<br>(0.35)                    | 31%<br>34%) <sup>2</sup> | 1.59                              | 55%   | 2.22<br>(1.15)                    | 34%<br>30%) <sup>2</sup> |
| RDX       | 0.25                              | 96%   | 0.63                              | 78%                      | 2.11                              | 72%   | 3.48<br>(0.80)                    | 67%<br>82%) <sup>2</sup> |
| TNB       | 0.28                              | 95%   | 0.98                              | 75%                      | 2.12                              | 74%   | 3.46<br>(0.76)                    | 53%<br>71%) <sup>2</sup> |
| 2,4-DNT   | 0.90                              | 84%   | 1.08                              | 78%                      | 2.10                              | 69%   | 3.52<br>(0.62)                    | 65%<br>79%) <sup>2</sup> |
| 2,4,6-TNT | 0.60                              | 97%   | 1.57                              | 74%                      | 1.68                              | 66%   | 3.18<br>(0.68)                    | 51%<br>60%) <sup>2</sup> |
| Tetryl    | 1.95                              | 88%   | 4.29                              | 51%                      | 2.60                              | 68%   | 4.13                              | 56%                      |
| DPA       | 0.50                              | 94%   | 2.44<br>(1.15)                    | 74%<br>84%) <sup>2</sup> | 2.14                              | 69%   | 3.71<br>(0.90)                    | 67%<br>82%) <sup>2</sup> |
| 2,6-DNT   | 2.85                              | 94%   | 6.46                              | 86%                      | 6.36                              | 49%   | 2.04                              | 79%                      |
| NG        | 9.36                              | 94%   | 21.7                              | 76%                      | 32.5                              | 44%   | 26.1                              | 72%                      |
| PETN      | 10.2                              | 86%   | 5.39                              | 83%                      | 20.6                              | 62%   | 10.0                              | 82%                      |

<sup>1</sup> Calculated from four days of target vs. found concentrations using the procedures specified in the 1980 USATHAMA QA Plan.

<sup>2</sup> Calculated from three days of target vs. found concentrations.

QUANTITATIVE METHOD FOR THE DETERMINATION OF DNP,  
RDX, TNB, 2,4-DNT, TNT    ETRYL, AND DPA ON SURFACES

1. Application

Method used to extract the following compounds from metal, brick, concrete, transite surfaces:

2,4,-dinitrophenol    DNP  
cyclotrimethylenetrinitramine    RDX  
1,3,5-trinitrobenzene    TNB  
2,4-dinitrotoluene    2,4-DNT  
2,4,6-trinitrotoluene    2,4,6-TNT  
2,4,6-trinitrophenylmethylnitramine    Tetryl  
diphenylamine    DPA

A. Tested Concentration Range

|         |                                                                 |
|---------|-----------------------------------------------------------------|
| DNP     | 0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$ |
| RDX     | 0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$ |
| TNB     | 0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$ |
| 2,4-DNT | 0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$ |
| TNT     | 0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$ |
| Tetryl  | 0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$ |
| DPA     | 0.25 $\mu\text{g}/\text{cm}^2$ - 5.00 $\mu\text{g}/\text{cm}^2$ |

B. Sensitivity

Instrument response for each analyte is given below:

| <u>Analyte</u> | <u>Concentration</u> | <u>Response</u> |
|----------------|----------------------|-----------------|
| DNP            | 25.04 ng/mL          | 3700 area units |
| RDX            | 24.86 ng/mL          | 1388 area units |
| TNB            | 24.97 ng/mL          | 4437 area units |
| 2,4-DNT        | 24.99 ng/mL          | 6090 area units |
| TNT            | 24.91 ng/mL          | 4104 area units |
| Tetryl         | 25.30 ng/mL          | 2895 area units |
| DPA            | 24.91 ng/mL          | 2101 area units |

C. Detection Limit

See Table II-1.

D. Interferences

Interferences present in some brick and transite samples were apparently random rather than systematic. For example, the HPLC analysis of extraction B of one blank (unspiked) brick surface (extracted Feb 9)

indicated the presence of a compound with a retention time of 1510 seconds. This compound interfered with the TNT. However, this interference was not observed in any other blank brick sample. Therefore, only the brick samples analyzed on Feb 9 were corrected for the interference.

#### E. Analysis Rate

Six samples can be extracted and prepared for analysis in three hours. Rate of analysis is given below, excluding calibration standards:

|         |                            |
|---------|----------------------------|
| DNP     | 8 samples in an 8 hour day |
| RDX     | 8 samples in an 8 hour day |
| TNB     | 8 samples in an 8 hour day |
| 2,4-DNT | 8 samples in an 8 hour day |
| TNT     | 8 samples in an 8 hour day |
| Tetryl  | 8 samples in an 8 hour day |
| DPA     | 8 samples in an 8 hour day |

#### 2. Chemistry

2,4,-dinitrophenol    C<sub>6</sub>H<sub>4</sub>N<sub>2</sub>O<sub>5</sub>  
CAS RN 51-28-5  
MP 112-114C

Cyclotrimethylenetrinitramine    C<sub>3</sub>H<sub>6</sub>N<sub>6</sub>O<sub>6</sub>  
CAS RN 121-82-4  
MP 205-206C

1,3,5-trinitrobenzene    C<sub>6</sub>H<sub>3</sub>N<sub>3</sub>O<sub>6</sub>  
CAS RN 99-35-4  
MP 122.5C

2,4-dinitrotoluene    C<sub>7</sub>H<sub>6</sub>N<sub>2</sub>O<sub>4</sub>  
CAS RN 121-14-2  
MP 71C

2,4,6-trinitrotoluene    C<sub>7</sub>H<sub>5</sub>N<sub>3</sub>O<sub>6</sub>  
CAS RN 118-96-7  
MP 80.1C

2,4,6-trinitrophenylmethylnitramine  
CAS RN 479-45-8  
MP 130C Explodes 187C

diphenylamine    C<sub>12</sub>H<sub>11</sub>N  
CAS RN 122-39-4  
MP 53-54C BP 302C

### 3. Apparatus

#### A. Instrumentation

Waters Associates Model 6000A Solvent Delivery System  
Waters Associates Model M-45 Solvent Delivery System  
Waters Associates Model 660 Solvent Programmer  
Waters Associates Model 440 Absorbance Detector  
Waters Associates Intelligent Sample Processor (WISP)  
Spectra-Physics Minigrator  
Hewlett Packard 7133A Recorder

#### B. Parameters

Column: Spherisorb ODS 5 $\mu$ , 250 x 4.6 mm ID  
Precolumn: Pellicular LC-18, 40 $\mu$ , 50 x 4.6 mm ID  
Solvent System: linear gradient  
    Initial: 30/70 CH<sub>3</sub>CN/0.08 M acetic acid adjusted  
            to pH 3.1 with NH<sub>4</sub>OH  
    Final: 50/50 CH<sub>3</sub>CN/0.08 M acetic acid adjusted  
            to pH 3.1 with NH<sub>4</sub>OH  
    Time: 35 minutes  
Detector: UV at 254 nm  
Flow Rate: 1.0 mL/min  
Attenuation: 0.01 AUFS  
Injection Volume: 100  $\mu$ L

#### C. Hardware/Glassware

Westinghouse Ultrasonic cleaner  
8 ounce jars with teflon lined caps  
25 mL graduated cylinders  
microliter syringes  
volumetric flasks - 50, 10, 5 mL  
vials - WISP and 14 mL, with teflon lined caps

#### D. Chemicals

acetonitrile, HPLC grade  
nitrogen  
acetic acid  
ammonium hydroxide  
Standard Analytical Reference Material for each analyte

#### 4. Standards

##### A. Calibration Standard

Stock Solution A: DNP, RDX, TNB, 2,4-DNT, TNT, TETRYL, DPA  
Prepare individual stock solutions of 5.0 mg/mL. Combine  
500  $\mu$ L individual stocks and dilute to 5 mL.

##### Stock Solution B:

Dilute 625  $\mu$ L Stock Solution A to 5 mL. Concentration is 62.5  
 $\mu$ g/mL of each analyte

Calibration standards prepared in 25 mL volumetrics adding H<sub>2</sub>O  
so final solution is 55% H<sub>2</sub>O/45% CH<sub>3</sub>CN:

| <u>Cal Std.</u> | <u><math>\mu</math>L Stock B Added</u> | <u>Concentration Each Analyte</u> |
|-----------------|----------------------------------------|-----------------------------------|
| 1               | 10                                     | 25.0 $\mu$ g/L                    |
| 2               | 20                                     | 50.0 $\mu$ g/L                    |
| 3               | 40                                     | 100.0 $\mu$ g/L                   |
| 4               | 80                                     | 200.0 $\mu$ g/L                   |
| 5               | 200                                    | 500.0 $\mu$ g/L                   |
| 6               | 400                                    | 1000.0 $\mu$ g/L                  |

##### B. Control Spikes

Spiking stock solutions were prepared using 0.5 mg/mL stock solution.

| <u><math>\mu</math>L Stock Soln/x mL CH<sub>3</sub>CN</u> | <u>Concentration Each Analyte</u> |
|-----------------------------------------------------------|-----------------------------------|
| 125 $\mu$ L/10 mL                                         | 6.2 mg/L                          |
| 250 $\mu$ L/10 mL                                         | 12.0 mg/L                         |
| 500 $\mu$ L/10 mL                                         | 25.0 mg/L                         |
| 500 $\mu$ L/5 mL                                          | 50.0 mg/L                         |
| 625 $\mu$ L/5 mL                                          | 62.5 mg/L                         |
| 1250 $\mu$ L/5 mL                                         | 125.0 mg/L                        |
| 2500 $\mu$ L/5 mL                                         | 250.0 mg/L                        |

400  $\mu$ L 6.2, 12.5, 25.0, 62.5, 125.0 mg/L Stocks spiked onto 10 cm<sup>2</sup>  
concrete and metal surfaces.

200  $\mu$ L 12.5, 25.0, 50.0, 125.0, 250.0 mg/L stocks spiked onto 10 cm<sup>2</sup>  
brick and transite surfaces.

Concentration of analytes on surface after spiking: 2.5  $\mu$ g, 5.0  $\mu$ g,  
10.0  $\mu$ g, 25.0  $\mu$ g, 50  $\mu$ g.



## 5. Procedure

### Extraction A

1. Spike 10 cm<sup>2</sup> surface sample with acetonitrile spike solution (volume dependent on surface type). Allow solvent to evaporate.
2. Transfer sample to 8 ounce jar and add 20 mL CH<sub>3</sub>CN. Cover jar with teflon lined cap.
3. Sonicate for 10 minutes.
4. Transfer extract A to 50 mL volumetric flask, add 27 mL 0.08 M acetic acid and bring to volume with CH<sub>3</sub>CN. Save surface sample for extraction B.

Extract A ready for analysis.

### Extraction B

1. Add 20 mL CH<sub>3</sub>CN to jar with surface sample. Sonicate for 10 minutes.
2. Transfer surface sample to a second jar. Add 20 mL CH<sub>3</sub>CN and sonicate for 10 minutes more.
3. Combine extracts from steps 1 and 2 and evaporate using nitrogen to less than 5 mL.
4. Transfer evaporated extract to 10 mL volumetric flask. Add 5.5 mL to 0.08 M acetic acid and bring to volume with CH<sub>3</sub>CN.

Extract B ready for analysis.

## 6. Calculations

Calculate found concentration for each analyte in each sample extract from daily calibration data.

Multiply found concentration by extract volume to find total µg in extract. Combine total µg in extracts A and B to find total µg on surface.

## 7. References

Lakings, D.B., Baker, R.J., and Crook, M.V., "Precision and Accuracy Assessment of the High Performance Liquid Chromatographic Analytical Technique for the Determination of Dinitrophenol (DNP); Cyclotrimethylene trinitramine (RDX); 1,3-Dinitrobenzene (DNB); 1,3,5-Trinitrobenzene (TNB); 2,4-Dinitrotoluene (2,4-DNT); Trinitrotoluene (TNT); 2,4,6-Trinitrophenyl-methylnitramine (Tetryl); and Diphenylamine (DPA)", Midwest Research Institute Technical Report No. 1, USATHAMA Contract No. DAAK11-81-C-0007, March, 1981.

QUANTITATIVE METHOD FOR THE DETERMINATION OF  
2,6-DNT AND NG ON SURFACES

1. Application

Method used to extract the following compounds from metal, brick, concrete, transite surfaces:

2,6-dinitrotoluene    2,6-DNT  
nitroglycerine        NG

A. Tested Concentration Range

2,6-DNT        1.00  $\mu\text{g}/\text{cm}^2$  to 20.00  $\mu\text{g}/\text{cm}^2$   
NG              12.50  $\mu\text{g}/\text{cm}^2$  to 125.00  $\mu\text{g}/\text{cm}^2$

B. Sensitivity

Instrument response for each analyte is given below:

| <u>Analyte</u> | <u>Concentration</u>         | <u>Response</u>   |
|----------------|------------------------------|-------------------|
| 2,6-DNT        | 0.10 $\mu\text{g}/\text{mL}$ | 208750 area units |
| NG             | 1.25 $\mu\text{g}/\text{mL}$ | 186010 area units |

C. Detection Limit

See Table II-1.

D. Interferences

No interferences were observed

E. Analysis Rate

Six samples can be extracted and prepared for analysis in three hours. Rate of analysis is given below excluding calibration standards:

2,6-DNT              16 Samples in an 8 hour day  
NG                    16 samples in an 8 hour day

2. Chemistry

2,6-dinitrotoluene    C7H6N2O4  
CAS RN 606-20-2  
MP 66C

Nitroglycerine        C3H5N3O9  
CAS RN 55-63-0  
MP Stable form 13.5C

### 3. Apparatus

#### A. Instrumentation

Beckman Model 110A Solvent Metering Pump  
Waters Associates Model 450 Variable Wavelength Detector  
Waters Associates Model U6K Injector  
Hewlett Packard 3390A Integrator/Recorder

#### B. Parameters

Column: Spherisorb ODS, 5 $\mu$ , 250 x 4.6 mm ID  
Precolumn: Pellicular LC-18, 40  $\mu$ , 50 x 4.6 mm ID  
Solvent System: 35/65 CH<sub>3</sub>CN/0.005 M t-butyl ammonium hydroxide,  
pH 6.5. adjusted with 1N H<sub>3</sub>PO<sub>4</sub>  
Detector: UV at 230 nm  
Flow Rate: 1.0 mL/min  
Attenuation: 0.01 AUFS  
Injection Volume: 100  $\mu$ L

#### C. Hardware/Glassware

Westinghouse Ultrasonic cleaner  
8 ounce jars with teflon lined caps  
25 mL graduated cylinders  
microliter syringes  
volumetric flasks - 50, 10, 5 mL  
vials - WISP and 14 mL, with teflon lined caps

#### D. Chemicals

Acetonitrile, HPLC grade  
nitrogen  
phosphoric acid  
t-butyl ammonium hydroxide

Standard Analytical Reference Material for each analyte.

### 4. Standards

#### A. Calibration Standards

Prepare individual stock solutions:  
5 mg/mL 2,6-DNT  
50 mg/mL NG

**Stock Solution A:**

Combine 200  $\mu$ L 2,6-DNT stock and 250  $\mu$ L NG stock and dilute to 10 mL  $\text{CH}_3\text{CN}$ . Concentration is 0.1 mg/mL 2,6-DNT and 1.25 mg/mL NG.

Calibration standards prepared in 10 mL volumetric flasks adding 50%  $\text{H}_2\text{O}$ /50%  $\text{CH}_3\text{CN}$ .

| <u>Cal Std.</u> | <u><math>\mu</math>L Stock A added</u> | <u>2,6-DNT</u> | <u>NG</u>       |
|-----------------|----------------------------------------|----------------|-----------------|
| 1               | 10                                     | 0.1 $\mu$ g/mL | 1.2 $\mu$ g/mL  |
| 2               | 20                                     | 0.2 $\mu$ g/mL | 2.5 $\mu$ g/mL  |
| 3               | 40                                     | 0.4 $\mu$ g/mL | 5.0 $\mu$ g/mL  |
| 4               | 80                                     | 0.8 $\mu$ g/mL | 10.0 $\mu$ g/mL |
| 5               | 200                                    | 2.0 $\mu$ g/mL | 25.0 $\mu$ g/mL |

**B. Control Spikes**

Spike Solutions prepared following chart below:

| <u>Spike Solution</u> | <u>Amount Stock</u>                                    | <u>Dilute with <math>\text{CH}_3\text{CN}</math> to</u> | <u>Concentration</u>                 |
|-----------------------|--------------------------------------------------------|---------------------------------------------------------|--------------------------------------|
| 1                     | 1 mL of 5 mg/mL 2,6-DNT                                | 5 mL                                                    | 1 mg/mL 2,6-DNT                      |
| 2                     | 0.5 mL of 5 mg/mL 2,6-DNT                              | 5 mL                                                    | 0.5 mg/mL 2,6-DNT                    |
| 3                     | 0.5 mL of 5 mg/mL 2,6-DNT<br>0.625 mL of 50 mg/mL NG   | 5 mL                                                    | 0.5 mg/mL 2,6-DNT<br>6.25 mg/mL NG   |
| 4                     | 0.250 mL of 5 mg/mL 2,6-DNT<br>0.312 mL of 50 mg/mL NG | 5 mL                                                    | 0.25 mg/mL 2,6-DNT<br>3.12 mg/mL NG  |
| 5                     | 0.200 mL of 5 mg/mL 2,6-DNT<br>0.250 mL of 50 mg/mL NG | 5 mL                                                    | 0.20 mg/mL 2,6-DNT<br>2.50 mg/mL NG  |
| 6                     | 0.200 mL of 5 mg/mL 2,6-DNT<br>0.250 mL of 50 mg/mL NG | 10 mL                                                   | 0.10 mg/mL 2,6-DNT<br>1.25 mg/mL NG  |
| 7                     | 0.100 mL of 5 mg/mL 2,6-DNT<br>0.125 mL of 50 mg/mL NG | 10 mL                                                   | 0.05 mg/mL 2,6-DNT<br>0.62 mg/mL NG  |
| 8                     | 0.025 mL of 5 mg/mL 2,6-DNT<br>0.031 mL of 50 mg/mL NG | 5 mL                                                    | 0.025 mg/mL 2,6-DNT<br>0.31 mg/mL NG |

400  $\mu$ L of spike solutions 8, 7, 6, 4, 2 spiked onto 10  $\text{cm}^2$  concrete and metal surfaces

200  $\mu$ L of spike solutions 7, 6, 5, 3, 1 spiked onto 10  $\text{cm}^2$  brick and transite surfaces

Concentration on surface after spiking:

NG - 125 µg, 250 µg, 500 µg, 1250 µg

2,6-DNT - 10 µg, 20 µg, 40 µg, 100 µg, 200 µg.

## 5. Procedure

### Extraction A

1. Spike 10 cm<sup>2</sup> surface sample with acetonitrile spike solution (volume dependent on surface type). Allow solvent to evaporate.
2. Transfer sample to 8 ounce jar and add 20 mL CH<sub>3</sub>CN. Cover jar with teflon lined cap.
3. Sonicate for 10 minutes.
4. Transfer extract A to 50 mL volumetric flask, add 25 mL H<sub>2</sub>O and bring to volume with CH<sub>3</sub>CN.

Save surface sample for Extraction B.

### Extraction B

1. Add 20 mL CH<sub>3</sub>CN to jar with surface sample. Sonicate for 10 minutes.
2. Transfer surface sample to a second jar. Add 20 mL CH<sub>3</sub>CN and sonicate for 10 minutes more.
3. Combine extracts from steps 1 and 2 and evaporate using nitrogen to less than 5 mL.
4. Transfer evaporated extract to 10 mL volumetric flask. Add 5.0 mL H<sub>2</sub>O and bring to volume with CH<sub>3</sub>CN. Extract B ready for analysis.

## 6. Calculations

Calculate found concentration for each analyte in each sample extract from daily calibration data.

Multiply found concentration by extract volume to find total µg in extract. Combine total µg in extracts A and B to find total µg on surface.

## 7. References

Lakings, D.B., Baker, R.J., and Crook, M.V., "Precision and Accuracy Assessment of the High Performance Liquid Chromatographic Analytical Technique for the Determination of Nitrobenzene (NB), 2,6-Dinitrotoluene (2,6-DNT), Nitroglycerin (NG), and Picric Acid (PA), Midwest Research Institute Technical Report No. 2, USATHAMA Contract No. DAAK11-81-C-0007, May, 1981.

QUANTITATIVE METHOD FOR THE DETERMINATION  
OF PETN ON SURFACES

1. Application

Method used to extract pentaerythrite tetranitrate (PETN) from metal, brick, concrete, transite surfaces.

A. Tested Concentration Range:

PETN      5.0  $\mu\text{g}/\text{cm}^2$  to 100.0  $\mu\text{g}/\text{cm}^2$

B. Sensitivity

Instrument response for PETN is given below:

| <u>Concentration</u>         | <u>Response</u>  |
|------------------------------|------------------|
| 0.50 $\mu\text{g}/\text{mL}$ | 98025 area units |

C. Detection Limit

See Table II-1.

D. Interferences

There were no interferences.

E. Analysis Rate

Six samples can be extracted and prepared for analysis in three hours. Rate of analysis is given below, excluding calibration standards:

PETN              32 samples in an 8 hour day

2. Chemistry

Pentaerythrite tetranitrate C5H8N4O12  
CAS RN 78-11-5  
MP 140-141 C

3. Apparatus

A. Instrumentation

Beckman Model 110A Solvent Metering Pump  
Waters Associates Model 450 Variable Wavelength Detector  
Waters Associates Model U6K Injector  
Hewlett Packard 3390A Integrator/Recorder

#### B. Parameters

Column: Spherisorb ODS, 5 $\mu$ , 250 x 4.6 mm ID  
Precolumn: Pellicular LC-18, 40 $\mu$ , 50 x 4.6 mm ID  
Solvent System: 65% CH<sub>3</sub>CN/35% H<sub>2</sub>O  
Detector: UV at 230 nm  
Flow Rate: 1.0 mL/min  
Attenuation: 0.01 AUFS  
Injection Volume: 100  $\mu$ L

#### C. Hardware/Glassware

Westinghouse Ultrasonic cleaner  
8 ounce jars with teflon lined caps  
25 mL graduated cylinders  
microliter syringes  
volumetric flasks - 50, 10, 5 mL  
vials - WISP and 14 mL, with teflon lined caps

#### D. Chemicals

Acetonitrile, HPLC grade  
Standard Analytical Reference Material for PETN

#### 4. Standards

##### A. Calibration Standards:

Prepare stock solution as follows:  
200  $\mu$ L of SARM (50 mg/mL) in 10 mL CH<sub>3</sub>CN = 1.0 mg/mL

Calibration Standards prepared in 10 mL volumetric flasks adding  
50% H<sub>2</sub>O/50% CH<sub>3</sub>CN.

| <u>Cal. Std</u> | <u><math>\mu</math>L Stock added</u> | <u>Concentration PETN</u> |
|-----------------|--------------------------------------|---------------------------|
| 1               | 5                                    | 0.5 $\mu$ g/mL            |
| 2               | 10                                   | 1.0 $\mu$ g/mL            |
| 3               | 20                                   | 2.0 $\mu$ g/mL            |
| 4               | 40                                   | 4.0 $\mu$ g/mL            |
| 5               | 100                                  | 10.0 $\mu$ g/mL           |
| 6               | 200                                  | 20.0 $\mu$ g/mL           |

## B. Control Spikes

Spike solutions prepared following chart below:

| <u>Spike Solution</u> | <u>Amount Stock</u>      | <u>Dilute with<br/>CH<sub>3</sub>CN to</u> | <u>Concentration</u> |
|-----------------------|--------------------------|--------------------------------------------|----------------------|
| 1                     | 12.5 $\mu$ L of 50 mg/mL | 5 mL                                       | 0.125 mg/mL          |
| 2                     | 50 $\mu$ L of 50 mg/mL   | 10 mL                                      | 0.25 mg/mL           |
| 3                     | 100 $\mu$ L of 50 mg/mL  | 10 mL                                      | 0.50 mg/mL           |
| 4                     | 100 $\mu$ L of 50 mg/mL  | 5 mL                                       | 1.0 mg/mL            |
| 5                     | 125 $\mu$ L of 50 mg/mL  | 5 mL                                       | 1.25 mg/mL           |
| 6                     | 500 $\mu$ L of 50 mg/mL  | 10 mL                                      | 2.5 mg/mL            |
| 7                     | 500 $\mu$ L of 50 mg/mL  | 5 mL                                       | 5.0 mg/mL            |

400  $\mu$ L spike solutions 1,2,3,5 and 6 spiked onto 10 cm<sup>2</sup> metal surfaces

200  $\mu$ L spike solutions 2,3,4,6,7 spiked onto 10 cm<sup>2</sup> concrete, transite, and brick surfaces.

Concentration of analytes on surface after spiking: 50  $\mu$ g, 100  $\mu$ g, 200  $\mu$ g, 500  $\mu$ g, 1000  $\mu$ g.

## 5. Procedure

1. Spike 10 cm<sup>2</sup> surface sample with acetonitrile spike solution (volume dependent on surface type). Allow solvent to evaporate.
2. Transfer sample to 8 ounce jar.
3. Add 12 mL CH<sub>3</sub>CN. Cover jar with teflon lined cap. Sonicate for 10 min.
4. Transfer extract to 50 mL volumetric flask.
5. Repeat steps 3 and 4 twice, adding extracts to same 50 mL volumetric flask.
6. Rinse jar with 17 mL H<sub>2</sub>O and transfer to volumetric flask.
7. Bring to volume with CH<sub>3</sub>CN.
8. Ready for HPLC analysis.

## 6. Calculations

Calculate found concentration for each analyte in each sample extract from daily calibration data.

Multiply found concentration by extract volume to find total  $\mu$ g in extract.

## 7. References

None



Table II-2. DNP on Metal - Target US. Found Concentrations

| 2,4-DINITROPHENOL (24DNP) |                                      |
|---------------------------|--------------------------------------|
| METAL SURFACE             |                                      |
| TARGET CONC. US           | FOUND CONC.                          |
| ug/10 sq cm               | ug/10 sq cm                          |
| 2.500                     | 1.920<br>2.270<br>2.260<br>2.470     |
| 5.000                     | 3.700<br>4.030<br>4.610<br>5.100     |
| 10.000                    | 9.060<br>8.020<br>10.040<br>11.220   |
| 25.000                    | 23.500<br>22.460<br>25.240<br>24.610 |
| 50.000                    | 46.870<br>47.850<br>48.840<br>48.220 |

Table II-3. DNP on Metal - Analysis of Target-  
Found Concentration Points

2,4-DINITROPHENOL (24DNP)  
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 17.6145 SD= 17.4133452982

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.226280339806

SLOPE= 0.964366504854

USE FOR ACCURACY

R= 0.998797261237

MEAN SQR DEV OF POINTS FROM REGRESSION= 0.769459159271

ST ERROR EST= 0.877188212

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 1.37295729236

X(D)= 3.30309425069

2,4-DINITROPHENOL (24DNP)  
METAL SURFACE  
FOUND CONC




Figure 11-1. DNP on Metal - Graph of Target-Found Concentration Points

Table II-4. DYP on Metal - Inaccuracy and Imprecision Data

**2,4-DINITROPHENOL (24DNP)  
METAL SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION**

| Target Conc<br>ug/10 sq cm | Mean Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|----------------------------|--------------------------------|-----------------------|------------------------|-------------|
| 2.500                      | 2.230                          | 0.228                 | -10.800                | 10.232      |
| 5.000                      | 4.360                          | 0.620                 | -12.900                | 14.229      |
| 10.000                     | 9.585                          | 1.367                 | -4.150                 | 14.261      |
| 25.000                     | 23.953                         | 1.228                 | -4.190                 | 5.126       |
| 50.000                     | 47.945                         | 0.825                 | -4.110                 | 1.720       |
| Means                      |                                | 0.854                 | -7.210                 | 9.114       |

2,4-DINITRUPHENOL (24DNP)

METAL SURFACE

MEAN INACCURACY

MEAN FOUND CONCENTRATION (REPORT)

-4.11

0

-12.80

10

20

30

40

50

TARGET CONC

0.869

VERTICAL AXIS TIC INTERVAL = 0.869

Figure II-2. DNP on Metal - Graph of Inaccuracy

Arthur D. Little, Inc.

2,4-DINITROPHENOL (24DNP)  
 METAL SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-3. DNP on Metal - Graph of Imprecision

Table II-5. DNP on Concrete - Target U.S. Found Concentrations

| 2,4-DINITROPHENOL (24DNP)<br>CONCRETE SURFACE<br>TARGET CONC. US FOUND CONC |                                      |
|-----------------------------------------------------------------------------|--------------------------------------|
| Target Conc<br>ug/10 sq cm                                                  | Found Conc<br>ug/10 sq cm            |
| 2.500                                                                       | 0.710<br>1.010<br>1.470<br>1.250     |
| 5.000                                                                       | 1.320<br>1.890<br>2.690<br>1.920     |
| 16.000                                                                      | 2.270<br>3.500<br>3.800<br>3.290     |
| 25.000                                                                      | 6.210<br>8.430<br>8.980<br>8.560     |
| 50.000                                                                      | 10.800<br>17.570<br>17.540<br>16.980 |

Table II-6. DNP on Concrete - Analysis of Target-  
Found Concentration Points

2,4-DINITROPHENOL (24DNP)  
CONCRETE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535072

FOUND CONC  
MEAN= 6.0095 SD= 5.74332424029

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.308970469256  
SLOPE= 0.308136731392  
USE FOR ACCURACY  
R= 0.967603816639  
MEAN SQR DEV OF POINTS FROM REGRESSION= 2.219410059  
ST ERROR EST= 1.48977141166  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18  
TWO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
Y(C)= 3.02503353561  
X(D)= 17.4060970424





Table II-7. DNP on Concrete - Inaccuracy and Precision Data

| 2,4-DINITROPHENOL (24DNP)<br>CONCRETE SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                              |                       |                        |             |  |
|---------------------------------------------------------------------------------------------------------------------------|------------------------------|-----------------------|------------------------|-------------|--|
| Mn Target Con<br>ug/10 sq cm                                                                                              | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |  |
| 2.500                                                                                                                     | 1.110                        | 0.326                 | -55.600                | 29.386      |  |
| 5.000                                                                                                                     | 1.955                        | 0.562                 | -60.900                | 20.768      |  |
| 10.000                                                                                                                    | 3.215                        | 0.664                 | -67.830                | 20.649      |  |
| 25.000                                                                                                                    | 8.045                        | 1.246                 | -67.820                | 15.483      |  |
| 50.000                                                                                                                    | 15.723                       | 3.293                 | -68.555                | 20.944      |  |
| Means                                                                                                                     |                              | 1.218                 | -64.145                | 23.046      |  |

2,4-DINITROPHENOL (24DNP)  
 CONCRETE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-5. DNP on Concrete - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-6. DNP on Concrete - Graph of Imprecision

Table II-8. DNP on Concrete (3 days) - Target vs.  
Found Concentrations

| 2,4-DINITROPHENOL (24DNP)<br>CONCRETE SURFACE |                            |
|-----------------------------------------------|----------------------------|
| TARGET CONC. US FOUND CONC                    |                            |
| Target Conc<br>ug/10 sq cm                    | Found Conc<br>ug/10 sq cm  |
| 2.500                                         | 1.010<br>1.470<br>1.250    |
| 5.000                                         | 1.890<br>2.690<br>1.920    |
| 10.000                                        | 3.500<br>3.800<br>3.290    |
| 25.000                                        | 8.430<br>8.980<br>8.560    |
| 50.000                                        | 17.570<br>17.540<br>16.980 |

Table II-9. DNP on Concrete (3 days) - Analysis of Target-Found Concentration Points

2,4-DINITROPHENOL (24DNP)  
CONCRETE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.1953683274

FOUND CONC  
MEAN= 6.592 SD= 6.17940727625

NØ. RUN: 1 TOTAL X-Y ALL RUNS 15 NØ. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.31678802589

SLOPE= 0.339200647249

USE FOR ACCURACY

R= 0.998781992786

MEAN SQR DEV OF POINTS FROM REGRESSION= 0.100113722929

ST ERROR EST= 0.316407326662

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 0.915196314813

X(D)= 3.50862633373

2,4-DINITROPHENOL (24DNP)  
CONCRETE SURFACE  
FOUND CONC

17.57

$y(c)$

$y(d)$

3

0.00

VERTICAL AXIS TIC INTERVAL= 1.757

TARGET CONC

50

40

30

20

10

0

Figure II-7. DNP on Concrete (3 days) - Graph of Target-Found Concentration Points

Table II-10. DNP on Concrete (3 days) - Inaccuracy and Imprecision Data

| 2,4-DINITROPHENOL (24DNP)<br>CONCRETE SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                              |                       |                        |             |  |
|---------------------------------------------------------------------------------------------------------------------------|------------------------------|-----------------------|------------------------|-------------|--|
| Mn Target Con<br>ug/10 sq cm                                                                                              | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |  |
| 2.500                                                                                                                     | 1.243                        | 0.230                 | -50.267                | 18.504      |  |
| 5.000                                                                                                                     | 2.167                        | 0.453                 | -56.667                | 20.929      |  |
| 10.000                                                                                                                    | 3.530                        | 0.256                 | -64.703                | 7.261       |  |
| 25.000                                                                                                                    | 8.657                        | 0.287                 | -65.373                | 3.321       |  |
| 50.000                                                                                                                    | 17.363                       | 0.332                 | -65.273                | 1.914       |  |
| Means                                                                                                                     |                              | 0.312                 | -60.456                | 10.386      |  |



2,4-DINITROPHENOL (24DNP)  
 CONCRETE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-8. DNP on Concrete (3 days) - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-9. DNP on Concrete (3 days) - Graph of Imprecision

Table II-11. DNP on Brick - Target vs. Found Concentrations

| 2,4-DINITROPHENOL (24DNP)<br>BRICK SURFACE |                                              |
|--------------------------------------------|----------------------------------------------|
| TARGET CONC.<br>ug/10 sq cm                | VS. FOUND CONC.<br>Found Conc<br>ug/10 sq cm |
| 2.500                                      | 1.000<br>1.360<br>0.810<br>1.110             |
| 5.000                                      | 1.690<br>4.000<br>2.780<br>3.620             |
| 10.000                                     | 5.560<br>3.500<br>4.890<br>6.800             |
| 25.000                                     | 13.450<br>10.200<br>11.320<br>14.100         |
| 50.000                                     | 28.950<br>34.090<br>21.730<br>26.490         |

Table II-12. DNP on Brick - Analysis of Target-Found Concentration Points

2,4-DINITROPHENOL (24DNP)  
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 9.8725 SD= 18.2455330148

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.352117718446  
SLOPE= 0.532682038835

USE FOR ACCURACY  
R= 0.972875761094

MEAN SQR DEV OF POINTS FROM REGRESSION= 5.92937890036  
ST ERROR EST= 2.43503160151

USE FOR PRECISION  
T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 4.08728778129

X(D)= 15.86533395681



Table II-13. DNP on Brick - Inaccuracy and Imprecision Data

**2,4-DINITROPHENOL (24DNP)**

**BRICK SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION**

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 2.500                        | 0.230                 | -57.200                | 21.461      |
| 5.000                        | 3.023                        | 1.024                 | -39.550                | 33.885      |
| 10.000                       | 5.188                        | 1.375                 | -48.125                | 26.514      |
| 25.000                       | 12.268                       | 1.819                 | -50.930                | 14.830      |
| 50.000                       | 27.815                       | 5.146                 | -44.370                | 18.501      |
| Means                        |                              | 1.919                 | -48.035                | 23.030      |

2,4-DINITROPHENOL (24DNP)  
 BRICK SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-11. DNP on Brick - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)  
 BRICK SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-12. DNP on Brick - Graph of Imprecision



Table II-14. DNP on Transite - Target vs. Found Concentrations

| 2,4-DINITROPHENOL (24DNP)   |                           |
|-----------------------------|---------------------------|
| TRANSITE SURFACE            |                           |
| TARGET CONC. VS. FOUND CONC |                           |
| Target Conc<br>ug/10 sq cm  | Found Conc<br>ug/10 sq cm |
| 2.500                       | 0.980                     |
|                             | 0.710                     |
|                             | 0.250                     |
|                             | 1.280                     |
| 5.000                       | 2.260                     |
|                             | 1.440                     |
|                             | 1.200                     |
|                             | 1.970                     |
| 10.000                      | 3.880                     |
|                             | 2.380                     |
|                             | 2.790                     |
|                             | 2.630                     |
| 25.000                      | 9.700                     |
|                             | 7.140                     |
|                             | 6.960                     |
|                             | 7.050                     |
| 50.000                      | 24.050                    |
|                             | 16.740                    |
|                             | 12.550                    |
|                             | 15.680                    |

Table II-15. DNP on Transite - Analysis of Target-  
Found Concentration Points

2,4-DINITROPHENOL (24DNP)  
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 6.002 SD= 6.55535747148

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.297536407767

SLOPE= 0.344839805825

USE FOR ACCURACY

R= 0.948720859866

MEAN SQR DEV OF POINTS FROM REGRESSION= 4.53277563377

ST ERROR EST= 2.12903161078

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C


Y(C)= 3.58398803391

X(D)= 22.2382830955

2,4-DINITROPHENOL (24DNP)  
TRANSITE SURFACE  
FOUND CONC

24.05

0



TARGET CONC

VERTICAL AXIS TIC INTERVAL = 2.405

Figure II-13. DNP on Transite - Graph of Target-Found Concentration Points

Table II-16. DNP on Transite - Inaccuracy and Imprecision Data

2,4-DINITROPHENOL (24DNP)

TRANSITE SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 0.805                        | 0.437                 | -67.800                | 54.304      |
| 5.000                        | 1.718                        | 0.484                 | -65.650                | 28.183      |
| 10.000                       | 2.920                        | 0.662                 | -70.800                | 22.667      |
| 25.000                       | 7.713                        | 1.327                 | -69.150                | 17.206      |
| 50.000                       | 17.235                       | 4.867                 | -65.490                | 28.205      |
| Means                        |                              | 1.555                 | -67.770                | 30.113      |

2,4-DINITROPHENOL (24DNP)  
TRANSITE SURFACE  
MEAN INACCURACY

-65.49

MEAN FOUND CONCENTRATION (REPORT)




Figure II-14. DNP on Transite - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-15. DNP on Transite - Graph of Imprecision

Table II-17. DNP on Transite (3 days) - Target vs. Found Concentrations

| 2,4-DINITROPHENOL (24DNP)  |                            |
|----------------------------|----------------------------|
| TRANSITE SURFACE           |                            |
| TARGET CONC. US FOUND CONC |                            |
| ug/10 sq cm                | ug/10 sq cm                |
| 2.500                      | 0.710<br>0.250<br>1.280    |
| 5.000                      | 1.440<br>1.200<br>1.970    |
| 10.000                     | 2.380<br>2.790<br>2.630    |
| 25.000                     | 7.140<br>6.960<br>7.030    |
| 50.000                     | 16.740<br>12.550<br>15.680 |

Table II-18. DNP on Transite (3 days) - Analysis of Target-Found Concentration Points

2,4-DINITROPHENOL (24DNP)  
TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.1953683274

FOUND CONC  
MEAN= 5.3846666667 SD= 5.52760328671

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.161968623674

SLOPE= 0.299817691478

USE FOR ACCURACY

R= 0.986918388416

MEAN SQ DEV OF POINTS FROM REGRESSION= 0.855263021743  
ST ERROR EST= 0.924804315378

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 1.5870829987

X(D)= 11.5032476202





Table II-19. DNP on Transite (3 days) - Inaccuracy and Imprecision Data

| 2,4-DINITROPHENOL (24DNP)<br>TRANSITE SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                              |                       |                        |             |
|---------------------------------------------------------------------------------------------------------------------------|------------------------------|-----------------------|------------------------|-------------|
| Mn Target Conc<br>ug/10 sq cm                                                                                             | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
| 2.500                                                                                                                     | 0.747                        | 0.516                 | -70.133                | 69.104      |
| 5.000                                                                                                                     | 1.537                        | 0.394                 | -69.267                | 25.640      |
| 10.000                                                                                                                    | 2.600                        | 0.207                 | -74.000                | 7.948       |
| 25.000                                                                                                                    | 7.050                        | 0.090                 | -71.800                | 1.277       |
| 50.000                                                                                                                    | 14.990                       | 2.179                 | -70.020                | 14.533      |
| Means                                                                                                                     |                              | 0.677                 | -71.044                | 23.700      |

ACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-17. DNP on Transit (3 days) - Graph of Inaccuracy

2,4-DINITROPHENOL (24DNP)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-18. DNP on Transite (3 days) - Graph of Imprecision

Table II-20. RDX on Metal - Target vs. Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

| METAL SURFACE            |                                       |
|--------------------------|---------------------------------------|
| TARGET CONC. ug/10 sq cm | VS FOUND CONC. Found Conc ug/10 sq cm |
| 2.500                    | 2.340                                 |
|                          | 2.040                                 |
|                          | 2.180                                 |
|                          | 2.620                                 |
| 5.000                    | 4.100                                 |
|                          | 4.250                                 |
|                          | 4.840                                 |
|                          | 5.070                                 |
| 25.000                   | 24.180                                |
|                          | 23.960                                |
|                          | 23.070                                |
|                          | 24.460                                |
| 50.000                   | 49.060                                |
|                          | 46.620                                |
|                          | 46.960                                |
|                          | 48.510                                |

Table II-21. RDX on Metal - Analysis of Target-Found Concentrations Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
METAL SURFACE

ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 20.625 SD= 19.6956001855

FOUND CONC  
MEAN= 19.64125 SD= 18.9047401022

N0. RUNS 1 TOTAL X-Y ALL RUNS 16 N0. CONCENTR 16  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.144172932331

SLOPE= 0.959293233083

USE FOR ACCURACY

R= 0.999424264881

MEAN SQR DEV OF POINTS FROM REGRESSION= 0.440790601532  
ST ERROR EST= 0.66392062894

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 1.10196465687

X(D)= 2.58815970399

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
METAL SURFACE  
FOUND CONC




Figure II-19. RDX on Metal - Graph of Target-  
Found Concentration Points

Table II-22. RDX on Metal - Inaccuracy and Imprecision Data

**CYCLOTRIMETHYLENETRINITRAMINE (RDX)**

**METAL SURFACE**

**STATISTICAL DATA USED TO DETERMINE PERCENT**

**INACCURACY AND IMPRECISION**

| Mn Target Conc<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                         | 2.295                        | 0.249                 | -8.200                 | 10.847      |
| 5.000                         | 4.565                        | 0.464                 | -8.700                 | 10.166      |
| 25.000                        | 23.918                       | 0.601                 | -4.330                 | 2.512       |
| 50.000                        | 47.788                       | 1.182                 | -4.425                 | 2.473       |
| Means                         |                              | 0.624                 | -6.414                 | 6.499       |






Figure 11-20. RDX on Metal - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 METAL SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-21. RDX on Metal - Graph of Imprecision

Table II-23. RDX on Concrete - Target vs. Found Concentrations

CYCLOTRIMETHYLENE TRINITRAMINE (RDX)

| CONCRETE SURFACE         |                         |
|--------------------------|-------------------------|
| TARGET CONC. us/10 sq cm | FOUND CONC. us/10 sq cm |
| 2.500                    | 1.610                   |
|                          | 1.740                   |
|                          | 1.950                   |
|                          | 2.090                   |
| 5.000                    | 2.860                   |
|                          | 3.710                   |
|                          | 4.320                   |
|                          | 4.210                   |
| 10.000                   | 6.730                   |
|                          | 7.650                   |
|                          | 7.520                   |
|                          | 7.930                   |
| 25.000                   | 19.370                  |
|                          | 20.410                  |
|                          | 20.410                  |
|                          | 17.490                  |
| 50.000                   | 37.810                  |
|                          | 41.120                  |
|                          | 41.240                  |
|                          | 35.160                  |

Table II-24. RDX on Concrete - Analysis of Target-  
Found Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
CONCRETE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 14.2665 SD= 14.1353562966

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.168798543689

SLOPE= 0.780286407767

USE FOR ACCURACY

R= 0.99553764773

MEAN SQR DEV OF POINTS FROM REGRESSION= 1.8713384477

ST ERROR EST= 1.36796580648

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 2.32519575231

x(d)= 6.34666977015

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
CONCRETE SURFACE  
FOUND CONC




Figure II-22. RDX on Concrete - Graph of Target-Found Concentration Points

Table II-25. Inaccuracy and Imprecision Data

**CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
CONCRETE SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION**

| Mn Target Conc<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                         | 1.848                        | 0.214                 | -26.100                | 11.570      |
| 5.000                         | 3.775                        | 0.665                 | -24.500                | 17.623      |
| 10.000                        | 7.450                        | 0.514                 | -25.425                | 6.896       |
| 25.000                        | 19.420                       | 1.377                 | -22.320                | 7.090       |
| 50.000                        | 38.833                       | 2.919                 | -22.335                | 7.517       |
| Means                         |                              | 1.138                 | -24.136                | 10.141      |




Figure 11-23. RDX on Concrete - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-24. RDX on Concrete - Graph of Imprecision



Table II-26. RDX on Brick - Target vs. Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

BRICK SURFACE

TARGET CONC. VS. FOUND CONC.

| Target Conc<br>ug/10 sq cm | Found Conc<br>ug/10 sq cm            |
|----------------------------|--------------------------------------|
| 2.500                      | 1.340<br>2.070<br>1.020<br>1.160     |
| 5.000                      | 2.450<br>7.380<br>4.120<br>5.010     |
| 10.000                     | 9.080<br>6.410<br>7.190<br>9.900     |
| 25.000                     | 21.180<br>14.890<br>14.480<br>20.990 |
| 50.000                     | 36.030<br>45.030<br>24.020<br>41.020 |

Table II-27. RDX on Brick - Analysis of Target-Found  
Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
BRICK SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 13.7385 SD= 13.6021253274

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.432779126214

SLOPE= 0.71922815534

USE FOR ACCURACY

R= 0.95362437934

MEAN SQR DEV OF POINTS FROM REGRESSION= 17.6939762945

ST ERROR EST= 4.20642084135

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(0) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U CB

y(c)= 8.10167664923

x(d)= 21.0600746217

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
BRICK SURFACE  
FOUND CONC




Figure 14-25. RDX on Brick - Graph of Target-Found Concentration Points

Table 11-28. RDX on Brick - Inaccuracy and Imprecision Data

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
BRICK SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.398                        | 0.467                 | -44.100                | 33.422      |
| 5.000                        | 4.740                        | 2.055                 | -5.200                 | 43.358      |
| 10.000                       | 8.145                        | 1.620                 | -18.550                | 19.894      |
| 25.000                       | 17.885                       | 3.700                 | -28.460                | 20.686      |
| 50.000                       | 36.525                       | 9.113                 | -26.950                | 24.951      |
| Means                        |                              | 3.391                 | -24.652                | 28.462      |

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 BRICK SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-26. RDX on Brick - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 BRICK SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-27. RDX on Brick - Graph of Imprecision

Table II-29. RDX on Transite - Target vs. Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

| TRANSITE SURFACE           |        |
|----------------------------|--------|
| TARGET CONC. US FOUND CONC |        |
| Target Conc Found Conc     |        |
| ug/10 sq cm ug/10 sq cm    |        |
| 2.500                      | 0.710  |
|                            | 1.640  |
|                            | 1.570  |
|                            | 1.410  |
| 5.000                      | 1.610  |
|                            | 3.900  |
|                            | 3.670  |
| 10.000                     | 3.480  |
|                            | 2.630  |
|                            | 7.590  |
| 25.000                     | 7.150  |
|                            | 7.840  |
|                            | 6.640  |
|                            | 15.830 |
| 50.000                     | 19.870 |
|                            | 18.270 |
|                            | 12.540 |
|                            | 39.670 |
|                            | 44.350 |
|                            | 38.070 |

Table 11-30. RDX on Transite - Analysis of Target-  
Found Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.3 SD= 18.0350535872

FOUND CONC  
MEAN= 11.922 SD= 13.6804011396

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.544051375404  
SLOPE= 0.673840614887

USE FOR ACCURACY

R= 0.888327791599

MEAN SQR DEV OF POINTS FROM REGRESSION= 41.6587590559  
ST ERROR EST= 6.45435969372

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 11.2231552201

x(d)= 34.829209243






Figure 11-28. RDX on Transit - Graph of Target-Found Concentration Points

Table II-31. RDX on Transite - Inaccuracy and Imprecision Data

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 TRANSITE SURFACE  
 STATISTICAL DATA USED TO DETERMINE PERCENT  
 INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.333                        | 0.426                 | -46.700                | 31.971      |
| 5.000                        | 3.165                        | 1.051                 | -36.700                | 33.200      |
| 10.000                       | 6.303                        | 2.465                 | -36.975                | 39.110      |
| 25.000                       | 15.153                       | 5.913                 | -39.390                | 39.024      |
| 50.000                       | 33.658                       | 14.328                | -32.685                | 42.571      |
| Means                        |                              | 4.837                 | -38.490                | 37.175      |

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 TRANSITE SURFACE  
 MEAN FOUND CONCENTRATION (REPORT)  
 MEAN INACCURACY

-32.69

-46.70

0

10

20

30

40

50

TARGET CONC

INTERVAL = 1.4015

VERTICAL AXIS TIC

Figure 11-29. RDX on Transite - Graph of Inaccuracy

Arthur D Little Inc.

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-30. RDX on Transite - Graph of Imprecision

Table II-32. RDX on Transite (3 days) - Target vs.  
Found Concentrations

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

| TRANSITE SURFACE<br>TARGET CONC. VS. FOUND CONC. |                           |
|--------------------------------------------------|---------------------------|
| Target Conc<br>ug/10 sq cm                       | Found Conc<br>ug/10 sq cm |
| 2.500                                            | 1.640                     |
|                                                  | 1.570                     |
|                                                  | 1.410                     |
| 5.000                                            | 3.900                     |
|                                                  | 3.670                     |
|                                                  | 3.480                     |
| 10.000                                           | 7.590                     |
|                                                  | 7.150                     |
|                                                  | 7.840                     |
| 25.000                                           | 15.830                    |
|                                                  | 19.870                    |
|                                                  | 18.270                    |
| 50.000                                           | 39.670                    |
|                                                  | 44.350                    |
|                                                  | 38.070                    |

Table II-33. RDX on Transite (3 days) - Analysis of Target-Found Concentration Points

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.1953683274

FOUND CONC  
MEAN= 14.287333333 SD= 14.9412293751

NØ. RUNS 1 TOTAL X-Y ALL RUNS 15 NØ. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.886491370011  
SLOPE= 0.815882416397

USE FOR ACCURACY

R= 0.993578286336

MEAN SQR DEV OF POINTS FROM REGRESSION= 3.07780840387  
ST ERROR EST= 1.75436837747

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 2.51147198116

X(D)= 8.04284789281

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
TRANSITE SURFACE  
FOUND CONC




Figure 11-31. RDX on Transite (3 days) - Graph of Target-Found Concentration Points

Table II-34. RDX on Transite (3 days) - Inaccuracy and Imprecision Data

CYCLOTRIMETHYLENETRINITRAMINE (RDX)

STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION

| Mn Target Con ug/10 sq cm | Mn Found Conc ug/10 sq cm | Standard Deviation | Mean Pct Inaccuracy | Imprecision |
|---------------------------|---------------------------|--------------------|---------------------|-------------|
| 2.500                     | 1.540                     | 0.118              | -38.400             | 7.656       |
| 5.000                     | 3.683                     | 0.210              | -26.333             | 5.710       |
| 10.000                    | 7.527                     | 0.349              | -24.733             | 4.641       |
| 25.000                    | 17.990                    | 2.035              | -28.040             | 11.309      |
| 50.000                    | 40.697                    | 3.263              | -18.607             | 8.019       |
| Means                     |                           | 1.195              | -27.223             | 7.467       |



CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 TRANSITE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-32. RDX on Transite (3 days) - Graph of Inaccuracy

CYCLOTRIMETHYLENETRINITRAMINE (RDX)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-33. RDX on Transite (3 days) - Graph of Imprecision

Table II-35. TNB on Metal - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

| METAL SURFACE               |                                            |
|-----------------------------|--------------------------------------------|
| TARGET CONC.<br>ug/10 sq cm | VS FOUND CONC<br>Found Conc<br>ug/10 sq cm |
| 2.500                       | 2.130<br>2.290<br>2.280<br>2.470           |
| 5.000                       | 3.690<br>4.460<br>4.740<br>4.820           |
| 10.000                      | 9.380<br>8.150<br>9.600<br>10.190          |
| 25.000                      | 23.560<br>23.480<br>25.550<br>24.340       |
| 50.000                      | 46.760<br>47.630<br>46.340<br>48.390       |

Table II-36. TNB on Metal - Analysis of Target-Found  
Concentration Points

1,3,5-TRINITROBENZENE (135TNB)  
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 17.5125 SD= 17.1872997405

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.102313915858  
SLOPE= 0.95215210356

USE FOR ACCURACY  
R= 0.999116467986

MEAN SQR DEV OF POINTS FROM REGRESSION= 0.550752890156  
ST ERROR EST= 0.742127273556

USE FOR PRECISION  
T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 1.25068860685

x(d)= 2.83188868416

1,3,5-TRINITROBENZENE (135TNB)  
METAL SURFACE  
FOUND CONC

48.39




Figure II-34. TNB on Metal - Graph of Target-Found Concentration Points

Table II-37. TNB on Metal - Inaccuracy and Imprecision Data

1,3,5-TRINITROBENZENE (135TNB)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 2.293                        | 0.139                 | -8.300                 | 6.069       |
| 5.000                        | 4.428                        | 0.515                 | -11.450                | 11.639      |
| 10.000                       | 9.330                        | 0.858                 | -6.700                 | 9.194       |
| 25.000                       | 24.233                       | 0.960                 | -3.070                 | 3.962       |
| 50.000                       | 47.280                       | 0.914                 | -5.440                 | 1.934       |
| Means                        |                              | 0.677                 | -6.992                 | 6.560       |

1,3,5-TRINITROBENZENE (135TNB)  
METAL SURFACE

MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-3.07




Figure 11-35. TNB on Metal - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)  
 METAL SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure 11-36. TNB on Metal - Graph of Imprecision



Table II-38. TNB on Concrete - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)  
CONCRETE SURFACE

| TARGET CONC. US FOUND CONC |                                      |
|----------------------------|--------------------------------------|
| Target Conc                | Found Conc                           |
| ug/10 sq cm                | ug/10 sq cm                          |
| 2.500                      | 1.950<br>1.920<br>2.090<br>1.800     |
| 5.000                      | 3.160<br>3.510<br>5.120<br>3.690     |
| 10.000                     | 5.880<br>7.170<br>8.570<br>7.620     |
| 25.000                     | 20.890<br>22.660<br>20.540<br>17.150 |
| 50.000                     | 38.570<br>42.370<br>34.610<br>34.000 |

Table II-39. TNB on Concrete - Analysis of Target-  
Found Concentration Points

1,3,5-TRINITROBENZENE (135TNB)  
CONCRETE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535072

FOUND CONC  
MEAN= 14.1635 SD= 13.7366717317

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.222028721683  
SLOPE= 0.753593042071  
USE FOR ACCURACY  
R= 0.989402029999  
MEAN SQR DEV OF POINTS FROM REGRESSION= 4.19942063895  
ST ERROR EST= 2.04924879869  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18

THO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
y(c)= 3.95809784769  
x(d)= 9.81755206131

1,3,5-TRINITROBENZENE (135TNB)  
CONCRETE SURFACE  
FOUND CONC




Figure 11-37. TNB on Concrete - Graph of Target-Found Concentration Points

Table II-40. TNB on Concrete - Inaccuracy and Imprecision Data

| 1,3,5-TRINITROBENZENE (135TNB)             |                           |                    |                     |             |  |
|--------------------------------------------|---------------------------|--------------------|---------------------|-------------|--|
| CONCRETE SURFACE                           |                           |                    |                     |             |  |
| STATISTICAL DATA USED TO DETERMINE PERCENT |                           |                    |                     |             |  |
| INACCURACY AND IMPRECISION                 |                           |                    |                     |             |  |
| Mn Target Con ug/10 sq cm                  | Mn Found Conc ug/10 sq cm | Standard Deviation | Mean Pct Inaccuracy | Imprecision |  |
| 2.500                                      | 1.940                     | 0.119              | -22.400             | 6.142       |  |
| 5.000                                      | 3.870                     | 0.062              | -22.600             | 22.271      |  |
| 10.000                                     | 7.310                     | 1.118              | -26.900             | 15.291      |  |
| 25.000                                     | 20.310                    | 2.302              | -18.760             | 11.334      |  |
| 50.000                                     | 37.388                    | 3.891              | -25.225             | 10.406      |  |
| Means                                      |                           | 1.658              | -23.177             | 13.089      |  |

# 1,3,5-TRINITROBENZENE (135TNB)

CONCRETE SURFACE

MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-38. TNB on Concrete - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-39. TNB on Concrete - Graph of Imprecision

Table II-41. TNB on Brick - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

BRICK SURFACE

TARGET CONC. VS. FOUND CONC.

| Target Conc<br>ug/10 sq cm | Found Conc<br>ug/10 sq cm            |
|----------------------------|--------------------------------------|
| 2.500                      | 1.390<br>2.180<br>1.270<br>2.340     |
| 5.000                      | 2.800<br>6.180<br>4.090<br>4.950     |
| 10.000                     | 8.250<br>6.340<br>7.600<br>9.480     |
| 25.000                     | 22.460<br>14.690<br>12.960<br>21.150 |
| 50.000                     | 39.220<br>43.160<br>23.840<br>43.340 |

Table II-42. TNB on Brick - Analysis of Target-Found  
(Concentration Points)

1,3,5-TRINITROBENZENE (135TNB)

BRICK SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.83505335872

FOUND CONC  
MEAN= 13.8845 SD= 13.9256762002

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.272975323625  
SLOPE= 0.735758090615

USE FOR ACCURACY

R= 0.952875566024

MEAN SQR DEV OF POINTS FROM REGRESSION= 18.8379829636

ST ERROR EST= 4.3402745263

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 8.18590696298

x(d)= 21.2428137955



A graph illustrating the relationship between variables  $y(c)$  and  $x(d)$ . The vertical axis ( $y(c)$ ) ranges from 0 to 43.34, and the horizontal axis ( $x(d)$ ) ranges from 0 to 50. Three parallel lines are plotted, representing different target concentrations. A specific point  $(c, d)$  is marked on the middle line. The vertical axis TIC interval is noted as 4.334.

Figure II-40. TNB on Brick - Graph of Target-Found Concentration Points

Table II-43. TNB on Brick - Inaccuracy and Imprecision  
Data

**1,3,5-TRINITROBENZENE (135TNB)**

**BRICK SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION**

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.795                        | 0.543                 | -28.200                | 30.257      |
| 5.000                        | 4.505                        | 1.424                 | -9.900                 | 31.608      |
| 10.000                       | 7.918                        | 1.309                 | -20.825                | 16.534      |
| 25.000                       | 17.815                       | 4.692                 | -28.740                | 26.335      |
| 50.000                       | 37.390                       | 9.231                 | -25.220                | 24.689      |
| Means                        |                              | 3.440                 | -22.577                | 25.885      |

1,3,5-TRINITROBENZENE (135TNB)  
 BRICK SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-41. TNB on Brick - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)  
 BRICK SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-42. TNB on Brick - Graph of Imprecision

Table II-44. TNB on Transite - Target vs. Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

| TRANSITE SURFACE           |                                      |
|----------------------------|--------------------------------------|
| TARGET CONC. VS FOUND CONC |                                      |
| Target Conc<br>ug/10 sq cm | Found Conc<br>ug/10 sq cm            |
| 2.500                      | 0.860<br>2.050<br>1.780<br>1.540     |
| 5.000                      | 1.600<br>3.760<br>4.620<br>3.730     |
| 10.000                     | 3.270<br>8.230<br>8.220<br>8.140     |
| 25.000                     | 6.990<br>17.090<br>20.220<br>18.330  |
| 50.000                     | 12.180<br>33.330<br>34.930<br>39.400 |

Table II-45. TNB on Transite - Analysis of Target-  
Found Concentration Points

1,3,5-TRINITROBENZENE (135TNB)  
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0350535072

FOUND CONC

MEAN= 11.5135 SD= 12.010847849

NB. RUNS 1 TOTAL X-Y ALL RUNS 20 NB. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.355752427104

SLOPE= 0.592310679612

USE FOR ACCURACY

R= 0.889392237864

MEAN SQR DEV OF POINTS FROM REGRESSION= 31.0226365031

ST ERROR EST= 5.64115560706

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 10.840373461

x(d)= 34.6225232479



Table II-46. TNB on Transite - Inaccuracy and Imprecision Data

1,3,5-TRINITROBENZENE (135TNB)

TRANSITE SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.558                        | 0.510                 | -37.700                | 32.715      |
| 5.000                        | 3.428                        | 1.286                 | -31.450                | 37.529      |
| 10.000                       | 6.965                        | 2.464                 | -30.350                | 35.372      |
| 25.000                       | 15.658                       | 5.920                 | -37.370                | 37.889      |
| 50.000                       | 29.960                       | 12.128                | -40.000                | 40.482      |
| Means                        |                              | 4.462                 | -35.390                | 36.781      |



1,3,5-TRINITROBENZENE (135TNB)  
 TRANSITE SURFACE  
 MEAN FOUND CONCENTRATION (REPORT)  
 MEAN INACCURACY




Figure II-44. TNB on Transite - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)  
 TRANSITE SURFACE  
 IMPRECISION NEAR FOUND CONCENTRATION (REPORT)




Figure II-45. TNB on Transite - Graph of Imprecision

Table II-47. TNB on Transite (3 days) - Target vs.  
Found Concentrations

1,3,5-TRINITROBENZENE (135TNB)

| TRANSITE SURFACE<br>TARGET CONC. VS. FOUND CONC. |                            |
|--------------------------------------------------|----------------------------|
| Target Conc<br>ug/10 sq cm                       | Found Conc<br>ug/10 sq cm  |
| 2.500                                            | 2.050<br>1.780<br>1.540    |
| 5.000                                            | 3.760<br>4.620<br>3.730    |
| 10.000                                           | 8.230<br>8.220<br>8.140    |
| 25.000                                           | 17.090<br>20.220<br>18.330 |
| 50.000                                           | 33.330<br>34.930<br>39.400 |

Table II-48. TNB on Transite (3 days) - Analysis of Target-Found Concentration Points

1,3,5-TRINITROBENZENE (135TNB)  
TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.1953683274

FOUND CONC

MEAN= 13.691333333 SD= 13.0087014102

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.541645091694

SLOPE= 0.710793959008

USE FOR ACCURACY

R= 0.994192846865

MEAN SQR DEV OF POINTS FROM REGRESSION= 2.11040858187

ST ERROR EST= 1.45275207171

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 3.28917417064

X(D)= 7.64791513283

1,3,5-TRINITROBENZENE (135TNB)  
TRANSITE SURFACE  
FOUND CONC





Figure II-46. TNB on Transite (3 days) - Graph of Target-Found Concentration Points

Table II-49. TNB on Transite (3 days) - Inaccuracy and Imprecision Data

| 1,3,5-TRINITROBENZENE (135TNB)             |                              |                       |                        |             |  |
|--------------------------------------------|------------------------------|-----------------------|------------------------|-------------|--|
| TRANSITE SURFACE                           |                              |                       |                        |             |  |
| STATISTICAL DATA USED TO DETERMINE PERCENT |                              |                       |                        |             |  |
| INACCURACY AND IMPRECISION                 |                              |                       |                        |             |  |
| Mn Target Con<br>ug/10 sq cm               | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |  |
| 2.500                                      | 1.790                        | 0.255                 | -28.400                | 14.254      |  |
| 5.000                                      | 4.037                        | 0.505                 | -19.267                | 12.520      |  |
| 10.000                                     | 8.197                        | 0.049                 | -18.033                | 0.602       |  |
| 25.000                                     | 18.547                       | 1.576                 | -25.813                | 8.499       |  |
| 50.000                                     | 35.887                       | 3.146                 | -28.227                | 8.767       |  |
| Means                                      |                              | 1.106                 | -23.948                | 8.928       |  |

1,3,5-TRINITROBENZENE (1351NB)  
 TRANSITE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION <REPORT>



VERTICAL AXIS TIC INTERVAL = 1.03666666667

Figure II-47. TNB on Transite (3 days) - Graph of Inaccuracy

1,3,5-TRINITROBENZENE (135TNB)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-48. TNB on Transite - Graph of Imprecision



Table II-50. 2,4-DNT on Metal - Target vs. Found Concentrations

| 2,4-DINITROTOLUENE (24DNT) |                                            |
|----------------------------|--------------------------------------------|
| METAL SURFACE              |                                            |
| TARGET CONC. us/10 sq cm   | US FOUND CONC<br>Found Conc<br>us/10 sq cm |
| 2.500                      | 2.120                                      |
|                            | 2.300                                      |
|                            | 2.480                                      |
|                            | 2.070                                      |
| 5.000                      | 3.710                                      |
|                            | 4.370                                      |
|                            | 4.170                                      |
|                            | 4.840                                      |
| 10.000                     | 9.020                                      |
|                            | 8.720                                      |
|                            | 10.040                                     |
|                            | 11.180                                     |
| 25.000                     | 21.000                                     |
|                            | 20.820                                     |
|                            | 25.960                                     |
|                            | 21.740                                     |
| 50.000                     | 39.160                                     |
|                            | 38.400                                     |
|                            | 45.230                                     |
|                            | 46.510                                     |

Table II-51. 2,4-DNT on Metal - Analysis of Target-  
Found Concentration Points

2,4-DINITROTOLUENE (24DNT)  
METAL SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 16.192 SD= 15.3408483262

NB. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.59741302589  
SLOPE= 0.842950647249  
USE FOR ACCURACY  
R= 0.990992138847  
MEAN SQR DEV OF POINTS FROM REGRESSION= 4.45523971932  
ST ERROR EST= 2.11974387819  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18

TWO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
Y(C)= 4.44559634128  
X(d)= 9.04498923486

2,4-DINITROTOLUENE (24DNT)  
METAL SURFACE  
FOUND CONC

46.51




Figure II-49. 2,4-DNT on Metal - Graph of Target-Found Concentration Points

Table II-52. 2,4-DNT on Metal - Inaccuracy and Imprecision Data

2,4-DINITROTOLUENE (24DNT)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 2.243                        | 0.187                 | -10.300                | 8.322       |
| 5.000                        | 4.273                        | 0.468                 | -14.530                | 10.965      |
| 10.000                       | 9.740                        | 1.114                 | -2.600                 | 11.436      |
| 25.000                       | 22.380                       | 2.420                 | -10.480                | 10.812      |
| 50.000                       | 42.325                       | 4.138                 | -15.350                | 9.777       |
| Means                        |                              | 1.665                 | -10.656                | 10.263      |

2,4-DINITROTOLUENE (24DNT)  
 METAL SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-50. 2,4-DNT on Metal - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)  
METAL SURFACE

IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure II-51. 2,4-DNT on Metal - Graph of Imprecision

Table II-53. 2,4-DNT on Concrete - Target vs. Found Concentrations

| 2,4-DINITROTOLUENE (24DNT)<br>CONCRETE SURFACE |                           |
|------------------------------------------------|---------------------------|
| Target Conc<br>ug/10 sq cm                     | Found Conc<br>ug/10 sq cm |
| 2.500                                          | 2.530                     |
|                                                | 2.140                     |
|                                                | 2.100                     |
|                                                | 2.090                     |
| 5.000                                          | 4.340                     |
|                                                | 4.510                     |
|                                                | 4.140                     |
|                                                | 3.830                     |
| 10.000                                         | 6.310                     |
|                                                | 8.360                     |
|                                                | 7.890                     |
|                                                | 8.160                     |
| 25.000                                         | 23.120                    |
|                                                | 20.680                    |
|                                                | 22.670                    |
|                                                | 17.900                    |
| 50.000                                         | 33.880                    |
|                                                | 41.990                    |
|                                                | 44.240                    |
|                                                | 36.070                    |

Table II-54. 2,4-DNT on Concrete - Analysis of  
Target Found Concentration Points

2,4-DINITROTOLUENE (24DNT)  
CONCRETE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 14.8475 SD= 14.301442281

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.364439724919  
SLOPE= 0.782868122977  
USE FOR ACCURACY  
R= 0.987247878379  
MEAN SQR DEV OF POINTS FROM REGRESSION= 5.47110764003  
ST ERROR EST= 2.33903989706  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18  
TWO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
Y(C)= 4.62883883032  
X(D)= 10.7804175652





Table II-55. 2,4-DNT on Concrete - Inaccuracy and Imprecision Data

| 2,4-DINITROTOLUENE (24DNT)<br>CONCRETE SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                              |                       |                        |             |
|----------------------------------------------------------------------------------------------------------------------------|------------------------------|-----------------------|------------------------|-------------|
| Mn Target Con<br>ug/10 sq cm                                                                                               | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
| 2.500                                                                                                                      | 2.215                        | 0.211                 | -11.400                | 9.531       |
| 5.000                                                                                                                      | 4.205                        | 0.292                 | -15.900                | 6.948       |
| 10.000                                                                                                                     | 7.680                        | 0.933                 | -23.200                | 12.154      |
| 25.000                                                                                                                     | 21.093                       | 2.378                 | -15.630                | 11.273      |
| 50.000                                                                                                                     | 39.045                       | 4.871                 | -21.910                | 12.476      |
| Means                                                                                                                      |                              | 1.737                 | -17.608                | 10.476      |

2,4-DINITROTOLUENE (24DNT)  
 CONCRETE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-53. 2,4-DNT on Concrete - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure II-54. 2,4-DNT on Concrete - Graph of Imprecision

Table II-56. 2,4-DNT on Brick - Target vs. Found Concentrations

| 2,4-DINITROTOLUENE (24DNT)<br>BRICK SURFACE |                                             |
|---------------------------------------------|---------------------------------------------|
| TARGET CONC.<br>ug/10 sq cm                 | US FOUND CONC.<br>Found Conc<br>ug/10 sq cm |
| 2.500                                       | 1.530<br>2.270<br>1.040<br>1.780            |
| 5.000                                       | 2.410<br>5.420<br>4.160<br>4.740            |
| 10.000                                      | 7.570<br>6.300<br>7.420<br>9.830            |
| 25.000                                      | 22.860<br>14.260<br>13.890<br>20.880        |
| 50.000                                      | 37.560<br>37.440<br>22.770<br>42.410        |

Table II-57. 2,4-DNT on Brick - Analysis of Target-  
Found Concentration Points

2,4-DINITROTOLUENE (24DNT)  
BRICK SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 13.327 SD= 13.1141246071

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.494765372168  
SLOPE= 0.693634304207  
USE FOR ACCURACY  
R= 0.953912839867  
MEAN SQR DEV OF POINTS FROM REGRESSION= 16.347255196  
ST ERROR EST= 4.04317390128  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18  
TWO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
Y(C)= 7.86604073574  
X(d)= 20.9893901903

2,4-DINITROTOLUENE (24DNT)  
BRICK SURFACE  
FOUND CONC

42.41




Figure II-55. 2,4-DNT on Brick - Graph of Target-  
Found Concentration Points

Table II-58. 2,4-DNT on Brick - Inaccuracy and Imprecision Data

| 2,4-DINITROTOLUENE (24DNT)<br>BRICK SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                              |                       |                        |             |
|-------------------------------------------------------------------------------------------------------------------------|------------------------------|-----------------------|------------------------|-------------|
| Mn Target Conc<br>ug/10 sq cm                                                                                           | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
| 2.500                                                                                                                   | 1.655                        | 0.512                 | -33.800                | 30.961      |
| 5.000                                                                                                                   | 4.183                        | 1.289                 | -16.350                | 30.819      |
| 10.000                                                                                                                  | 7.780                        | 1.479                 | -22.200                | 19.016      |
| 25.000                                                                                                                  | 17.973                       | 4.575                 | -28.110                | 25.455      |
| 50.000                                                                                                                  | 35.045                       | 8.505                 | -29.910                | 24.267      |
| Means                                                                                                                   |                              | 3.272                 | -26.074                | 26.104      |



2,4-DINITROTOLUENE (24DNT)  
 BRICK SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-56. 2,4-DNT on Brick - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)

BRICK SURFACE

IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-57. 2,4-DNT on Brick - Graph of Imprecision

Table II-59. 2,4-DNT on Transite - Target vs. Found Concentrations

| 2,4-DINITROTOLUENE (24DNT) |                                       |
|----------------------------|---------------------------------------|
| TRANSITE SURFACE           |                                       |
| TARGET CONC. ug/10 sq cm   | US FOUND CONC. Found Conc ug/10 sq cm |
| 2.500                      | 0.700                                 |
|                            | 2.220                                 |
|                            | 2.040                                 |
|                            | 1.660                                 |
| 5.000                      | 1.540                                 |
|                            | 3.940                                 |
|                            | 4.250                                 |
|                            | 3.050                                 |
| 10.000                     | 2.850                                 |
|                            | 7.620                                 |
|                            | 7.190                                 |
|                            | 8.020                                 |
| 25.000                     | 7.950                                 |
|                            | 16.380                                |
|                            | 18.930                                |
|                            | 21.220                                |
| 50.000                     | 11.620                                |
|                            | 37.580                                |
|                            | 39.370                                |
|                            | 41.520                                |

Table II-60. 2,4-DNT on Transite - Analysis of Target-  
Found Concentration Points

2,4-DINITROTOLUENE (24DNT)  
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 11.9825 SD= 13.2377478109

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.0521541262134  
SLOPE= 0.65052184466  
USE FOR ACCURACY  
R= 0.886268306037  
MEAN SQR DEV OF POINTS FROM REGRESSION= 39.6820662231  
ST ERROR EST= 6.29937030385  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73486096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 11.4324849867

x(d)= 35.2284612374




Table II-61. 2,4-DNT on Transite - Inaccuracy and Imprecision Data

| 2,4-DINITROTOLUENE (24DNT)<br>TRANSITE SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |               |                    |                     |             |
|----------------------------------------------------------------------------------------------------------------------------|---------------|--------------------|---------------------|-------------|
| Mn Target Con                                                                                                              | Mn Found Conc | Standard Deviation | Mean Pct Inaccuracy | Imprecision |
| 2.500                                                                                                                      | 1.655         | 0.678              | -33.800             | 48.973      |
| 5.000                                                                                                                      | 3.195         | 1.215              | -36.100             | 38.026      |
| 10.000                                                                                                                     | 6.400         | 2.404              | -35.800             | 37.446      |
| 25.000                                                                                                                     | 16.120        | 5.794              | -35.520             | 35.945      |
| 50.000                                                                                                                     | 32.523        | 14.028             | -34.955             | 43.133      |
| Means                                                                                                                      |               | 4.824              | -35.235             | 39.104      |

2,4-DINITROTOLUENE (24DNT)  
 TRANSITE SURFACE  
 MEAN INACCURACY

MEAN FOUND CONCENTRATION (REPORT)

-33.800

0

-36.10

VERTICAL AXIS TIC INTERVAL = 0.23 TARGET CONC

8

9

10

20

30

40

50

Figure II-59. 2,4-DNT on Transite - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-60. 2,4-DNT on Transite - Graph of Imprecision



Table II-62. 2,4-DNT on Transite (3 days) - Target vs.  
Found Concentrations

| 2,4-DINITROTOLUENE (24DNT)<br>TRANSITE SURFACE<br>TARGET CONC. VS. FOUND CONC. |                            |
|--------------------------------------------------------------------------------|----------------------------|
| Target Conc<br>ug/10 sq cm                                                     | Found Conc<br>ug/10 sq cm  |
| 2.500                                                                          | 2.220<br>2.040<br>1.660    |
| 5.000                                                                          | 3.940<br>4.250<br>3.050    |
| 10.000                                                                         | 7.620<br>7.190<br>8.020    |
| 25.000                                                                         | 16.380<br>18.930<br>21.220 |
| 50.000                                                                         | 37.580<br>39.370<br>41.520 |

Table 11-63. 2,4-DNT on Transite (3 days) - Analysis  
of Target-Found Concentration Points

2,4-DINITROTOLUENE (24DNT)  
TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.1953683274

FOUND CONC  
MEAN= 14.332666667 SD= 14.4185805061

NØ. RUNS 1 TOTAL X-Y ALL RUNS 15 NØ. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.270698489752

SLOPE= 0.789371089536

USE FOR ACCURACY

R= 0.996138122967

MEAN SQR DEV OF POINTS FROM REGRESSION= 1.72591231227

ST ERROR EST= 1.31373981909

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77893170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 2.21392254463

X(D)= 6.23754930521

2,4-DINITROTOLUENE (24DNT)  
TRANSITE SURFACE  
FOUND CONC




Figure 11-61. 2,4-DNT on Transite (3 days) - Graph of Target-found Concentration Points.

Table 11-64. 2,4-DNT on Transite (3 days) - Inaccuracy and Imprecision Data

2,4-DINITROTOLUENE (24DNT)  
TRANSITE SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Target Con-<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                         | 1.973                        | 0.286                 | -21.067                | 14.488      |
| 5.000                         | 3.747                        | 0.623                 | -25.067                | 16.626      |
| 10.000                        | 7.610                        | 0.415                 | -23.900                | 5.455       |
| 25.000                        | 18.843                       | 2.421                 | -24.627                | 12.849      |
| 50.000                        | 39.490                       | 1.973                 | -21.020                | 4.996       |
| Means                         |                              | 1.144                 | -23.136                | 10.883      |

2,4-DINITROTOLUENE (24DNT)  
 TRANSITE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-62. 2,4-DNT on Transite (3 days) - Graph of Inaccuracy

2,4-DINITROTOLUENE (24DNT)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-63. 2,4-DNT on Transite (3 days) - Graph of Imprecision

Table 11-65. TNT on Metal - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

| METAL SURFACE               |                                             |
|-----------------------------|---------------------------------------------|
| TARGET CONC.<br>ug/10 sq cm | US FOUND CONC.<br>Found Conc<br>ug/10 sq cm |
| 2.500                       | 1.920                                       |
|                             | 2.140                                       |
|                             | 2.120                                       |
| 5.000                       | 1.810                                       |
|                             | 3.360                                       |
|                             | 4.510                                       |
| 10.000                      | 4.020                                       |
|                             | 3.910                                       |
| 25.000                      | 9.560                                       |
|                             | 9.090                                       |
|                             | 8.070                                       |
| 50.000                      | 8.680                                       |
|                             | 23.630                                      |
|                             | 24.690                                      |
|                             | 21.180                                      |
|                             | 20.570                                      |
|                             | 48.040                                      |
|                             | 50.620                                      |
|                             | 43.920                                      |
|                             | 50.080                                      |

2,4,6-TRINITROTOLUENE (246TNT)  
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 17.096 SD= 17.6143310075

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.90170105178  
SLOPE= 0.972848705502  
USE FOR ACCURACY  
R= 0.996085433419  
MEAN SQR DEV OF POINTS FROM REGRESSION= 2.55903491889  
ST ERROR EST= 1.59969838373  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F. = 18

TWO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE  
MEASURED 1 TIME(S))  
Y(C)= 2.01477436407  
X(D)= 5.95490116603



2,4,6-TRINITROTOLUENE (246TNT)  
METAL SURFACE  
FOUND CONC




Figure 11-64. TNT on Metal - Graph of Target-Found Concentration Points

Table 11-67. TNT on Metal - Inaccuracy and Imprecision Data

2,4,6-TRINITROTOLUENE (246TNT)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.998                        | 0.160                 | -20.100                | 7.993       |
| 5.000                        | 3.950                        | 0.472                 | -21.000                | 11.948      |
| 10.000                       | 8.850                        | 0.632                 | -11.500                | 7.143       |
| 25.000                       | 22.518                       | 1.961                 | -9.930                 | 8.710       |
| 50.000                       | 48.165                       | 3.040                 | -3.670                 | 6.312       |
| Means                        |                              | 1.253                 | -13.240                | 8.421       |

2,4,6-TRINITROTOLUENE (246TNT)  
 METAL SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-65. TNT on Metal - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)  
 METAL SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-66. TNT on Metal - Graph of Imprecision

Table II-68. TNT on Concrete - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

| CONCRETE SURFACE            |                              |
|-----------------------------|------------------------------|
| TARGET CONC.<br>ug/10 sq cm | VS FOUND CONC<br>ug/10 sq cm |
| 2.500                       | 1.540                        |
|                             | 2.020                        |
|                             | 1.310                        |
|                             | 1.790                        |
| 5.000                       | 3.950                        |
|                             | 4.430                        |
|                             | 3.090                        |
|                             | 3.200                        |
| 10.000                      | 6.960                        |
|                             | 7.630                        |
|                             | 5.250                        |
|                             | 6.690                        |
| 25.000                      | 19.840                       |
|                             | 20.130                       |
|                             | 14.290                       |
|                             | 17.430                       |
| 50.000                      | 40.730                       |
|                             | 43.930                       |
|                             | 27.160                       |
|                             | 35.640                       |

Table II-69. TNT on Concrete Analysis of Target-  
Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)  
CONCRETE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 13.3505 SD= 13.7346130078

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.36244538835

SLOPE= 0.741240291262

USE FOR ACCURACY

R= 0.973329817619

MEAN SQR DEV OF POINTS FROM REGRESSION= 10.4794771204

ST ERROR EST= 3.23720205122

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 5.53942979924

x(d)= 15.7269052076

2,4,6-TRINITROTOLUENE (246TNT)  
 CONCRETE SURFACE  
 FOUND CONC




Figure 11-67. TNT on Concrete - Graph of Target-Found Concentration Points

Table II-70. TNT on Concrete - Inaccuracy and Imprecision Data

2,4,6-TRINITROTOLUENE (246TNT)

CONCRETE SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.665                        | 0.307                 | -33.400                | 18.456      |
| 5.000                        | 3.668                        | 0.636                 | -26.650                | 17.340      |
| 10.000                       | 6.633                        | 1.003                 | -33.675                | 15.120      |
| 25.000                       | 17.923                       | 2.707                 | -28.310                | 15.105      |
| 50.000                       | 36.865                       | 7.315                 | -26.270                | 19.843      |
| Means                        |                              | 2.394                 | -29.661                | 17.173      |



2,4,6-TRINITROTOLUENE (246TNT)  
 CONCRETE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)  
 -26.27




Figure 11-68. TNT on Concrete - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-69. TNT on Concrete - Graph of Imprecision

Table 11-71. TNT on Brick - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

BRICK SURFACE

TARGET CONC. VS. FOUND CONC.

| Target Conc<br>ug/10 sq cm | Found Conc<br>ug/10 sq cm |
|----------------------------|---------------------------|
| 2.500                      | 1.470                     |
|                            | 1.960                     |
|                            | 0.890                     |
|                            | 1.800                     |
| 5.000                      | 2.320                     |
|                            | 5.930                     |
|                            | 4.100                     |
|                            | 4.180                     |
| 10.000                     | 8.080                     |
|                            | 6.280                     |
|                            | 7.310                     |
|                            | 9.400                     |
| 25.000                     | 20.760                    |
|                            | 15.210                    |
|                            | 13.750                    |
|                            | 20.760                    |
| 50.000                     | 36.650                    |
|                            | 34.530                    |
|                            | 24.570                    |
|                            | 38.460                    |

Table 11-72. TNT on Brick - Analysis of Target-  
Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)  
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0350535872

FOUND CONC

MEAN= 12.9205 SD= 12.3402646372

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.645345873786

SLOPE= 0.66352184466

USE FOR ACCURACY

R= 0.969724100462

MEAN SQR DEV OF POINTS FROM REGRESSION= 9.50589122507

ST ERROR EST= 3.09610904606

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 6.28998864488

X(D)= 16.80061912189

[illegible]

### ENT on Brick - Graph Concentration Points

Table 11-73. TNT on Brick - Inaccuracy and Imprecision Data

**2,4,6-TRINITROTOLUENE (246TNT)**

**BRICK SURFACE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION**

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.530                        | 0.473                 | -38.800                | 38.911      |
| 5.000                        | 4.133                        | 1.474                 | -17.350                | 35.672      |
| 10.000                       | 7.768                        | 1.315                 | -22.325                | 16.925      |
| 25.000                       | 17.620                       | 3.674                 | -29.520                | 20.854      |
| 50.000                       | 33.553                       | 6.200                 | -32.895                | 18.478      |
| Means                        |                              | 2.627                 | -28.178                | 24.568      |

2,4,6-TRINITROTOLUENE (246TNT)  
 BRICK SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-71. TNT on Brick - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)  
 BRICK SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-72. TNT on Brick - Graph of Imprecision



Table U-74. TNT on Transite - Target vs. Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

| TRANSITE SURFACE         |                                      |
|--------------------------|--------------------------------------|
| TARGET CONC. ug/10 sq cm | US FOUND CONC Found Conc ug/10 sq cm |
| 2.500                    | 1.160                                |
|                          | 1.770                                |
|                          | 1.770                                |
|                          | 0.900                                |
| 5.000                    | 1.560                                |
|                          | 3.900                                |
|                          | 3.510                                |
|                          | 2.520                                |
| 10.000                   | 2.820                                |
|                          | 7.770                                |
|                          | 6.400                                |
|                          | 5.880                                |
| 25.000                   | 6.450                                |
|                          | 17.370                               |
|                          | 16.330                               |
|                          | 15.800                               |
| 50.000                   | 11.720                               |
|                          | 28.790                               |
|                          | 29.470                               |
|                          | 32.170                               |

Table 11-75. TNT on Transit - Analysis of Target-Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)  
TRANSITE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 9.903 SD= 10.1474591583

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.496468446602

SLOPE= 0.508461165049

USE FOR ACCURACY

R= 0.903686747162

MEAN SQR DEV OF POINTS FROM REGRESSION= 19.9286214266  
ST ERROR EST= 4.46414845482

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 8.63523970102

X(D)= 31.8194111464

[illegible]

Arthur D. Little, Inc.

Table 11-76. TNT on Transite - Inaccuracy and Imprecision Data


2,4,6-TRINITROTOLUENE (246TNT)  
 TRANSITE SURFACE  
 STATISTICAL DATA USED TO DETERMINE PERCENT  
 INACCURACY AND IMPRECISION

| Mn Targt Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-----------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                       | 1.400                        | 0.440                 | -44.000                | 31.445      |
| 5.000                       | 2.873                        | 1.050                 | -42.550                | 36.562      |
| 10.000                      | 5.718                        | 2.090                 | -42.825                | 36.549      |
| 25.000                      | 13.988                       | 5.067                 | -44.050                | 36.226      |
| 50.000                      | 25.538                       | 9.327                 | -48.925                | 36.521      |
| Means                       |                              | 3.595                 | -44.470                | 35.461      |

2,4,6-TRINITROTOLUENE (246TNT)  
TRANSITE SURFACE

MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-42.55



VERTICAL AXIS TIC INTERVAL = 0.6375

Figure 11-74. TNT on Transite - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-75. TNT on Transite - Graph of Imprecision

Table II-77. TNT on Transite (3 days) - Target vs.  
Found Concentrations

2,4,6-TRINITROTOLUENE (246TNT)

TRANSITE SURFACE

| TARGET CONC. VS. FOUND CONC. |                           |
|------------------------------|---------------------------|
| Target Conc<br>ug/10 sq cm   | Found Conc<br>ug/10 sq cm |
| 2.500                        | 1.770                     |
|                              | 1.770                     |
|                              | 0.900                     |
| 5.000                        | 3.900                     |
|                              | 3.510                     |
|                              | 2.520                     |
| 10.000                       | 7.770                     |
|                              | 6.400                     |
|                              | 5.880                     |
| 25.000                       | 17.370                    |
|                              | 16.330                    |
|                              | 15.800                    |
| 50.000                       | 28.790                    |
|                              | 29.470                    |
|                              | 32.170                    |

Table 11-78. TNT on Transite (3 days) - Analysis of  
Target-Found Concentration Points

2,4,6-TRINITROTOLUENE (246TNT)

TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.1953683274

FOUND CONC

MEAN= 11.623333333 SD= 11.0214567787

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.468352211435

SLOPE= 0.602971952535

USE FOR ACCURACY

R= 0.995448876481

MEAN SQR DEV OF POINTS FROM REGRESSION= 1.18801497594

ST ERROR EST= 1.08996099744

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 2.52974974353

X(D)= 6.77062537754



32.17




Figure 11-7b. TNT on Transit (3 days) - Graph of Target-found Concentration Points

Table 11-79. TNT on Transite (3 days) - Inaccuracy and Imprecision Data

2,4,6-TRINITROTOLUENE (246TNT)  
 TRANSITE SURFACE  
 STATISTICAL DATA USED TO DETERMINE PERCENT  
 INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 1.480                        | 0.502                 | -40.800                | 33.939      |
| 5.000                        | 3.310                        | 0.711                 | -33.800                | 21.493      |
| 10.000                       | 6.683                        | 0.976                 | -33.167                | 14.689      |
| 25.000                       | 16.500                       | 0.799                 | -34.000                | 4.841       |
| 50.000                       | 30.143                       | 1.788                 | -39.713                | 5.931       |
| Means                        |                              | 0.955                 | -36.296                | 16.162      |

2,4,6-TRINITROTOLUENE (246TNT)  
 TRANSITE SURFACE  
 MEAN INACCURACY




Figure 11-77. TNT on Transite (3 days) - Graph of Inaccuracy

2,4,6-TRINITROTOLUENE (246TNT)  
TRANSITE SURFACE  
IMPRECISION MEAN FOUND CONCENTRATION <REPORT>




Figure 11-78. TNT on Transite (3 days) - Graph of Imprecision

Table 11-80. Tetryl on Metal - Target vs. Found Concentrations

**2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)**

| <b>METAL SURFACE</b>            |                                              |
|---------------------------------|----------------------------------------------|
| <b>TARGET CONC. ug/10 sq cm</b> | <b>US FOUND CONC. Found Conc ug/10 sq cm</b> |
| <b>2.500</b>                    | <b>1.820</b>                                 |
|                                 | <b>0.810</b>                                 |
|                                 | <b>0.000</b>                                 |
|                                 | <b>0.000</b>                                 |
| <b>5.000</b>                    | <b>2.660</b>                                 |
|                                 | <b>2.430</b>                                 |
|                                 | <b>0.110</b>                                 |
|                                 | <b>0.250</b>                                 |
| <b>10.000</b>                   | <b>8.520</b>                                 |
|                                 | <b>8.140</b>                                 |
|                                 | <b>2.620</b>                                 |
|                                 | <b>2.070</b>                                 |
| <b>25.000</b>                   | <b>22.940</b>                                |
|                                 | <b>19.600</b>                                |
|                                 | <b>8.610</b>                                 |
|                                 | <b>11.160</b>                                |
| <b>50.000</b>                   | <b>48.900</b>                                |
|                                 | <b>46.740</b>                                |
|                                 | <b>41.180</b>                                |
|                                 | <b>32.350</b>                                |

Table 11-81. Tetryl on Metal - Analysis of Target-  
Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0330535872

FOUND CONC

MEAN= 13.0455 SD= 16.5283732046

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -3.22953074434

SLOPE= 0.879731391586

USE FOR ACCURACY

R= 0.959925250553

MEAN SQR DEV OF POINTS FROM REGRESSION= 22.6491360617

ST ERROR EST= 4.7591108478

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 5.44699848777

x(d)= 19.476685492

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
METAL SURFACE  
FOUND CONC




Figure 11-79. Tetryl on Metal - Graph of Target-Found Concentration Points

Table II-82. Tetryl on Metal - Inaccuracy and Imprecision Data

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

METAL SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT

INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 0.658                        | 0.864                 | -73.700                | 131.401     |
| 5.000                        | 1.363                        | 1.370                 | -72.750                | 100.340     |
| 10.000                       | 5.338                        | 3.466                 | -46.625                | 64.941      |
| 25.000                       | 15.578                       | 6.793                 | -37.690                | 43.610      |
| 50.000                       | 42.293                       | 7.383                 | -15.415                | 17.457      |
| Means                        |                              | 3.975                 | -49.236                | 71.590      |



2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 METAL SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-15.42




Figure 11-80. Tetryl on Metal - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 METAL SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-81. Tetryl on Metal - Graph of Imprecision

Table II-83. Tetryl on Concrete - Target vs. Found Concentrations

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

| TETRYL CONC. VS. FOUND CONC.<br>CONCRETE SURFACE |                           |
|--------------------------------------------------|---------------------------|
| Target Conc<br>ug/10 sq cm                       | Found Conc<br>ug/10 sq cm |
| 2.500                                            | 0.230                     |
|                                                  | 0.520                     |
|                                                  | 0.250                     |
|                                                  | 0.350                     |
| 5.000                                            | 1.500                     |
|                                                  | 1.790                     |
|                                                  | 0.350                     |
|                                                  | 1.100                     |
| 10.000                                           | 2.880                     |
|                                                  | 3.960                     |
|                                                  | 1.280                     |
|                                                  | 1.250                     |
| 25.000                                           | 15.930                    |
|                                                  | 14.650                    |
|                                                  | 3.460                     |
|                                                  | 8.500                     |
| 50.000                                           | 28.080                    |
|                                                  | 40.160                    |
|                                                  | 9.780                     |
|                                                  | 18.020                    |

Table 11-84. Tetryl on Concrete - Analysis of  
Target-Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
CONCRETE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 7.702 SD= 10.8150614277

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -1.68730380259  
SLOPE= 0.507529935275

USE FOR ACCURACY  
R= 0.846350216412

MEAN SQR DEV OF POINTS FROM REGRESSION= 35.0255619146  
ST ERROR EST= 5.91823976488

USE FOR PRECISION  
T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 9.10248066865

x(d)= 42.9096479565



Table 11-35. Tetryl on Concrete - Inaccuracy and Imprecision Data

**2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)**

**CONCRETE SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION**

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                        | 0.338                        | 0.133                 | -86.500                | 39.262      |
| 5.000                        | 1.185                        | 0.624                 | -76.300                | 52.694      |
| 10.000                       | 2.343                        | 1.320                 | -76.575                | 56.353      |
| 25.000                       | 10.635                       | 5.779                 | -57.460                | 54.341      |
| 50.000                       | 24.010                       | 13.112                | -51.980                | 54.610      |
| Means                        |                              | 4.194                 | -69.763                | 51.452      |

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 CONCRETE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-83. Tetryl on Concrete - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-84. Tetryl on Concrete - Graph of Imprecision



Table 11-86. Tetryl on Brick - Target vs. Found Concentrations

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

| TARGET CONC. US FOUND CONC |                                      |
|----------------------------|--------------------------------------|
| Target Conc<br>ug/10 sq cm | Found Conc<br>ug/10 sq cm            |
| 2.500                      | 2.210<br>1.280<br>0.620<br>0.060     |
| 5.000                      | 2.300<br>6.050<br>3.910<br>0.600     |
| 10.000                     | 9.340<br>6.520<br>8.060<br>2.560     |
| 25.000                     | 24.170<br>17.540<br>17.820<br>11.820 |
| 50.000                     | 41.510<br>41.820<br>30.920<br>20.310 |

Table 11-87. Tetryl on Brick - Analysis of Target-  
Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
BRICK SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 12.471 SD= 13.2601609006

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.212591019417

SLOPE= 0.685599514563

USE FOR ACCURACY

R= 0.931916945925

MEAN SQR DEV OF POINTS FROM REGRESSION= 24.441622711

ST ERROR EST= 4.94384695465

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 8.80073814918

x(d)= 26.0154833108

41.82




Figure 11-35. Tetryl on Brick - Graph of Target-Found Concentration Points

Table 11-88. Tetryl on Brick - Inaccuracy and Imprecision Data

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

BRICK SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Targt Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-----------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                       | 1.043                        | 0.924                 | -58.300                | 88.667      |
| 5.000                       | 3.215                        | 2.323                 | -35.700                | 72.270      |
| 10.000                      | 6.620                        | 2.942                 | -33.800                | 44.441      |
| 25.000                      | 17.838                       | 5.046                 | -28.650                | 28.291      |
| 50.000                      | 33.640                       | 10.230                | -32.720                | 30.409      |
| Means                       |                              | 4.293                 | -37.834                | 52.816      |

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 BRICK SURFACE  
 MEAN INACCURACY




Figure 11-86. Tetryl on Brick - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 BRICK SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-S7. Tetryl on Brick - Graph of Imprecision

Table 11-89. Tetryl on Transite - Target vs. Found Concentrations

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)

| TRANSITE SURFACE |                                      |
|------------------|--------------------------------------|
| TARGET CONC. US  | FOUND CONC                           |
| ug/10 sq cm      | ug/10 sq cm                          |
| 2.500            | 0.950<br>1.340<br>0.140<br>0.170     |
| 5.000            | 1.460<br>3.080<br>0.540<br>0.450     |
| 10.000           | 2.670<br>7.850<br>2.640<br>0.900     |
| 25.000           | 6.880<br>22.300<br>9.000<br>4.300    |
| 50.000           | 12.710<br>40.460<br>21.230<br>35.070 |

Table II-90. Tetryl on Transite - Analysis of  
Target-Found Concentration Points

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
TRANSITE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 8.711 SD= 11.9174652462

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -1.73395631068  
SLOPE= 0.56459223301

USE FOR ACCURACY

R= 0.854414170035

MEAN SQR DEV OF POINTS FROM REGRESSION= 40.4738695793  
ST ERROR EST= 6.36190769968

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 9.86469728103

X(D)= 41.354710617



2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 TRANSITE SURFACE  
 FOUND CONC




Figure 11-88. Tetryl on Transite - Graph of Target-Found Concentration Points

Table 11-91. Tetryl on Transite - Inaccuracy and Imprecision Data

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
TRANSITE SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Target Conc<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 2.500                         | 0.650                        | 0.593                 | -74.000                | 91.302      |
| 5.000                         | 1.383                        | 1.220                 | -72.350                | 88.262      |
| 10.000                        | 3.515                        | 3.006                 | -64.850                | 85.522      |
| 25.000                        | 10.640                       | 8.015                 | -57.440                | 75.329      |
| 50.000                        | 27.368                       | 12.692                | -45.265                | 46.376      |
| Means                         |                              | 5.105                 | -62.781                | 77.358      |

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 TRANSITE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-89. Tetryl on Transite - Graph of Inaccuracy

2,4,6-TRINITROPHENYLMETHYLNITRAMINE (TETRYL)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-90. Tetryl on Transite - Graph of Imprecision

Table II-92. DPA on Metal - Target vs. Found Concentrations

| DIPHENYLAMINE (DPA)               |                                      |
|-----------------------------------|--------------------------------------|
| METAL SURFACE                     |                                      |
| TARGET CONC. Target Conc<br>ng/ml | VS FOUND CONC<br>Found Conc<br>ng/ml |
| 2.500                             | 2.050                                |
|                                   | 2.610                                |
|                                   | 2.560                                |
|                                   | 2.410                                |
| 5.000                             | 3.770                                |
|                                   | 5.070                                |
|                                   | 4.980                                |
|                                   | 4.720                                |
| 10.000                            | 9.520                                |
|                                   | 10.140                               |
|                                   | 10.270                               |
|                                   | 9.910                                |
| 25.000                            | 25.240                               |
|                                   | 25.420                               |
|                                   | 24.400                               |
|                                   | 21.800                               |
| 50.000                            | 49.370                               |
|                                   | 47.940                               |
|                                   | 47.220                               |
|                                   | 43.510                               |

Table II-93. DPA on Metal - Analysis of Target-Found  
Concentration Points

DIPHENYLAMINE (DPA)  
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 18.5 SD= 18.0350535872

FOUND CONC

MEAN= 17.6455 SD= 17.002404621

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.253584142395

SLOPE= 0.940103559871

USE FOR ACCURACY

R= 0.997201187585

MEAN SQR DEV OF POINTS FROM REGRESSION= 1.70567937344

ST ERROR EST= 1.30601660534

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

y(c)= 2.63463632131

x(d)= 5.03549506469




Figure 11-91. DPA on Metal - Graph of Target-Found Concentration Points

Table II-94. DPA on Metal - Inaccuracy and Imprecision Data

| DIPHENYLAMINE (DPA)<br>METAL SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                        |                       |                        |             |
|------------------------------------------------------------------------------------------------------------------|------------------------|-----------------------|------------------------|-------------|
| Mn Target Con<br>ng/ml                                                                                           | Mn Found Conc<br>ng/ml | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
| 2.500                                                                                                            | 2.400                  | 0.253                 | -3.700                 | 10.510      |
| 5.000                                                                                                            | 4.635                  | 0.595                 | -7.300                 | 12.847      |
| 10.000                                                                                                           | 9.950                  | 0.329                 | -0.400                 | 3.303       |
| 25.000                                                                                                           | 24.215                 | 1.670                 | -3.140                 | 6.898       |
| 50.000                                                                                                           | 47.010                 | 2.499                 | -5.980                 | 5.315       |
| Means                                                                                                            |                        | 1.069                 | -4.104                 | 7.774       |






Figure II-92. DPA on Metal - Graph of Inaccuracy

DIPHENYLAMINE (DPA)  
METAL SURFACE  
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-93. DPA on Metal - Graph of Imprecision

Table 11-95. DPA on Concrete - Target vs. Found Concentrations

| DIPHENYLAMINE (DPA)<br>CONCRETE SURFACE<br>TARGET CONC. VS. FOUND CONC. |                     |
|-------------------------------------------------------------------------|---------------------|
| Target Conc<br>ng/ml                                                    | Found Conc<br>ng/ml |
| 2.500                                                                   | 2.010               |
|                                                                         | 2.200               |
|                                                                         | 1.720               |
|                                                                         | 1.560               |
| 5.000                                                                   | 4.940               |
|                                                                         | 4.130               |
|                                                                         | 4.140               |
|                                                                         | 2.950               |
| 10.000                                                                  | 5.300               |
|                                                                         | 8.470               |
|                                                                         | 7.470               |
|                                                                         | 4.880               |
| 25.000                                                                  | 19.400              |
|                                                                         | 21.100              |
|                                                                         | 19.580              |
|                                                                         | 9.750               |
| 50.000                                                                  | 35.410              |
|                                                                         | 43.400              |
|                                                                         | 47.480              |
|                                                                         | 22.730              |

Table 11-96. DPA on Concrete - Analysis of Target-  
Found Concentration Points

DIPHENYLAMINE (DPA)  
CONCRETE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535072

FOUND CONC  
MEAN= 13.431 SD= 14.2563485027

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.303319174757  
SLOPE= 0.742395631068  
USE FOR ACCURACY  
R= 0.939170713081  
MEAN SQR DEV OF POINTS FROM REGRESSION= 25.3061729268  
ST ERROR EST= 5.03052412049  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18  
TWO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE  
MEASURED 1 TIME(S))  
Y(C)= 8.86803467264  
X(D)= 24.4266262877

DIPHENYLAMINE (DPA)  
CONCRETE SURFACE  
FOUND CONC




Figure 11-94. DPA on Concrete - Graph of Target-  
Found Concentration Points

Table II-97. DPA on Concrete - Inaccuracy and Imprecision Data

| DIPHENYLAMINE (DPA)<br>CONCRETE SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                        |                       |                        |             |
|---------------------------------------------------------------------------------------------------------------------|------------------------|-----------------------|------------------------|-------------|
| Mn Target Conc<br>ng/ml                                                                                             | Mn Found Conc<br>ng/ml | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
| 2.500                                                                                                               | 1.873                  | 0.287                 | -25.100                | 15.326      |
| 5.000                                                                                                               | 4.040                  | 0.820                 | -19.200                | 20.292      |
| 10.000                                                                                                              | 6.530                  | 1.721                 | -34.700                | 26.351      |
| 25.000                                                                                                              | 17.450                 | 5.195                 | -30.170                | 29.756      |
| 50.000                                                                                                              | 37.255                 | 10.904                | -25.490                | 29.269      |
| Means                                                                                                               |                        | 3.785                 | -26.932                | 24.199      |

DIPHENYLAMINE (DPA)  
 CONCRETE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-19.20

-34.70

VERTICAL AXIS TIC INTERVAL = 1.55 TARGET CONC




Figure 11-95. DPA on Concrete - Graph of Inaccuracy

DIPHENYLAMINE (DPA)  
CONCRETE SURFACE  
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-96. DPA on Concrete - Graph of Imprecision



Table II-98. BPA on Concrete (3 days) - Target vs. Found Concentrations

| DIPHENYLAMINE (DPA)<br>CONCRETE SURFACE |                     |
|-----------------------------------------|---------------------|
| TARGET CONC. VS FOUND CONC              |                     |
| Target Conc<br>ng/ml                    | Found Conc<br>ng/ml |
| 2.500                                   | 2.010               |
|                                         | 2.200               |
|                                         | 1.720               |
| 5.000                                   | 4.940               |
|                                         | 4.130               |
|                                         | 4.140               |
| 10.000                                  | 5.300               |
|                                         | 8.470               |
|                                         | 7.470               |
| 25.000                                  | 19.400              |
|                                         | 21.100              |
|                                         | 19.580              |
| 50.000                                  | 35.410              |
|                                         | 43.400              |
|                                         | 47.480              |

Table 11-99. BPA on Concrete (3 days) - Analysis of Target-Found Concentration Points

DIPHENYLAMINE (DPA)  
CONCRETE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.1953683274

FOUND CONC  
MEAN= 15.1166666667 SD= 15.5726285992

NO. RUNS 1 TOTAL X-Y ALL RUNS 15 NO. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.509147788565

SLOPE= 0.844638619202

USE FOR ACCURACY

R= 0.986893039028

MEAN SQR DEV OF POINTS FROM REGRESSION= 6.80118484773  
ST ERROR EST= 2.60790813637

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F. = 13

TWO TAIL P LEVEL IS .1

t= 1.77093170942

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 4.42307951945

X(D)= 11.514520374



Table 11-100. DPA on Concrete (3 days) - Inaccuracy and Imprecision Data

**DIPHENYLAMINE (DPA)**

**CONCRETE SURFACE**

**STATISTICAL DATA USED TO DETERMINE PERCENT**

**INACCURACY AND IMPRECISION**

| Mn Target Conc<br>ng/ml | Mn Found Conc<br>ng/ml | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-------------------------|------------------------|-----------------------|------------------------|-------------|
| 2.500                   | 1.977                  | 0.242                 | -20.933                | 12.229      |
| 5.000                   | 4.403                  | 0.465                 | -11.933                | 10.556      |
| 10.000                  | 7.080                  | 1.621                 | -29.200                | 22.890      |
| 25.000                  | 20.027                 | 0.934                 | -19.893                | 4.663       |
| 50.000                  | 42.097                 | 6.140                 | -15.807                | 14.585      |
| Means                   |                        | 1.880                 | -19.553                | 12.984      |

DIPHENYLAMINE (DPA)  
 CONCRETE SURFACE  
 MEAN INACCURACY  
 -11.93  
 MEAN FOUND CONCENTRATION (REPORT)




Figure 11-98. DPA on Concrete (3 days) - Graph of Inaccuracy

DIPHENYLAMINE (DPA)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-99. BPA on Concrete (3 days) - Graph of Imprecision

Table II-101. DPA on Brick - Target vs. Found Concentrations

| DIPHENYLAMINE (DPA)<br>BRICK SURFACE<br>TARGET CONC. VS. FOUND CONC. |                     |
|----------------------------------------------------------------------|---------------------|
| Target Conc<br>ng/ml                                                 | Found Conc<br>ng/ml |
| 2.500                                                                | 1.920               |
|                                                                      | 2.180               |
|                                                                      | 1.710               |
|                                                                      | 2.200               |
| 5.000                                                                | 2.260               |
|                                                                      | 4.950               |
|                                                                      | 4.550               |
|                                                                      | 4.750               |
| 10.000                                                               | 7.900               |
|                                                                      | 7.060               |
|                                                                      | 7.890               |
|                                                                      | 8.260               |
| 25.000                                                               | 22.880              |
|                                                                      | 15.150              |
|                                                                      | 13.180              |
|                                                                      | 16.430              |
| 50.000                                                               | 35.990              |
|                                                                      | 44.510              |
|                                                                      | 23.090              |
|                                                                      | 38.060              |

Table II-102. DPA on Brick - Analysis of Target-  
Found Concentration Points

DIPHENYLAMINE (DPA)  
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 13.246 SD= 13.1338996494

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.41527710356  
SLOPE= 0.693552588997

USE FOR ACCURACY

R= 0.952364373267

MEAN SQR DEV OF POINTS FROM REGRESSION= 16.9340657172  
ST ERROR EST= 4.11510215148

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE

MEASURED 1 TIME(S))

Y(C)= 7.91768779587

X(D)= 21.366954558






Table 11-103. DPA on Brick - Inaccuracy and Imprecision Data

DIPHENYLAMINE (DPA)

BRICK SURFACE

STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION

| Mn Target Con<br>ng/ml | Mn Found Conc<br>ng/ml | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------|------------------------|-----------------------|------------------------|-------------|
| 2.500                  | 2.003                  | 0.233                 | -19.900                | 11.636      |
| 5.000                  | 4.128                  | 1.256                 | -17.450                | 36.422      |
| 10.000                 | 7.778                  | 0.508                 | -22.225                | 6.536       |
| 25.000                 | 16.910                 | 4.198                 | -32.360                | 24.828      |
| 50.000                 | 35.413                 | 8.981                 | -29.175                | 25.360      |
| Means                  |                        | 3.835                 | -24.222                | 19.756      |

DIPHENYLAMINE (DPA)  
 BRICK SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-101. DPA on Brick - Graph of Inaccuracy

DIPHENYLAMINE (DPA)  
BRICK SURFACE  
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)

30.42




Figure II-102. DPA on Brick - Graph of Imprecision

Table II-104. DPA on Transite - Target vs. Found Concentrations

| DIPHENYLAMINE (DPA)<br>TRANSITE SURFACE<br>TARGET CONC. VS. FOUND CONC. |                     |
|-------------------------------------------------------------------------|---------------------|
| Target Conc<br>ng/ml                                                    | Found Conc<br>ng/ml |
| 2.500                                                                   | 0.810               |
|                                                                         | 1.490               |
|                                                                         | 1.880               |
|                                                                         | 1.740               |
| 5.000                                                                   | 1.320               |
|                                                                         | 4.640               |
|                                                                         | 4.740               |
|                                                                         | 3.910               |
| 10.000                                                                  | 2.840               |
|                                                                         | 9.790               |
|                                                                         | 9.030               |
|                                                                         | 7.970               |
| 25.000                                                                  | 7.550               |
|                                                                         | 18.970              |
|                                                                         | 21.120              |
|                                                                         | 17.340              |
| 50.000                                                                  | 11.650              |
|                                                                         | 45.170              |
|                                                                         | 37.090              |
|                                                                         | 41.580              |

Table II-105. DPA on Transite - Analysis of Target-  
Found Concentration Points

DIPHENYLAMINE (DPA)  
TRANSITE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.0350535872

FOUND CONC  
MEAN= 12.5315 SD= 13.8086106908

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.11885315534  
SLOPE= 0.670953883495  
USE FOR ACCURACY  
R= 0.876314751306  
MEAN SQR DEV OF POINTS FROM REGRESSION= 46.7094406587  
ST ERROR EST= 6.83443052922

USE FOR PRECISION  
T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE  
MEASURED 1 TIME(S))

y(c)= 12.5789824038

x(d)= 37.1490430551

DIPHENYLAMINE (DPA)  
TRANSITE SURFACE  
FOUND CONC

45.17




Figure 11-103. DPA on Transite - Graph of Target-Found Concentration Points


Table II-106. DPA on Transite - Inaccuracy and Imprecision Data

| DIPHENYLAMINE (DPA)<br>TRANSITE SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                        |                       |                        |             |  |
|---------------------------------------------------------------------------------------------------------------------|------------------------|-----------------------|------------------------|-------------|--|
| Mn Target Con<br>ng/ml                                                                                              | Mn Found Conc<br>ng/ml | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |  |
| 2.500                                                                                                               | 1.480                  | 0.475                 | -40.800                | 32.080      |  |
| 5.000                                                                                                               | 3.653                  | 1.590                 | -26.950                | 43.762      |  |
| 10.000                                                                                                              | 7.408                  | 3.135                 | -25.925                | 42.324      |  |
| 25.000                                                                                                              | 16.245                 | 6.000                 | -35.020                | 36.933      |  |
| 50.000                                                                                                              | 33.873                 | 15.179                | -32.255                | 44.813      |  |
| Means                                                                                                               |                        | 5.278                 | -32.190                | 39.984      |  |



DIPHENYLAMINE (DPA)  
 TRANSITE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-25.93



VERTICAL AXIS TIC INTERVAL = 1.4875  
 -40.80

Figure 11-104. DPA on Transite - Graph of Inaccuracy

DIPHENYLAMINE (DPA)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-105. DPA on Transite - Graph of Imprecision

Table II-107. DPA on Transite (3 days) - Target vs. Found Concentrations

| DIPHENYLAMINE (DPA)<br>TRANSITE SURFACE<br>TARGET CONC. VS. FOUND CONC. |                     |
|-------------------------------------------------------------------------|---------------------|
| Target Conc<br>ng/ml                                                    | Found Conc<br>ng/ml |
| 2.500                                                                   | 1.490               |
|                                                                         | 1.880               |
|                                                                         | 1.740               |
| 5.000                                                                   | 4.640               |
|                                                                         | 4.740               |
|                                                                         | 3.910               |
| 10.000                                                                  | 9.790               |
|                                                                         | 9.030               |
|                                                                         | 7.970               |
| 25.000                                                                  | 18.970              |
|                                                                         | 21.120              |
|                                                                         | 17.340              |
| 50.000                                                                  | 45.170              |
|                                                                         | 37.090              |
|                                                                         | 41.580              |

Table II-108. DPA on Transite (3 days) - Target vs.  
Found Concentrations

DIPHENYLAMINE (DPA)

TRANSITE SURFACE

ANALYSIS OF 15 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 18.5 SD= 18.1953683274

FOUND CONC  
MEAN= 15.097333333 SD= 14.9797851723

N0. RUNS 1 TOTAL X-Y ALL RUNS 15 N0. CONCENTR 15  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -0.011459007551  
SLOPE= 0.816691477886  
USE FOR ACCURACY  
R= 0.992803695584  
MEAN SQR DEV OF POINTS FROM REGRESSION= 3.84924285745  
ST ERROR EST= 1.96194873976  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 13  
TWO TAIL P LEVEL IS .1  
t= 1.77093178942  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
(EACH TARGET CONC CONSIDERED INDEP SAMPLE  
MEASURED 1 TIME(S))  
y(c)= 3.69909237891  
x(d)= 8.97736979202



Table II-109. DPA on Transite (3 days) - Inaccuracy and Imprecision Data


**DIPHENYLAMINE (DPA)**

**TRANSITE SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION**

| Mn Target Con<br>ng/ml | Mn Found Conc<br>ng/ml | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------|------------------------|-----------------------|------------------------|-------------|
| 2.500                  | 1.703                  | 0.198                 | -31.867                | 11.599      |
| 5.000                  | 4.430                  | 0.453                 | -11.400                | 10.228      |
| 10.000                 | 8.930                  | 0.914                 | -10.700                | 10.236      |
| 25.000                 | 19.143                 | 1.896                 | -23.427                | 9.984       |
| 50.000                 | 41.280                 | 4.848                 | -17.440                | 9.807       |
| Means                  |                        | 1.502                 | -18.967                | 10.355      |

DIPHENYLAMINE (DPA)  
 TRANSITE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)

-10.70



-31.87

TARGET CONC

VERTICAL AXIS TIC INTERVAL = 2.11666666667

Figure 11-107. DPA on Transite (3 days) - Graph of Inaccuracy




Figure 11-108. DPA on Transite (3 days) - Graph of Imprecision



Table II-110. 2,6-DNT on Metal - Target vs. Found Concentrations

| 2,6-DINITROTOLUENE (26DNT)<br>METAL SURFACE |                                              |
|---------------------------------------------|----------------------------------------------|
| TARGET CONC.<br>ug/10 sq cm                 | VS. FOUND CONC.<br>Found Conc<br>ug/10 sq cm |
| 10.000                                      | 9.550                                        |
|                                             | 7.500                                        |
|                                             | 11.350                                       |
|                                             | 10.500                                       |
| 20.000                                      | 19.000                                       |
|                                             | 21.000                                       |
|                                             | 22.500                                       |
|                                             | 19.300                                       |
| 40.000                                      | 29.000                                       |
|                                             | 39.900                                       |
|                                             | 40.000                                       |
|                                             | 37.100                                       |
| 100.000                                     | 86.700                                       |
|                                             | 97.000                                       |
|                                             | 95.900                                       |
|                                             | 93.600                                       |
| 200.000                                     | 168.000                                      |
|                                             | 208.000                                      |
|                                             | 192.000                                      |
|                                             | 189.000                                      |

Table II-111 - 2,6-DNT on Metal - Analysis of Target-  
Found Concentration Points

2,6-DINITROTOLUENE (26DNT)  
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 74 SD= 72.140214349

FOUND CONC

MEAN= 69.845 SD= 68.3900825739

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.0850323624595

SLOPE= 0.942702265372

USE FOR ACCURACY

R= 0.994394814712

MEAN SQR DEV OF POINTS FROM REGRESSION= 55.1910252609

ST ERROR EST= 7.42906624421

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 13.6292662533

x(d)= 28.5092514656

2,6-DINITROTOLUENE (26DNT)  
METAL SURFACE  
FOUND CONC




Figure II-109. 2,6-DNT on Metal - Graph of Target-  
Found Concentration Points

TABLE II-112. 2,6-DNT on Metal - Inaccuracy and Imprecision Data

2,6-DINITROTOLUENE (26DNT)  
METAL SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Target Conc<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 10.000                        | 9.725                        | 1.656                 | -2.750                 | 17.024      |
| 20.000                        | 20.450                       | 1.626                 | 2.250                  | 7.950       |
| 40.000                        | 36.500                       | 5.178                 | -8.750                 | 14.185      |
| 100.000                       | 93.300                       | 4.622                 | -6.700                 | 4.954       |
| 200.000                       | 189.250                      | 16.439                | -5.375                 | 1.687       |
| Means                         |                              | 5.904                 | -4.265                 | 10.560      |

2,6-DINITROTOLUENE (26DNT)  
 METAL SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-110. 2,6-DNT on Metal - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)  
 METAL SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-111. 2,6-DNT on Metal - Graph of Imprecision

Table II-113. 2,6-DNT on Concrete - Target vs.  
Found Concentrations

| 2,6-DINITROTOLUENE (26DNT)  |                                             |
|-----------------------------|---------------------------------------------|
| CONCRETE                    |                                             |
| TARGET CONC.<br>ug/10 sq cm | VS FOUND CONC.<br>Found Conc<br>ug/10 sq cm |
| 10.000                      | 15.000                                      |
|                             | 10.200                                      |
|                             | 9.000                                       |
|                             | 10.200                                      |
| 20.000                      | 10.000                                      |
|                             | 16.600                                      |
|                             | 19.000                                      |
|                             | 20.400                                      |
| 40.000                      | 37.600                                      |
|                             | 35.000                                      |
|                             | 37.400                                      |
|                             | 35.700                                      |
| 100.000                     | 76.900                                      |
|                             | 94.600                                      |
|                             | 97.200                                      |
|                             | 104.000                                     |
| 200.000                     | 120.000                                     |
|                             | 189.000                                     |
|                             | 181.000                                     |
|                             | 196.000                                     |

Table II-114. 2,6-DNT on Concrete - Analysis of Target-Found Concentration Points

2,6-DINITROTOLUENE (26DNT)  
CONCRETE SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 74 SD= 72.140214349

FOUND CONC  
MEAN= 65.53 SD= 63.8466262535

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 1.8729368932

SLOPE= 0.860230582524

USE FOR ACCURACY

R= 0.971973340682

MEAN SQR DEV OF POINTS FROM REGRESSION= 237.809596818

ST ERROR EST= 15.4210763832

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 29.9877306664

x(d)= 64.5501819952



2,6-DINITROTOLUENE (26DNT)  
CONCRETE SURFACE  
FOUND CONC




Figure 11-112. 2,6-DNT on Concrete - Graph of Target-  
Found Concentration Points

Table II-115. 2,6-DNT on Concrete - Inaccuracy and Imprecision Data

**2,6-DINITROTOLUENE (26DNT)**

**CONCRETE STATISTICAL DATA USED TO DETERMINE PERCENT INACCURACY AND IMPRECISION**

| Mn Target Conc<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|-------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 10.000                        | 9.850                        | 0.574                 | -1.500                 | 5.832       |
| 20.000                        | 16.300                       | 4.609                 | -17.500                | 27.931      |
| 40.000                        | 36.625                       | 1.014                 | -8.438                 | 2.770       |
| 100.000                       | 93.175                       | 11.551                | -6.825                 | 12.397      |
| 200.000                       | 171.300                      | 34.876                | -14.250                | 20.336      |
| Means                         |                              | 10.525                | -9.702                 | 13.853      |

2,6-DINITROTOLUENE (26DNT)  
 CONCRETE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-113. 2,6-DNT on Concrete - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)  
CONCRETE  
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)

27.93




Figure 11-114. 2,6-DNT on Concrete - Graph of Imprecision

Table II-114. 2,6-DNT on Brick - Target vs.  
Found Concentrations

| 2,6-DINITROTOLUENE (26DNT)<br>BRICK SURFACE |                                             |
|---------------------------------------------|---------------------------------------------|
| TARGET CONC.<br>ug/10 sq cm                 | VS FOUND CONC.<br>Found Conc<br>ug/10 sq cm |
| 10.000                                      | 4.000<br>5.300<br>5.600<br>4.300            |
| 20.000                                      | 8.000<br>9.100<br>8.680<br>7.100            |
| 40.000                                      | 17.300<br>24.800<br>25.200<br>16.700        |
| 100.000                                     | 46.200<br>62.400<br>57.400<br>45.400        |
| 200.000                                     | 89.800<br>94.300<br>123.000<br>80.820       |

TABLE II-117. 2,6-DNT on Brick - Analysis of Target-  
Found Concentration Points

2,6-DINITROTOLUENE (26DNT)  
BRICK SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 74 SD= 72.140214349

FOUND CONC  
MEAN= 36.77 SD= 36.3865408894

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 0.463001618123

SLOPE= 0.490635113269

USE FOR ACCURACY

R= 0.972736659578

MEAN SQR DEV OF POINTS FROM REGRESSION= 75.1641619382  
ST ERROR EST= 8.66972675107

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 16.2691343172

x(d)= 63.6301961361

2,6-DINITROTOLUENE (26DNT)  
BRICK SURFACE  
FOUND CONC




Figure 11-115. 2,6-DNT on Brick - Graph of Target-Found Concentration Points

TABLE 11-118. 2,6-DNT on Brick - Inaccuracy and Imprecision Data

2,6-DINITROTOLUENE (26DNT)  
BRICK SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 10.000                       | 4.800                        | 0.770                 | -52.000                | 16.040      |
| 20.000                       | 0.220                        | 0.073                 | -50.900                | 10.626      |
| 40.000                       | 21.000                       | 4.620                 | -47.500                | 22.039      |
| 100.000                      | 52.850                       | 8.399                 | -47.150                | 15.092      |
| 200.000                      | 96.980                       | 10.229                | -51.510                | 10.797      |
| Means                        |                              | 6.580                 | -51.412                | 16.680      |



2,6-DINITROTOLUENE (26DNT)  
BRICK SURFACE  
MEAN INACCURACY

-47.15

MEAN FOUND CONCENTRATION (REPORT)




Figure 11-146. 2,6-DNT on Brick - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)  
BRICK SURFACE  
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-117. 2,6-DNT on Brick - Graph of Imprecision

TABLE II-119. 2,6-DNT on Transite - Target vs.  
Found Concentrations

| 2,6-DINITROTOLUENE (26DNT)   |                           |
|------------------------------|---------------------------|
| TRANSITE SURFACE             |                           |
| TARGET CONC. VS. FOUND CONC. |                           |
| Target Conc<br>ug/10 sq cm   | Found Conc<br>ug/10 sq cm |
| 10.000                       | 7.500                     |
|                              | 8.700                     |
|                              | 7.500                     |
|                              | 7.800                     |
| 20.000                       | 14.700                    |
|                              | 16.300                    |
|                              | 17.800                    |
|                              | 19.900                    |
| 40.000                       | 27.700                    |
|                              | 36.000                    |
|                              | 36.600                    |
|                              | 33.200                    |
| 100.000                      | 89.000                    |
|                              | 87.300                    |
|                              | 81.900                    |
|                              | 89.700                    |
| 200.000                      | 152.000                   |
|                              | 156.000                   |
|                              | 164.000                   |
|                              | 160.000                   |

TABLE II-120. 2,6-DNT on Transite - Analysis of  
Target-Found Concentration Points

2,6-DINITROTOLUENE (26DNT)  
TRANSITE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 74 SD= 72.140214349

FOUND CONC  
MEAN= 60.68 SD= 57.3720507785

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 1.99743122977

SLOPE= 0.793007686084

USE FOR ACCURACY

R= 0.99713612601

MEAN SQR DEV OF POINTS FROM REGRESSION= 19.8720941201  
ST ERROR EST= 4.45781158418

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 10.1246494767

X(D)= 20.3745583959

2,6-DINITROTOLUENE (26DNT)  
TRANSITE SURFACE  
FOUND CONC





Figure 11-118. 2,6-DNT on Transite - Graph of Target-Found Concentration Points

TABLE II-121. 2,6-DNT on Transite - Inaccuracy and Imprecision Data

2,6-DINITROTOLUENE (26DNT)  
 TRANSITE SURFACE  
 STATISTICAL DATA USED TO DETERMINE PERCENT  
 INACCURACY AND IMPRECISION

| Mn Target Con ug/10 sq cm | Mn Found Conc ug/10 sq cm | Standard Deviation | Mean Pct Inaccuracy | Imprecision |
|---------------------------|---------------------------|--------------------|---------------------|-------------|
| 10.000                    | 7.875                     | 0.560              | -21.250             | 7.211       |
| 20.000                    | 17.175                    | 2.214              | -14.125             | 12.892      |
| 40.000                    | 33.375                    | 4.063              | -16.563             | 12.174      |
| 100.000                   | 86.975                    | 3.530              | -13.025             | 4.059       |
| 200.000                   | 158.000                   | 5.164              | -21.000             | 3.260       |
| Means                     |                           | 3.108              | -17.193             | 7.921       |

2,6-DINITROTOLUENE (26DNT)  
 TRANSITE SURFACE  
 MEAN FOUND CONCENTRATION (REPORT)  
 MEAN INACCURACY  
 -13.03



VERTICAL AXIS TIC INTERVAL = 0.8225

Figure II-119. 2,6-DNT on Transite - Graph of Inaccuracy

2,6-DINITROTOLUENE (26DNT)  
 TRANSITE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-120. 2,6-DNT on Transite - Graph of Imprecision



TABLE II-122. NG on Metal - Target vs. Found Concentrations

| NITROGLYCERIN<br>METAL SURFACE<br>TARGET CONC. VS. FOUND CONC. |                                              |
|----------------------------------------------------------------|----------------------------------------------|
| Target Conc<br>ug/10 sq cm                                     | Found Conc<br>ug/10 sq cm                    |
| 125.000                                                        | 122.000<br>99.500<br>102.000<br>119.000      |
| 250.000                                                        | 205.000<br>210.000<br>206.000<br>220.000     |
| 500.000                                                        | 415.000<br>423.000<br>402.000<br>443.000     |
| 1250.000                                                       | 1150.000<br>1200.000<br>1150.000<br>1140.000 |

TABLE II-123. NG on Metal - Analysis of Target-  
Found Concentration Points

NITROGLYCERINE (NG)  
METAL SURFACE

ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC

MEAN= 531.25 SD= 450.693909433

FOUND CONC

MEAN= 475.40625 SD= 424.537778308

NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -24.3256419256

SLOPE= 0.940671794872

USE FOR ACCURACY

R= 0.998580335082

MEAN SQR DEV OF POINTS FROM REGRESSION= 547.954235351  
ST ERROR EST= 23.4084223166

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 19.9865337733

x(d)= 93.6030226134




Figure 11-121. NG on Metal - Graph of Target-Found Concentration Points

TABLE II-124. NG on Metal - Inaccuracy and Imprecision Data

| NITROGLYCERIN (NG)<br>METAL SURFACE<br>STATISTICAL DATA USED TO DETERMINE PERCENT<br>INACCURACY AND IMPRECISION |                              |                       |                        |             |
|-----------------------------------------------------------------------------------------------------------------|------------------------------|-----------------------|------------------------|-------------|
| Mn Target Con<br>ug/10 sq cm                                                                                    | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
| 125.000                                                                                                         | 110.625                      | 11.514                | -11.500                | 10.408      |
| 250.000                                                                                                         | 210.250                      | 6.050                 | -15.900                | 3.258       |
| 500.000                                                                                                         | 420.750                      | 17.173                | -15.050                | 4.082       |
| 1250.000                                                                                                        | 1160.000                     | 27.000                | -7.200                 | 2.334       |
| Means                                                                                                           |                              | 15.654                | -12.613                | 5.020       |

NITROGLYCERIN (NG)  
 METAL SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure 11-122. NG on Metal - Graph of Inaccuracy

NITROGLYCERIN (NG)  
METAL SURFACE

IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-123. NG on Metal - Graph of Imprecision

TABLE IT-125. NG on Concrete Target vs. Found Concentrations

| NITROGLYCERIN (NG)<br>CONCRETE SURFACE<br>TARGET CONC. VS FOUND CONC. |                           |
|-----------------------------------------------------------------------|---------------------------|
| Target Conc<br>ug/10 sq cm                                            | Found Conc<br>ug/10 sq cm |
| 125.000                                                               | 93.300                    |
|                                                                       | 86.600                    |
|                                                                       | 95.900                    |
|                                                                       | 99.500                    |
| 250.000                                                               | 181.000                   |
|                                                                       | 142.000                   |
|                                                                       | 193.000                   |
|                                                                       | 196.000                   |
| 500.000                                                               | 396.000                   |
|                                                                       | 392.000                   |
|                                                                       | 436.000                   |
|                                                                       | 429.000                   |
| 1250.000                                                              | 866.000                   |
|                                                                       | 883.000                   |
|                                                                       | 1007.000                  |
|                                                                       | 1014.000                  |

TABLE II-126. NG on Concrete - Analysis of Target-  
Found Concentration Points

**NITROGLYCERINE (NG)**

**CONCRETE SURFACE**

**ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS**

**TARGET CONC**

MEAN= 531.25 SD= 450.693909433

**FOUND CONC**

MEAN= 406.89375 SD= 343.646924364

NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 4.94128205129

SLOPE= 0.756616410256

USE FOR ACCURACY

R= 0.992304553607

MEAN SQR DEV OF POINTS FROM REGRESSION= 1939.89261997

ST ERROR EST= 44.0442121052

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 88.3170343321

X(C)= 217.473873833



NITROGLYCERINE (NG)  
CONCRETE SURFACE  
FOUND CONC

1014.00




Figure II-124. NG on Concrete - Graph of Target-Found Concentration Points

TABLE II-127. NG on Concrete - Inaccuracy and Imprecision Data

**NITROGLYCERIN (NG)  
CONCRETE SURFACE  
STATISTICAL DATA USED TO DETERMINE PERCENT  
INACCURACY AND IMPRECISION**

| Mn Target Conc<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Fct<br>Inaccuracy | Imprecision |
|-------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 125.000                       | 93.825                       | 5.446                 | -24.940                | 5.885       |
| 250.000                       | 178.000                      | 24.860                | -28.800                | 13.966      |
| 500.000                       | 413.250                      | 22.470                | -17.350                | 5.437       |
| 1250.000                      | 942.500                      | 78.878                | -24.600                | 8.369       |
| Means                         |                              | 32.913                | -23.923                | 8.394       |

NITROGLYCERIN (NG)  
 CONCRETE SURFACE  
 MEAN INACCURACY MEAN FOUND CONCENTRATION (REPORT)




Figure II-125. NG on Concrete - Graph of Inaccuracy

NITROGLYCERIN (NG)  
 CONCRETE SURFACE  
 IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure II-126. NG on Concrete - Graph of Imprecision

TABLE II-128. NG on Brick - Target vs. Found Concentrations

| NITROGLYCERIN<br>BRICK SURFACE<br>TARGET CONC. VS. FOUND CONC. |                           |
|----------------------------------------------------------------|---------------------------|
| Target Conc<br>ug/10 sq cm                                     | Found Conc<br>ug/10 sq cm |
| 125.000                                                        | 33.800                    |
|                                                                | 56.600                    |
|                                                                | 58.800                    |
|                                                                | 42.800                    |
| 250.000                                                        | 85.900                    |
|                                                                | 87.500                    |
|                                                                | 88.800                    |
|                                                                | 77.400                    |
| 500.000                                                        | 222.000                   |
|                                                                | 282.000                   |
|                                                                | 242.000                   |
|                                                                | 202.000                   |
| 1250.000                                                       | 537.000                   |
|                                                                | 614.000                   |
|                                                                | 560.000                   |
|                                                                | 453.000                   |

TABLE II-129. NG on Brick - Analysis of Target-  
Found Concentration Points

**NITROGLYCERINE (NG)**  
**BRICK SURFACE**  
**ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS**

**TARGET CONC**  
**MEAN= 531.25 SD= 450.693909433**


**FOUND CONC**  
**MEAN= 227.725 SD= 203.495737875**

**NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16**  
**MEASURES (Y'S) EACH TARGET CONC 1**

**INTERCEPT= -8.04974358974**  
**SLOPE= 0.443811282051**  
**USE FOR ACCURACY**  
**R= 0.982934796801**  
**MEAN SQR DEV OF POINTS FROM REGRESSION= 1501.39082051**  
**ST ERROR EST= 38.7477847175**  
**USE FOR PRECISION**  
**T FOR CONFIDENCE BAND**  
**D.F. = 14**  
**TWO TAIL P LEVEL IS .1**  
**t= 1.7613101065**  
**X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C**  
**Y(C)= 65.29986478**  
**X(D)= 324.819610735**

NITROGLYCERINE (NG)  
BRICK SURFACE  
FOUND CONC

614.00



VERTICAL AXIS TIC INTERVAL= 61.4  
TARGET CONC

Figure II-127. NG on Brick - Graph of Target-Pound Concentration Points

TABLE II-130. NG on Brick - Inaccuracy and Imprecision Data

| NITROGLYCERIN (NG)                         |                              |                       |                        |             |  |
|--------------------------------------------|------------------------------|-----------------------|------------------------|-------------|--|
| BRICK SURFACE                              |                              |                       |                        |             |  |
| STATISTICAL DATA USED TO DETERMINE PERCENT |                              |                       |                        |             |  |
| INACCURACY AND IMPRECISION                 |                              |                       |                        |             |  |
| Mn Target Con<br>ug/10 sq cm               | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |  |
| 125.000                                    | 40.000                       | 11.822                | -61.600                | 24.629      |  |
| 250.000                                    | 84.900                       | 5.139                 | -66.040                | 6.053       |  |
| 500.000                                    | 237.000                      | 34.157                | -52.600                | 14.412      |  |
| 1250.000                                   | 541.000                      | 66.958                | -56.720                | 12.377      |  |
| Means                                      |                              | 29.519                | -59.240                | 14.368      |  |






Figure 11-128. NG on Brick - Graph of Inaccuracy




Figure II-129. NC on Brick - Graph of Imprecision

TABLE II-131. NC on Transite - Target vs. Found Concentrations

| NITROGLYCERIN (NG)<br>TRANSITE SURFACE<br>TARGET CONC. VS FOUND CONC. |                           |
|-----------------------------------------------------------------------|---------------------------|
| Target Conc<br>ug/10 sq cm                                            | Found Conc<br>ug/10 sq cm |
| 125.000                                                               | 83.600                    |
|                                                                       | 76.000                    |
|                                                                       | 91.650                    |
|                                                                       | 85.500                    |
| 250.000                                                               | 164.000                   |
|                                                                       | 138.000                   |
|                                                                       | 174.000                   |
|                                                                       | 138.000                   |
| 500.000                                                               | 325.000                   |
|                                                                       | 442.000                   |
|                                                                       | 396.000                   |
|                                                                       | 369.000                   |
| 1250.000                                                              | 938.000                   |
|                                                                       | 952.000                   |
|                                                                       | 764.000                   |
|                                                                       | 902.000                   |

TABLE II-132. NG on Transite - Analysis of Target-  
Found Concentration Points

NITROGLYCERINE (NG)  
TRANSITE SURFACE

ANALYSIS OF 16 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 531.25 SD= 450.693909433

FOUND CONC  
MEAN= 377.421875 SD= 328.83558048

NO. RUNS 1 TOTAL X-Y ALL RUNS 16 NO. CONCENTR 16  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -3.90769230769  
SLOPE= 0.721561538462

USE FOR ACCURACY

R= 0.988954389276

MEAN SQR DEV OF POINTS FROM REGRESSION= 2545.27895262

ST ERROR EST= 50.4507577011

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 14

TWO TAIL P LEVEL IS .1

t= 1.7613101065

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

Y(C)= 89.5956595952

X(D)= 260.722616729

NITROGLYCERINE (NG)  
TRANSITE SURFACE  
FOUND CONC

952.00




Figure 11-130. NG on Transite - Graph of Target-Found Concentration Points

TABLE II-133. NG on Transite - Inaccuracy and Imprecision Data

**NITROGLYCERIN (NG)**

**TRANSITE SURFACE**

**STATISTICAL DATA USED TO DETERMINE PERCENT**

**INACCURACY AND IMPRECISION**

| Mn Target Con<br>ug/10 sq cm | Mn Found Conc<br>ug/10 sq cm | Standard<br>Deviation | Mean Pct<br>Inaccuracy | Imprecision |
|------------------------------|------------------------------|-----------------------|------------------------|-------------|
| 125.000                      | 84.188                       | 6.450                 | -32.650                | 7.661       |
| 250.000                      | 153.500                      | 18.358                | -38.600                | 11.959      |
| 500.000                      | 383.000                      | 49.024                | -23.400                | 12.800      |
| 1250.000                     | 889.000                      | 85.953                | -28.880                | 9.669       |
| Means                        |                              | 39.946                | -30.883                | 10.522      |




Figure II-131. NG on Transite Graph of Inaccuracy

NITROGLYCERIN (NG)  
TRANSITE SURFACE  
IMPRECISION MEAN FOUND CONCENTRATION (REPORT)




Figure 11-132. NG on Transite - Graph of Imprecision



TABLE II-134. PETN on Metal - Target vs. Found Concentrations

PENTERYTHRITE TETRAHITRATE (PETN)

| PENTERYTHRITE TETRAHITRATE<br>METAL SURFACE |                                             |
|---------------------------------------------|---------------------------------------------|
| TARGET CONC.<br>UG/ 10 SQ CM                | US FOUND CONC<br>Found Conc<br>UG/ 10 SQ CM |
| 50.000                                      | 30.200                                      |
| 50.000                                      | 32.400                                      |
| 50.000                                      | 47.400                                      |
| 50.000                                      | 47.800                                      |
| 100.000                                     | 78.000                                      |
| 100.000                                     | 73.500                                      |
| 100.000                                     | 90.800                                      |
| 100.000                                     | 79.600                                      |
| 200.000                                     | 199.000                                     |
| 200.000                                     | 159.000                                     |
| 200.000                                     | 187.000                                     |
| 200.000                                     | 189.000                                     |
| 500.000                                     | 450.000                                     |
| 500.000                                     | 406.000                                     |

Table 11-134. PETN on Metal - Target vs. Found Concentrations (Continued)

PENTAERYTHRITE TETRANITRATE (PETN)

| METAL SURFACE | TARGET CONC. UG/ 10 SQ CM | US FOUND CONC Found Conc UG/ 10 SQ CM |
|---------------|---------------------------|---------------------------------------|
|               | 349.000                   | 308.000                               |
|               | 349.000                   | 363.000                               |
|               | 1000.000                  | 902.000                               |
|               | 1000.000                  | 859.000                               |
|               | 697.000                   | 612.000                               |
|               | 1388.000                  | 1150.000                              |

TABLE II-135. PETN on Metal - Analysis of Target-  
Found Concentration Points

PENTAERYTHRITE TETRANITRATE (PETN)  
METAL SURFACE

ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 359.15 SD= 382.937438205

FOUND CONC  
MEAN= 313.185 SD= 328.387893672

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 5.97887894415

SLOPE= 0.855369959782

USE FOR ACCURACY

R= 0.997458150645

MEAN SQR DEV OF POINTS FROM REGRESSION= 577.940153813

ST ERROR EST= 24.0403858915

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

THO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 49.6274158122

x(d)= 101.531895278

PENTAERYTHRITE TETRANITRATE (PETN)  
METAL SURFACE  
FOUND CONC




Figure 11-133. PETN on Metal - Graph of Target-Found Concentration Points

TABLE II-136 PETN on Concrete - Target vs. Found Concentrations

(PETN)

| PENTAERYTHRITE TETRANITRATE<br>CONCRETE SURFACE |                                              |
|-------------------------------------------------|----------------------------------------------|
| TARGET CONC.<br>UG/10 SQ CM                     | VS. FOUND CONC.<br>Found Conc<br>UG/10 SQ CM |
| 50.000                                          | 30.600                                       |
| 50.000                                          | 44.400                                       |
| 50.000                                          | 45.800                                       |
| 50.000                                          | 44.200                                       |
| 100.000                                         | 97.000                                       |
| 100.000                                         | 86.100                                       |
| 100.000                                         | 84.700                                       |
| 100.000                                         | 86.900                                       |
| 200.000                                         | 167.000                                      |
| 200.000                                         | 194.000                                      |
| 200.000                                         | 181.000                                      |
| 200.000                                         | 182.000                                      |
| 349.000                                         | 323.000                                      |
| 349.000                                         | 292.000                                      |

Table II-136. PETN on Concrete - Target vs. Found Concentrations (Continued)

PENTERYTHRITE TETRANITRATE (PETN)

| PENTERYTHRITE TETRANITRATE<br>CONCRETE SURFACE<br>TARGET CONC. VS. FOUND CONC. |                            |
|--------------------------------------------------------------------------------|----------------------------|
| Target Conc<br>UG/ 10 SQ CM                                                    | Found Conc<br>UG/ 10 SQ CM |
| 349.000                                                                        | 300.000                    |
| 349.000                                                                        | 304.000                    |
| 697.000                                                                        | 583.000                    |
| 697.000                                                                        | 610.000                    |
| 697.000                                                                        | 603.000                    |
| 1388.000                                                                       | 1140.000                   |

TABLE II-137 PEIN on Concrete - Analysis of  
Target-Found Concentration Points

PENTAERYTHRITE TETRANITRATE (PETN)  
CONCRETE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 313.75 SD= 333.434578922

FOUND CONC  
MEAN= 270.335 SD= 276.82953766

N0. RUNS 1 TOTAL X-Y ALL RUNS 20 N0. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= 10.0941892935

SLOPE= 0.829452783128

USE FOR ACCURACY

R= 0.999056104403

MEAN SQR DEV OF POINTS FROM REGRESSION= 152.635268141

ST ERROR EST= 12.3545646682

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 32.5286725913

x(d)= 53.9135179075

PENTAERYTHRITE TETRANITRATE (PETN)  
 CONCRETE SURFACE  
 FOUND CONC

1149.00




Figure II-134. PETN on Concrete - Graph of Target-  
 Found Concentration Points



TABLE II-138 PETN on Brick - Target vs. Found Concentrations

(PETN)

| PENTAERYTHRITE TETRANITRATE<br>BRICK SURFACE<br>TARGET CONC. VS. FOUND CONC.<br>Target Conc Found Conc<br>UG/ 10 SQ CM UG/ 10 SQ CM |         |
|-------------------------------------------------------------------------------------------------------------------------------------|---------|
| 50.000                                                                                                                              | 39.700  |
| 50.000                                                                                                                              | 25.000  |
| 50.000                                                                                                                              | 26.400  |
| 50.000                                                                                                                              | 22.200  |
| 100.000                                                                                                                             | 70.400  |
| 100.000                                                                                                                             | 47.200  |
| 100.000                                                                                                                             | 52.800  |
| 100.000                                                                                                                             | 47.200  |
| 200.000                                                                                                                             | 112.000 |
| 200.000                                                                                                                             | 100.000 |
| 200.000                                                                                                                             | 97.100  |
| 200.000                                                                                                                             | 101.000 |
| 500.000                                                                                                                             | 364.000 |
| 349.000                                                                                                                             | 184.000 |

Table II-138. PETN on Brick - Target vs. Found Concentrations (Continued)

PENTAERYTHRITE TETRANITRATE (PETN)

| BRICK SURFACE             |                                        |
|---------------------------|----------------------------------------|
| TARGET CONC. UG/ 10 SQ CM | US FOUND CONC. Found Conc UG/ 10 SQ CM |
| 349.000                   | 168.000                                |
| 349.000                   | 180.000                                |
| 1000.000                  | 677.000                                |
| 697.000                   | 425.000                                |
| 697.000                   | 357.000                                |
| 697.000                   | 345.000                                |

TABLE II-139 PETN on Brick - Analysis of Target-  
Found Concentration Point

PENTAERYTHRITE TETRANITRATE (PETN)  
BRICK SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 301.9 SD= 277.200079745

FOUND CONC  
MEAN= 172.05 SD= 174.438034542

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -14.2223284125  
SLOPE= 0.617000094112  
USE FOR ACCURACY  
R= 0.980476968454  
MEAN SQR DEV OF POINTS FROM REGRESSION= 1241.08232912  
ST ERROR EST= 35.2403537031  
USE FOR PRECISION  
T FOR CONFIDENCE BAND  
D.F.= 18  
TWO TAIL P LEVEL IS .1  
t= 1.73406096408  
X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C  
Y(C)= 50.2303186945  
X(D)= 206.251664127

677.80




Figure 11-135. PETM on Brick - Graph of Target-Found Concentration Points

TABLE 11-140 PETN on Transite - Target vs. Found Concentrations

PENTERYTHRITE TETRANITRATE (PETN)

| TRANSITE SURFACE<br>TARGET CONC. US FOUND CONC.<br>Target Conc Found Conc<br>UG/ 10 SQ CM UG/ 10 SQ CM |         |
|--------------------------------------------------------------------------------------------------------|---------|
| 50.000                                                                                                 | 30.400  |
| 50.000                                                                                                 | 38.600  |
| 50.000                                                                                                 | 38.100  |
| 50.000                                                                                                 | 48.800  |
| 100.000                                                                                                | 61.300  |
| 100.000                                                                                                | 88.900  |
| 100.000                                                                                                | 77.200  |
| 100.000                                                                                                | 88.900  |
| 200.000                                                                                                | 110.000 |
| 200.000                                                                                                | 167.000 |
| 200.000                                                                                                | 152.000 |
| 200.000                                                                                                | 164.000 |
| 500.000                                                                                                | 360.000 |
| 349.000                                                                                                | 239.000 |

Table II-140. PETN on Transit - Target vs. Found Concentrations (Continued)

PENTAERYTHRITE TETRANITRATE (PETN)

| TRANSIT SURFACE<br>TARGET CONC. VS. FOUND CONC. |                            |
|-------------------------------------------------|----------------------------|
| Target Conc<br>UG/ 10 SQ CM                     | Found Conc<br>UG/ 10 SQ CM |
| 349.000                                         | 288.000                    |
| 349.000                                         | 299.000                    |
| 1000.000                                        | 801.000                    |
| 697.000                                         | 525.000                    |
| 1000.000                                        | 804.000                    |
| 1388.000                                        | 1158.000                   |

TABLE 11-141 PETN on Transite - Analysis of Target-  
Found Concentration Points

PENTAERYTHRITE TETRANITRATE (PETN)  
TRANSITE SURFACE  
ANALYSIS OF 20 TARGET CONC-FOUND CONC POINTS

TARGET CONC  
MEAN= 351.6 SD= 381.500141409

FOUND CONC  
MEAN= 276.96 SD= 312.148352719

NO. RUNS 1 TOTAL X-Y ALL RUNS 20 NO. CONCENTR 20  
MEASURES (Y'S) EACH TARGET CONC 1

INTERCEPT= -10.0036586066

SLOPE= 0.816165104114

USE FOR ACCURACY

R= 0.997497183375

MEAN SQR DEV OF POINTS FROM REGRESSION= 514.183809356

ST ERROR EST= 22.6756214767

USE FOR PRECISION

T FOR CONFIDENCE BAND

D.F.= 18

TWO TAIL P LEVEL IS .1

t= 1.73406096408

X(D) FOR CALIBRATION CURVE OR UNKNOWN SAMPLE? C/U C

y(c)= 31.1370879006

x(d)= 100.308711628

PENTAERYTHRITE TETRANITRATE (PETN)  
 TRANSITE SURFACE  
 FOUND CONC




Figure II-136. PETN on Transite - Graph of Target-Found Concentration Points