

Optically-Based Diagnostics for Gas-Phase Laser Development

**W. Terry Rawlins, Seonkyung Lee,
Kristin L. Galbally-Kinney, William J. Kessler,
Adam J. Hicks, Ian M. Konen, Emily P. Plumb, and Steven J. Davis**
Physical Sciences Inc.
Andover, MA 01810

**XVIII International Symposium on Gas Flow and Chemical Lasers and
High Power Lasers 2010**
Sofia, Bulgaria

August 2010

Acknowledgement of Support and Disclaimer

This material is based upon work supported by Air Force Office of Scientific Research under Contract Number FA9550-07-1-0575. Any opinions, findings and conclusions or recommendations expressed in this material are those of the author(s) and do not necessarily reflect the views of Air Force Office of Scientific Research.

Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE AUG 2010		2. REPORT TYPE		3. DATES COVERED 00-00-2010 to 00-00-2010	
4. TITLE AND SUBTITLE Optically-Based Diagnostics for Gas-Phase Laser Development				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Physical Sciences Inc,20 New England Business Center,Andover,MA,01810				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution unlimited					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT Same as Report (SAR)	18. NUMBER OF PAGES 28	19a. NAME OF RESPONSIBLE PERSON
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified			

Overview


VG10-174-1

- **Theme: multi-species diagnostics for **absolute** concentrations are essential for effective development of high-energy gas lasers**
 - Precursor production, loss, optimization
 - Transference from subscale reactors
 - Scaling of gain, power, efficiency
- **Current applications**
 - Electric Oxygen Iodine Laser (EOIL): precursor kinetics and gain dynamics
 - Related systems: COIL, micro-COIL
 - Alkali laser systems: DPAL, XPAL gain, multi-photon effects
- **Outline of presentation:**
 - Brief overview of diagnostics and apparatus
 - Absolute emission spectrometry
 - Near-infrared spectroscopy: $\text{O}_2(a^1\Delta_g)$, $\text{I}(^2\text{P}_{1/2})$
 - Air afterglow photometry: $\text{O}(^3\text{P})$
 - Ultrasensitive absorption photometry: I_2 , O_3
 - High-resolution absorption/gain spectroscopy: atomic iodine, alkali metals

Role of Optical Diagnostics in High Energy Gas Laser Development




VG10-174-2

- **Chemically rich, energetic, reacting flow with competing phenomena**
 - Multispecies detection required
 - COIL, EOIL: $O_2(a,b)$, I^* , I , I_2 , gain, T , O , O_3
 - DPAL, XPAL: ground-state M , numerous M^* , M_2^* , MX^* , gain
- **Objective: detect key species concentrations vs. flow time**
 - Vary operating conditions systematically
 - Quantify species production and loss rates
 - Relate to system design requirements
- **Put a “cage” around the model:**




Electric Oxygen Iodine Laser (EOIL)





VG10-174-3



Ionization



Dilution in He required




$$\frac{[I^*]}{[I]} \rightarrow K_{EQ}(T) \frac{[O_2(a)]}{[O_2]}$$

Hybrid EOIL: Catalytically enhanced $O_2(a)$

PSI Microwave Discharge Flow Reactors (2450 MHz)

VG10-174-4


Low-Pressure Reactor: Active-O₂ Kinetics



- 50-100 W, 10-100 Td
- External cavity discharge
- 1-10 Torr, 1-5 mmole/s
- 5-cm i.d. Pyrex flow tube
- Flow T ~ 300-350 K

EOIL Subsonic/Supersonic Reactor

- 1-5 kW, 10-50 Td
- Coaxial MIDJet discharge
- 30-70 Torr, 40-100 mmole/s
- M ~ 2 supersonic cavity
- Lasing typically 100-150 mW
-- $M^2 = 1.08 \pm 0.01$




J-5985


Near-IR Absolute Emission: O₂(a) and I*

VG10-174-5

- InGaAs array – monochromator: observe entire band
- **Concentration = Intensity ÷ Einstein coefficient**





- $A_{00} = 2.20 \times 10^{-4} \text{ s}^{-1} (\pm 10\%)$
- Detection limit $\sim 5 \times 10^{12} \text{ cm}^{-3}$ (5 cm path)



- $A = 8.0 \text{ s}^{-1} (\pm 20\%)$
- Detection limit $\sim 10^8 \text{ cm}^{-3}$ (5 cm path)

NIR Photometric Calibration: $O_2(a)$ and I^* Emission

VG10-174-6



- Collimated field of view: **no reflective surfaces**
 - Eliminate stray light, e.g. discharge emission
- Etendue ($A\Omega$) for calibration is identical to that for volume emission
- Spectral responsivity = (signal) / (Planck function)
- Blackbody calibrations 800-1000° C agree within 1%



Principle of Absolute Calibration Method

- **Instrumental signal for both blackbody and $O_2(a \rightarrow X)$ emission:**
 - $S(\lambda) = F(\lambda) T(\lambda) A\Omega \delta\lambda I(\lambda)$ $F = \text{spectral responsivity}$
 - Use the same $\{A\Omega \delta\lambda\}$ for calibration and gas emission measurements
- **Determine F by measuring blackbody spectrum (Planck function):**
 - $F(\lambda) = S_{BB}(\lambda) / N(\lambda, T_{BB})$ $T_{BB} \sim 1000 \text{ C}$
 - Area of source > area of fov
 - Accuracy <1%
- **Gas radiance in photons/cm²-s-sr-nm:**
 - $I_{aX}(\lambda) = S_{aX}(\lambda) / \{F(\lambda) T(\lambda)\}$
 - Correct for spectral baseline/background
- **Determine $O_2(a)$ concentration from spectrally integrated intensity:**
 - $[O_2(a)] = (4\pi/l) \int I_{aX}(\lambda) d\lambda / A_{aX}$ $A_{aX} = \text{Einstein coefficient}$

Spectral Fitting Analysis: $[O_2(a,v)]$, T_{rot}

- $O_2(a^1\Delta_g \rightarrow X^3\Sigma_g^-)$ spectroscopy:
 - Magnetic dipole transition
 - Hund's coupling case (b)
 - Bose-Einstein statistics ($^{16}O_2$)
 - → 9 rotational branches
- Our procedure: determine line strengths for (0,0) from line-by-line compilation (HITRAN)
 - Boltzmann rotational temperature
 - Shift (0,0) envelope to band centers for (1,1), (2,2), etc.
- Convolve with instrument scan function
 - Triangular slit function for grating monochromator (0.3 nm FWHM)
- Linear least squares solutions are $\{[O_2(a,v)]\}$

VG10-174-8



Estimation of A_{11} , A_{22} , A_{33} Values

VG10-174-9

- Scale values from A_{00} via Franck-Condon factors, transition moment vs. r-centroid:

$$A_{v',v''} = (64\pi^4/3h) (v^3 q_{v',v''}) (R(r))^2$$

- Franck-Condon factors, r-centroid values from Krupenie (1972)

- Estimate scaling of $(R(r))^2$ from $A_{00}:A_{01}:A_{10}$


- Literature: A_{00}/A_{01} is either ~ 50 or ~ 80
- **PSI measurement:** $A_{00}/A_{01} = 52 \pm 6$
- Badger et al. (1965): $A_{00}/A_{10} > 200$
- Solution: $(R(r))^2$ varies slightly with v'

- Solutions for $A_{00} = 2.20 \times 10^{-4} \text{ s}^{-1}$:

$$A_{11} = 2.17 \times 10^{-4} \text{ s}^{-1}$$

$$A_{22} = 2.12 \times 10^{-4} \text{ s}^{-1}$$



$$A_{33} = 2.06 \times 10^{-4} \text{ s}^{-1}$$



F-5589

BUT WE DO **NOT** OBSERVE $O_2(a, v>0)$!


VG10-174-10



- **Typical fits: $[(v=1)]/[(v=0)] \sim 3 \text{ to } 7 \%$**
 - Tends to track with temperature of discharge, i.e. thermal populations only
- **True for large range of conditions:**
 - 50 W – 2 kW discharges
 - 0.5 – 50 Torr
 - Cl_2 /BHP generators, energy pooling conditions
- **Slanger, Copeland 2003: $O_2(a, v)$ exchange with $O_2(X, v=0)$ is fast**
- **Implications for COIL I_2 dissociation mechanism?**

Spectral Fitting vs. Integration

VG10-174-11





- **Philosophy:** spectral fitting confirms band shape, T_{rot} , T_{vib} , no other radiators; then integration gives accurate values

Air Afterglow: $O + NO \rightarrow NO_2 + h\nu$

VG10-174-12



- Fiber-coupled photomultiplier
 - 580 nm filter, collimated field-of-view
- Calibrate with blackbody and titration reaction



J. Phys. Chem. 90, 320-325 (1986)

- Calibrations allow measurement of absolute emission rates for known $[O]$, $[NO]$
- $\rightarrow$ Determination of $k(580 \text{ nm})$
- Scale to other wavelengths via relative intensity measurements

Air Afterglow: Determination of [O]




VG10-174-13






For $[NO] < [N]$: $N_2(B \rightarrow A)$ emission
 For $[NO] > [N]$: $O + NO$ emission
 slope $\div [O]$ = calibration factor
 correct for $O+NO+M$ reaction

Blackbody calibration: $\rightarrow k_{580}$



Measurements: Discharge Production of O, O₂(a)


VG10-174-14



- Observed [O]/[O₂(a)] is ~ 1 or less
- Discharge models predict [O] >> [O₂(a)]
 - Electron-impact O₂ dissociation cross sections are too large
 - Possible O loss on hot walls

Ultrasensitive Dual-Beam Absorption: O₃ and I₂

VG10-174-15




$$\text{Absorbance} = \ln(I_0/I) = \sigma N \ell$$

Detection limit $\sim 10^{-5}$


$$\sigma(\text{I}_2) = 1.64 \times 10^{-18} \text{ cm}^2, 488 \text{ nm}$$

$$\sigma(\text{O}_3) = 1.14 \times 10^{-17} \text{ cm}^2, 254 \text{ nm}$$



O₃ Formation in Active-Oxygen Flow

VG10-174-16




- Surprise: O and O₃ are in ~ steady state!
- Requires O₃ conversion to O



J. Chem. Phys. 87, 5209-5221 (1987)


Small Signal Gain: Atomic Iodine

VG10-174-17




- Probe transmission on (3,4) line
- $G/\sigma(T) = [I^*] - [I]/2$
- $[I^*]$ from IR emission
 $\rightarrow [I], [I^*]/[I], ([I^*]+[I])/2$

- Scanning tunable diode laser
- Balanced ratiometric detection
- Detection limit $\sim 10^{-5} \text{ %/cm}$
- Doppler width $\rightarrow$ temperature
- Method widely used for COIL, EOIL development



Observations of I^*/I Behavior

VG10-174-18




← Subsonic flow reactor

- I^*/I ratio is limited by I^* loss reaction
- Catalytic environment enhances attainable I^*/I
- Chem. Phys. Lett. **469**, 68-70 (2009)
Proc. SPIE 7196-04 (2009)
Proc. SPIE 7581-06 (2010)

Supersonic flow reactor →

- I^*/I ratio is limited by I^* loss reaction
- Addition of NO enhances I^*/I
- I^*/I continues to increase past optimum gain point
- Proc. SPIE 6874-10 (2008)
Proc. SPIE 7581-03 (2010)
J. Appl. Phys. D **43** 025208 (2010)




J-5998

Application to Alkali and Alkali-Exciplex Systems

VG10-174-19


Cs-Ar Potential Energy Diagram



- DPAL: pump D_2 , lase D_1 (C_2H_6 promotes spin-orbit transfer)
- XPAL: pump broadband exciplex $X \rightarrow B$, lase on D_1 or D_2

DPAL/XPAL Gain Measurement Test Bed (Diode laser scanning D_1 line)

VG10-174-20




- Direct probe of population inversion dynamics
- Aids in design of optical resonators
- Portable: take to other facilities
- **Can extend to spatial imaging of gain**
 - Expect significant spatial effects in power scaling
 - Valuable tool for scaling DPAL to high powers




Optical Layout for DPAL/XPAL Gain Measurements

VG10-174-21



K-4551



K-4552

Computed D₁ Absorption Spectra: Cs



Collisional Broadening Effect

VG10-174-22

Cs $^2S_{1/2} - ^2P_{1/2}$, 894 nm

Low Pressure, Doppler broadening

High Pressure, collisional broadening





- Collisional broadening greatly expands required scan range
- High optical thickness at elevated temperatures

Absorption/Gain Spectra: $\text{Cs}(^2\text{S}_{1/2}, F''=4 \rightarrow ^2\text{P}_{1/2}, F')$, 894 nm

500 Torr Kr + 75 Torr C_2H_6 , 338 K

VG10-174-23



Pump Laser: $^2\text{S}_{1/2} \rightarrow ^2\text{P}_{3/2}$, 852 nm



- Continuing work: investigate absorption and gain dynamics for DPAL, XPAL configurations: Cs, Rb, K

State-Selected Absorption and Saturation

VG10-174-24




3-D Image of D₁ Gain, Absorption

Cs + 500 Torr Kr + 75 Torr C₂H₆

VG10-174-25

- Probe beam diameter > pump beam diameter
- Sample 9 combinations: probe{on peak, off peak, blocked} x pump{on peak, off peak, blocked}



K-4596

- Gain profile follows Gaussian profile of pump beam

Conclusions

- **Multispecies diagnostic suite**
 - Absolute emission spectrometry
 - Ultrasensitive absorption photometry
 - Scanning TDL absorption/gain spectroscopy
- **Extreme sensitivity enables subscale operation at low species concentrations**
 - Simplify chemistry, focus on primary reaction steps
 - Transfer to large scale systems: establish models for scaling
- **EOIL, catalytic EOIL, COIL, micro-COIL: operational parametrics**
 - $O_2(a)$ yield vs. small-signal gain
 - I_2 dissociation: $[I_2]$, $[I^*] + [I]$
 - O, O_3 effects
- **DPAL, XPAL: power scaling phenomena**
 - Gain vs. pump power, spatial effects at high optical depth
 - General emission spectroscopy: multi-photon effects vs. pump power

Acknowledgements

VG10-174-27

**Air Force Office of Scientific Research
High Energy Laser Joint Technology Office
Air Force Research Laboratory
Defense Advanced Research Projects Agency**

PSI

**Daniel Maser
David Oakes
Lawrence Piper**

MIT

Carol Livermore

UIUC/CUA

**David Carroll
Wayne Solomon
Joe Verdeyen
Gary Eden
Joe Zimmerman
Andrew Palla
Brian Woodard
Gabriel Benavides**

AFRL

**Tim Madden
David Hostutler**

Emory University

Michael Heaven