

REPORT DOCUMENTATION PAGE

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, gathering existing data needed, and completing and reviewing this collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a valid OMB control number. PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.

1. REPORT DATE (DD-MM-YYYY) 20-03-2001		2. REPORT TYPE Final		3. DATES COVERED (From - To) 06-01-00 to 11-30-00	
4. TITLE AND SUBTITLE Activity Monitoring Equipment				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER F49620-00-1-0322	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S) Czeisler, Charles A., Ph.D., M.D.				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Brigham and Women's Hospital 75 Francis Street Boston, MA 02115				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES) Willard D.Larkin, Ph.D. AFOSR/NL 801 North Randolph, Room 732 Arlington, VA 22203-1977					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT To evaluate the impact of sleep and waking alertness on the operational success of Air Force personnel, a system of individual monitoring is required. Monitoring is essential for evaluating the operational readiness of critical personnel, predicting individual differences in performance during sustained operations and jet-lag/shift work operations, and for use as a selection tool to identify individuals who may be better suited for particular missions. The current research was designed to implement the use of an unobtrusive ambulatory monitor (Actiwatch-L, Mini-Mitter Corporation, Bend, Oregon) to record ambient light level and wrist actigraphy of all subjects prior to the sleep deprivation jet-lag protocol; and to develop and conduct near-real-time analysis of these sleep and light exposure data. Data collected from these individuals will be analyzed using our most recent mathematical model that predicts waking alertness and neurobehavioral performance as well as circadian phase from activity (sleep-wake) and light history.					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT SAR	18. NUMBER OF PAGES 4	19a. NAME OF RESPONSIBLE PERSON Willard Larkin/AFOSR
a. REPORT U	b. ABSTRACT U	c. THIS PAGE U			19b. TELEPHONE NUMBER (include area code) (703) 696-7793

20010426 047

AIR FORCE OFFICE OF SCIENTIFIC RESEARCH (AFOSR)
NOTICE OF TRANSMITTAL DTIC. THIS TECHNICAL REPORT
HAS BEEN REVIEWED AND IS APPROVED FOR PUBLIC RELEASE
LAW AFR 190-12. DISTRIBUTION IS UNLIMITED.

PROGRESS REPORT

Submitted on 28 February, 2001

to

AFOSR

110 Duncan Avenue, Suite B 115

Bolling AFB

Washington, DC 20332-0001

by

Department of Medicine

Brigham and Women's Hospital

75 Francis Street

Boston, Massachusetts 02115

ACTIVITY MONITORING EQUIPMENT

Grant Number: F49620-00-1-0322

Performance Period: 01 Jun 00 – 30 Nov 00

Principal Investigator:

Charles A. Czeisler, Ph.D., M.D.
Professor of Medicine
Harvard Medical School
Chief, Division of Sleep Medicine
Department of Medicine
Brigham and Women's Hospital
221 Longwood Avenue
Boston, MA 02115

Telephone: (617) 732-4013
Facsimile: (617) 732-4015

Performance Sites:

Brigham and Women's Hospital
75 Francis Street (and)
221 Longwood Avenue
Boston, MA 02115

Name of Institution:

Brigham and Women's Hospital
75 Francis Street
Boston, MA 02115

2. Objectives:

Sleep is a basic biological need that is vital to the operational success of Air Force personnel. Laboratory-based studies indicate that chronically - inadequate sleep and sustained wakefulness, such as have been reported to occur in aircrew members, lead to profound decrements in human performance, including slowing of reaction time, impairment of short-term memory, degraded neurobehavioral performance, and lapses of attention, each of which could jeopardize mission success.

Given the impact of sleep and waking alertness on the operational success of Air Force personnel, a system of individual monitoring may be an essential component for: a) evaluating the operational readiness of critical personnel, b) predicting individual differences in performance during sustained operations and jet-lag/shift work operations, c) use as a selection tool to identify individuals who may be better suited for particular missions.

The timing of both the light-dark cycle, which is the most potent circadian synchronizer in the environment, and the sleep-wake cycle of personnel prior to deployment, can be monitored continuously and non-invasively on a routine operational basis using technology that has been developed over the past decade. The routine monitoring of wrist actigraphy and light exposure in personnel prior to critical missions would be a significant first step in assessing the problems associated with sleep-related neurobehavioral performance deficits in a meaningful way.

The Acitwatch-L light/activity monitors allow us to record wake-sleep activity patterns and light-exposure history in all individuals prior to admission into the laboratory for our AFOSR funded study and during the entire 31 days that they are in the laboratory our goal is to then use these subjects' previous light exposure and sleep-wake history to predict individual differences in a adjustment to shift work and SUSOPS.

The Acitwatch-L recorder (Mini-Mitter, Inc., Bend, OR) is a small wrist worn device (17 grams) that measures activity and ambient light exposure. The Acitwatch-L is made of titanium and a reinforced plastic case; is waterproof; powered by a 3V, 150 mAmp-hr Lithium Manganese battery that has a lifetime of 180 days; data is preserved if the battery expires. Once initialized, the Acitwatch-L can record data for 22 days (2 min sampling rate) and it can be programmed to become active at a designated time. Data is downloaded from the Acitwatch-L via the Acitwatch Reader. This interface reader is attached to a computer and the Acitwatch-L transfers data to the reader telemetrically (no cable connection required).

Data files in ASCII format allow raw data to be imported into spreadsheet and database software. Actiware-Sleep scoring software calculates sleep summary data such as estimates of sleep efficiency, sleep fragmentation, total sleep time and the number of wake and sleep bouts. The software also analyzes light exposure data in lux (e.g. maximum, average exposure).

The Acitwatch-L recorder senses motion (an accelerometer). The degree and speed of each motion produces an electrical current that varies in magnitude and is stored as an activity count.

The software allows for choosing bin length for both collection and analysis of data. The scoring software allows the sensitivity of the watch to be modified for analysis.

Specifically, we proposed to:

- 1) Implement the use of an unobtrusive ambulatory monitor (Actiwatch-L, Mini-Mitter Corporation, Bend, Oregon) to record ambient light level and wrist actigraphy of all subjects prior to the sleep deprivation jet-lag protocol.
- 2) Develop and conduct near-real-time analysis of these sleep and light exposure data collected from these individuals using our most recent mathematical model which predicts waking alertness and neurobehavioral performance as well as circadian phase from activity (sleep-wake) and light history.

3. Status of Effort:

We purchased 25 Actiwatch-L actigraphy and light exposure-recording units. They were delivered on 08/30/00 and were subsequently inspected and tested before use in our studies.

The equipment purchased accelerated our currently funded AFOSR research grant F49620-95-1-0388 entitled "Homeostatic and Circadian Regulation of Wakefulness During Jet Lag and Sleep Deprivation: Effect of Wake-Promoting Measures", by significantly upgrading our older equipment, and therefore allowing us to better achieve the objectives outlined within the scope of that project.

4. Accomplishments/New Findings:

The new Actiwatch-L recorders worked very effectively in measuring sleep-wake activity and light exposure patterns in research subjects. Work is ongoing to develop and optimize a transformation function so that data collected by the wrist worn Actiwatch-L can be directly input into the mathematical model of the effects of light exposure of varying duration and intensity on the human circadian system and performance. We have currently determined that the Actiwatch-L readings can be transformed to the equivalent of IL-1400 meter readings by a simple transformation function. Wrist position is a significant factor in the Actiwatch-L light readings that we are investigating further.

5. Personnel Report:

There are no budgetary costs for personnel.

6. Publications:

None

7. Interactions/Transitions:

A. Participation/presentations:

None

B. Consultative and advisory functions:

None

C. Transitions:

None

8. New discoveries, inventions, or patent disclosures:

None

9. Honors/Awards:

None